Instruction Manual
Manuale di istruzioni
Manuel d'instructions
Manual de instrucciones
Bedienungsanleitung
指导手册

AREX 6 DIGITAL PRO Heating Magnetic Stirrer

F20500460, F20510460

AREX 6 DIGITAL Heating Magnetic Stirrer

F20500461, F20510461

General Information / Informazioni Generali / Informations Générales / Información General / Allgemeine Hinweise / 一般信息

Before using the unit, please read the following instruction manual carefully.
Prima dell'utilizzo dello strumento si raccomanda di leggere attentamente il seguente manuale operativo.
Avant d'utiliser l'instrument, il est recommandé de lire attentivement le présent manuel d'instructions.
Antes de utilizar el instrumento, le recomendamos que lea con atención el siguiente manual de funcionamiento.
Bitte lesen Sie vor Inbetriebnahme des Geräts diese Bedienungsanleitung sorgfältig durch
在使用本装置之前，请仔细阅读以下使用说明书。

Caution, hot surface! / Attenzione, superficie calda! / Attention, surface chaude! / Prudencia, superficie caliente! / Vorsicht, heiße Oberfläche! / 注意，热的表面!

Do not dispose of this equipment as urban waste, in accordance with EEC directive 2002/96/CE.
Non smaltire l'apparecchiatura come rifiuto urbano, secondo quanto previsto dalla Direttiva 2002/96/CE.
Ne pas recycler l'appareil comme déchet solide urbain, conformément à la Directive 2002/96/CE.
No tirar el aparato en los desechos urbanos, como exige la Directiva 2002/96/CE.
Dieses Gerät unterliegt der Richtlinie 2002/96/EG und darf nicht mit dem normalen Hausmüll entsorgt werden.
根据 EEC 指令 2002/96/CE，请不要将本设备作为城市垃圾处理。

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This unit must be used for laboratory applications indoor only. The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions. If the product is used in a not specified way by the manufacturer or with non specified accessories, product's safety may be compromised.

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This unit has been designed and manufactured in compliance with the following standards:

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Electrical equipment for laboratory use

UL 61010-1

General requirement - Canadian electrical code

CAN/CSA-C22.2 No.61010-1

VELP reserves the right to modify the characteristics of its products with the aim to constantly improving their quality.

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The heated solution may release toxic, dangerous or poisonous gases. Adequate safety measures must be taken, in accordance with the safety regulations in force, including the presence of hood and personal protective equipment (masks, gloves, goggles, etc.).

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Do not use with explosive or dangerous materials for which the equipment is not designed. The stirrer must not be used in explosive atmospheres, in bain-marie or to stir combustible liquids that have a low combustion temperature. The product is intended for use with very small quantities of flammable liquids or flammable liquids that have a fire point higher than 625°C and a flash point higher than 600°C.

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It is responsibility of the user appropriately decontaminate the instrument in case of dangerous substances fall on or in it accordingly to the safety datasheet of substances used and to the current laboratories safety standards. It is not possible to decontaminate the product under steam.

It is also responsibility of the user to use substances for cleaning or decontaminating which do not react with internal parts of the instrument or with the material contained in it. In case of doubts on the compatibility of a cleaning solution, contact the manufacturer or local distributor.

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
1. Introduction

EN

The AREX 6 Digital/Digital PRO heating magnetic stirrer is used in the laboratory for general use and for all those applications that require a precise regulation for stirring speed and heating liquids.

- 1 Aluminum alloy plate with ceramic coating
- 2 Main switch
- 3 Speed control knob
- 4 Temperature control knob
- 5 Display

- A Probe led
- B Heating ON led
- C Temperature set point led
- D Timer led
- E Autoreverse led
- F Timer/autoreverse/lock button

The Aluminum alloy heating plate with ceramic coating (1) ensures:

- Optimum heat distribution and a high specific power thanks to the circular configuration
- Temperature homogeneity
- High resistance to thermal stress and thermal shock
- High resistance to chemical agents and surface abrasion

Magnetic stirring is generated by an AlNiCo magnet, driven by brushless motor which offers a virtually unlimited duration. Size and shape of the magnetic stir bar determines the stirring efficiency at any given speed. The stirring bar code which satisfies most applications is Ø8x40mm.

NOTE: The vessel must be made of a suitable material to withstand the foreseen temperature.

NOTE: Using the heating plate at high temperatures may cause discoloring. This does not alter the thermal, mechanical and chemical resistance of the plate in any way.

2. Assembly and installation

- Unpacking
 - Check the integrity of the unit after unpacking.
- The box includes
 - AREX 6 Digital/Digital PRO heating magnetic stirrer
 - Power supply cord
 - Instruction manual
- First installation
 - Place the unit on non-flammable surface
 - Make sure that the values on the rating plate, correspond to those of the power supply
 - Move the main switch (2) to the OFF position
 - Ensure that the socket provided with grounding is accordant to current safety norms and easy to reach. Use only the cable provided with the instrument.
 - Insert the mains power cable into the socket

NOTE: the main cable must remain far away from the hot plate. It can be substituted only by main cables with same features (T=90°C, connector C15).

3. Operating controls

Commissioning	➤ Switch on the instrument through the main switch (2) ➤ Display shows software version, the last set point values and OFF on both displays (if "Start mode" is set on Stop, see chapter 5)
Stirring	➤ Adjust motor speed set point by turning the speed control knob (3) ➤ Click the knob to start stirring ➤ The speed increases until set point achievement ➤ A microprocessor ensures constant speed even when the viscosity changes (counter-reaction) ➤ Switch off the stirring by clicking the knob
Heating	➤ Adjust heating plate temperature set point by turning the temperature control knob (4) ➤ During all the time of set point adjusting, the temperature set point led (C) is switched on ➤ Click the knob to start heating ➤ Display shows real plate temperature for 3 seconds and temperature set point for 1 second alternatively ➤ When temperature set point is visualized, the relative led (C) turns on ➤ While the instrument is heating, the Heating ON led (B) turns on ➤ Switch off the heating by clicking the knob ➤ When heating function is off and the heating plate temperature exceeds 50 °C, the message "Hot" is shown on the temperature display until temperature falls below 50 °C NOTE: this warning is not active if the instrument is not powered. NOTE: in case of black-out, once the power is back the device restarts in the set restarting mode (see chapter 5 "Strt mode")
Timer	➤ To access the timer function, click the Timer/autoreverse/lock button (F) ➤ On the displays t l N E H H.N N are visualized ➤ Set the timer turning the speed control knob until the desired value and click to confirm ➤ Timer led (D) turns on ➤ Timer countdown starts when one of the main functions is on or when temperature reaches set point value (see chapter 5, "Time Strt") ➤ When countdown is active, timer led flashes ➤ Esc timer menu: a) click twice timer/autoreverse button b) wait 5 seconds without touching any knob ➤ When timer countdown ends, E n d E t l N e is shown until operator action. It's possible stirring continues if the parameter "time Sped" is ON (see chapter 5)
Autoreverse	➤ To access autoreverse function click twice the Timer/autoreverse/lock button (F) ➤ On the displays A.rEu N N. s are visualized ➤ Set the autoreverse turning the speed control knob until the desired value and click to confirm ➤ The autoreverse led (E) turns on ➤ Autoreverse starts when stirring is on
Lock	➤ Holding the Timer/autoreverse/lock button (F) for 3 seconds, the instrument will lock it's settings during operations. ➤ Unlock the control panel by holding the Timer/autoreverse/lock button (F) for 3 seconds. ➤ If other buttons are clicked while the instrument is locked, the LEDs around the lock button will blink for few seconds.

4. External thermometers

VTF (only for AREX 6 Digital Pro)	➤ Move the main switch to off position ➤ Screw the threaded support rod into its seat on the back of the instrument ➤ Fasten the VTF thermoregulator onto the support rod ➤ Place the temperature probe in the flask ➤ Plug the VTF into the dedicated socket on the back of the instrument ➤ Switching on the instrument through the main switch, the probe led (A) turns on and display shows OFF ➤ Select operating temperature on VTF and after clicking temperature control knob, a dashed line (----) appears on the display NOTE: while VTF is inserted, it is not possible to set any temperature turning the knob on the instrument. The plate can heat until the maximum temperature set in the full scale (see chapter 5, "TEMP LIMIT"). It's possible to set the timer on the VTF only.
External probe	➤ Move the main switch to off position ➤ Screw the threaded support rod into its seat on the back of the instrument ➤ Fasten the clamp onto the support rod ➤ Place the temperature probe in the flask ➤ Plug the probe into the dedicated socket on the back of the instrument ➤ Switching on the instrument through the main switch, the probe led (A) turns on and display shows OFF ➤ Select operating temperature using the temperature control knob ➤ Click the knob to start heating ➤ The display shows the real temperature obtained by the probe. Each 3 seconds the temperature set point is shown on the display and the temperature set point led (C) turns on

5. Setting mode

Press both knobs for 3 seconds to enter into the setting mode when heating and stirring are off.

Clicking the Timer/autoreverse button (F) it is possible to pass from one of the following parameter to the next ones:

Parameter shown		Default value	Range	Description
Display 1	Display 2			
TENP	LINt	370	OFF - 50 - 370	Set full scale temperature: ➤ It limits the maximum value of temperature set point for the heating plate* (from 50 to 370°C step 10°C) ➤ If "OFF" is set, it is not possible to perform any action with the temperature control knob
SPEd	LINt	1700	100 - 1700	Set full scale speed: ➤ It limits the maximum value of speed set point for the motor (from 100 to 1700rpm, step 100rpm)
Strt	NOdE	StOP	StOP - run	Set starting mode: ➤ Stop: when the instrument is switched on through the main switch, OFF is visualized on both displays ➤ Run: when the instrument is switched on, it restarts to work with the last set points set
Ctrl	tyPE	FlNE	FlNE - FAST	Set control temperature with external probe: ➤ Fine: slow temperature set point reaching. Overshoot optimized ➤ Fast: fast temperature set point reaching. High overshoot
PrOb	AlAr	On	On - OFF	Set external probe safety: ➤ If "on" it gives AL6 alarm when a slow temperature increase is detected by the external probe. This function is active in case of: - Probe temperature < 50°C - T set point - T probe > 5°C ➤ If "on" it gives AL7 alarm when a fast temperature decrease is detected by the external probe.
PrOb	CAL	0.0	-10.0 ÷ 10.0	Alignment external probe: ➤ It allows external probe temperature alignment to a reference thermometer
tINE	Strt	SnAP	SnAP - SEtP	Set start timer countdown: ➤ Snap: timer countdown starts as soon as one main function (Temperature/Speed) is on ➤ SetP: timer countdown starts when temperature set point is reached
tINE	SPEd	OFF	On - OFF	Continue stirring after timer end: ➤ If "on", stirring continues after timer end
ANb.t	CAL	0.0	-10.0 ÷ 10.0	Heating plate probe calibration: ➤ It allows plate ambient temperature alignment to a reference thermometer.
COun	NOtO			Running time. Days (24H): ➤ Motor and plate working time are displayed in hours until 9999 hours (around 416 days). After, they are shown in days starting from 417 and a dot turns on as following, to differentiate days from hours . 4 1 7
COun	HEAt			
rSEt		nO	YES - nO	Reset: ➤ It allows to restore default setting parameters

Click the stirring control knob to visualize the set value. Turn the knob to set the new value and click to confirm.

To esc setting mode wait 10 seconds without touching any knob or push both knobs at the same time.

* The temperature control function of the heating plate can also be used as a safety thermostat. In this case the maximum temperature of the heating plate will not exceed the full scale temperature set on the AREX 6 Digital/Digital PRO meaning that a longer heating time is required in order to reach the temperature set with VTF or external probe inserted, but reducing temperature oscillation at set point value.

6. Error messages

When the display shows an error message, the stirring and heating functions stop automatically. To remove the error message, disconnect the instrument from the power supply.

Error code	Cause
AL1	Overtemperature (T > 430 °C)
AL2	Excessive heating time
AL3	The stirring system doesn't run correctly
AL4	Overtemperature of the external probe (T > 310 °C) – (with probe inserted)
AL5	Overtemperature of the safety probe (T > 430 °C)
AL6	Slow temperature increase read by the external probe
AL7	Fast temperature decrease read by the external probe

If an error message appears on the display, please contact VELP Scientifica's technical service department.

7. Maintenance and cleaning

Maintenance	➤ No routine or extraordinary maintenance is necessary ; ➤ Repairs must be carried out by authorized Velp personnel only ; ➤ Instrument must be transported in its original packaging any indications present on the original packaging must be followed (e.g. palletized); ➤ It is the responsibility of the user, to properly decontaminate the unit in case of hazardous substances remaining on the surface or interior of the device. If in doubt about the compatibility of a cleaning or decontamination product, contact the manufacturer or distributor.
Cleaning	➤ Disconnect the unit from the power supply and use a cloth dampened with an non-inflammable non-aggressive detergent.

8. Technical data

	Models	F20500460, F20500461	F20510460, F20510461
General features	Power supply	230 V / 50-60 Hz (+/-10%)	115V / 60 Hz (+/-10%)
	Dimensions (WxHxD)	160x105x280 mm (6.5x4.1x11 in)	160x105x280 mm (6.5x4.1x11 in)
	Weight	2.6 Kg (5.7 lb)	2.6 Kg (5.7 lb)
	Power input	630W, 2,7A	630W, 5,5A
	Construction material (body)	Aluminium body –Technopolymer enclosure	Aluminium body –Technopolymer enclosure
	Working in continuous	Admitted	Admitted
	Settable restart modality	Stop or work	Stop or work
	Maximum load on the plate	25kg	25kg
	Noisiness	<< 80 dBa	<< 80 dBa
	Environmental temperature admitted	+5...+40 °C	+5...+40 °C
	Storage temperature admitted	-10...+60 °C	-10...+60 °C
	Max humidity	80%	80%
	Level of electrical protection CEI EN60529	IP 42	IP 42
	Overvoltage category	II	II
	Pollution degree CEI EN61010-1	2	2
Heating plate	Max altitude	2000 m	2000 m
	Heating output	600 W	600 W
	Heating plate dimensions	Ø 135mm	Ø 135mm
	Programmable temperature range	0 - 370 °C	0 - 370 °C
	Type of temperature control	Digital	Digital
	Temperature selection	1 °C	1 °C
	Reading temperature resolution	1 °C	1 °C
	Hot plate alarm	Over 50°C	Over 50°C
	Overtemperature alarm	Over 430°C	Over 430°C
	Construction material (plate)	Aluminium	Aluminium
	Full scale temperature range	From 50°C to 370°C	From 50°C to 370°C
	Safety circuit	Separated with dedicated probe	Separated with dedicated probe
External Probe	Type	Pt100 Class A – Ø 3mm	Pt100 Class A – Ø 3mm
	Temperature control type	2 modalities: Fine, Fast	2 modalities: Fine, Fast
	Programmable temperature range	0 – 300°C	0 – 300°C
	Temperature selection	1°C	1°C
	Reading temperature resolution	1°C	1°C
	Accuracy	+/- 1°C*	+/- 1°C*
	Calibration	Da -10.0 a +10.0°C	Da -10.0 a +10.0°C
	Position alarm	Not inserted in the sample	Not inserted in the sample
	Overtemperature alarm	310°C	310°C
Stir	Electrical data	5 Vdc – 1W (max)	5 Vdc – 1W (max)
	Stirring capacity	20 l di H ₂ O	20 l di H ₂ O
	Programmable speed range	30 – 1700 rpm	30 – 1700 rpm
	Motor type	BLDC	BLDC
	Speed selection	5 rpm step	5 rpm step
	Autoreverse	From 5s to 99m:59s	From 5s to 99m:59s
	Autoreverse selection	1 s	1 s
	Full scale speed range	From 100 to 1700 rpm	From 100 to 1700 rpm
	Stirring alarm	Motor fault	Motor fault
Timer	Motor rating output	10 W	10 W
	Timer	From 1min to 99h:59min	From 1min to 99h:59min
	Timer selection	1 min	1 min
Counters	2 modalities	Immediate start or at temperature set point	Immediate start or at temperature set point
	Motor counter	Hours of work	Hours of work
Counters	Heating plate counter	Hours of work	Hours of work

VTF	Temperature accuracy	+/-0,5°C*	+/-0,5°C*
	Electrical data	12Vdc-1.2W	12Vdc-1.2W

* in following conditions: 800ml water in 1 liter glass beaker (diameter 105mm), stirring bar 8x40mm, 600rpm, 50°C.

9. Accessories / Spare parts

A00001056	Magnetic stir bar Ø 6x35mm	A00000325	MultiAluBlockTM, 4 pos. Ø28 x h 30mm
A00001057	Magnetic stir bar, Ø6x20 mm	A00000326	MultiAluBlockTM, 4 pos. Ø28 x h 24mm
A00001061	Magnetic stir bar, Ø10x60 mm	A00000327	MultiAluBlockTM, 4 pos. Ø21 x h 31 mm
A00001062	Magnetic stir bar Ø 3x6mm	A00000328	MultiAluBlockTM, 8 pos. Ø17 x h 26mm
A00001063	Magnetic stir bar Ø 4.5x12mm	A00000329	MultiAluBlockTM, 11 pos. Ø15 x h 20mm
A00000336	Magnetic cross shape stir bar, Ø10x5 mm	A00000337	MultiAluBlockTM, 11 pos. Ø12 x h 14 mm
A00000352	Magnetic cross shape stir bar, Ø20x8 mm	A00000338	MonoAluBlock, 17 pos. Ø28 x h 24 mm
A00000354	Magnetic disc stir bar, Ø10x6 mm	A00000339	MonoAluBlock, 17 pos. Ø28 x h 30 mm
A00000355	Magnetic disc stir bar, Ø20x10 mm	A00000340	MonoAluBlock, 17 pos. Ø28 x h 43 mm
A00000356	Magnetic stir bar, Ø8x40 mm	A00000341	MonoAluBlock, 40 pos. Ø12 x h 14 mm
A00000357	Magnetic octagonal stir bar, Ø10x13 mm	A00000351	Handle for AluBlock removal
A00000330	Hemispheric bowl for 50ml flasks	A00000342	PTFE Safety cover for bowl 50 ml *
A00000331	Hemispheric bowl for 100ml flasks	A00000343	PTFE Safety cover for bowl 100 ml *
A00000332	Hemispheric bowl for 250ml flasks	A00000344	PTFE Safety cover for bowl 250 ml *
A00000333	Hemispheric bowl for 500ml flasks	A00000345	PTFE Safety cover for bowl 500 ml *
A00000334	Hemispheric bowl for 1000ml flasks	A00000346	PTFE Safety cover for bowl 1000 ml *
A00000280	Clamp with probe support	A00001069	Support rod
A00000268	Temp. Probe AISI 316 Ti Ø3mm	A00000335	Protective cover ARE/AREX 6 series
A00000349	Temp. probe glass-coated Ø5mm		
A00000323	MultiAluBlockTM Base Ø135 mm	10000239	Bumpon 13Dx5H embedded
A00000324	MultiAluBlockTM, 4 pos. Ø28 x h 43mm	10007175	Knob 35D blue

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We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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EN 61326-1	Electrical equipment for measurement, control and laboratory use - EMC requirements
2015/863/EU (RoHS III)	Restriction of the use of certain hazardous substances
2012/19/EU (WEEE)	Waste of electric and electronic equipment
EN 61010-1	Safety requirements for electrical equipment for measurement, control, and laboratory use.
EN 61010-2-051	Particular requirements for laboratory equipment for mixing and stirring
EN 61010-2-010	Particular requirements for laboratory equipment for the heating of material

and satisfies the essential requirements of the following directives:

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11. Declaration of conformity ^{UK}CA

We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following regulations:

S.I. 2016/1101 Electrical Equipment (Safety) Regulations 2016

S.I. 2016/1091 Electromagnetic Compatibility Regulations 2016

according to the relevant designated standards:

EN 61010-1 Safety requirements for electrical equipment for measurement, control, and laboratory use
Part 1: General requirements

EN 61010-2-010 Particular requirements for laboratory equipment for the heating of material

EN 61010-2-051 Particular requirements for laboratory equipment for mixing and stirring

EN 61326-1 Electrical equipment for measurement, control and laboratory use - EMC requirements -
Part 1: General requirements

and satisfies the essential requirements of regulations:

S.I. 2008/1597 Supply of Machinery (Safety) Regulations 2008

S.I. 2012/3032 Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic
Equipment Regulations 2012

S.I. 2013/3113 Waste Electrical and Electronic Equipment Regulations 2013

plus modifications.

Thank you for having chosen VELP!

Established in 1983, VELP is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. VELP works according to **ISO 9001**, **ISO14001** and **OHSAS 18001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

Analytical instruments

Elemental Analyzers
Digestion Units
Distillation Units
Solvent Extractors
Fiber Analyzers
Dietary Fiber Analyzers
Oxidation Stability Reactor
Consumables

Laboratory Equipment

Magnetic Stirrers
Heating Magnetic Stirrers
Heating Plates
Overhead stirrers
Vortex mixers
Homogenizers
COD Thermoreactors
BOD and Respirometers
Cooled Incubators
Flocculators
Overhead Shakers
Turbidimeter
Radiation Detector
Open Circulating Baths
Pumps

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