Instruction Manual
Manuale di istruzioni
Manuel d'instructions
Manual de instrucciones
Bedienungsanleitung
指导手册

ARE 6 Heating Magnetic Stirrer

F20500462, F20510462

AREX 6 Heating Magnetic Stirrer

F20500463, F20510463

**General Information / Informazioni Generali / Informations Générales / Información General /
Allgemeine Hinweise / 一般信息**

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This unit has been designed and manufactured in compliance with the following standards:

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Electrical equipment for laboratory use

UL 61010-1

General requirement - Canadian electrical code

CAN/CSA-C22.2 No.61010-1

VELP reserves the right to modify the characteristics of its products with the aim to constantly improving their quality.

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The plug disconnects the instrument. Therefore, place the instrument where it can be quickly disconnected.

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The heated solution may release toxic, dangerous or poisonous gases. Adequate safety measures must be taken, in accordance with the safety regulations in force, including the presence of hood and personal protective equipment (masks, gloves, goggles, etc.).

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Beware of the effect of the magnetic field on cardiac pacemakers and data media.

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Position the instrument on a flat surface, with a distance from the wall of 30 cm (at least).

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Do not use with explosive or dangerous materials for which the equipment is not designed. The stirrer must not be used in explosive atmospheres, in bain-marie or to stir combustible liquids that have a low combustion temperature. The product is intended for use with very small quantities of flammable liquids or flammable liquids that have a fire point higher than 625°C and a flash point higher than 600°C.

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It is responsibility of the user appropriately decontaminate the instrument in case of dangerous substances fall on or in it accordingly to the safety datasheet of substances used and to the current laboratories safety standards. It is not possible to decontaminate the product under steam.

It is also responsibility of the user to use substances for cleaning or decontaminating which do not react with internal parts of the instrument or with the material contained in it. In case of doubts on the compatibility of a cleaning solution, contact the manufacturer or local distributor.

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1. Introduction

EN

The ARE/AREX 6 heating magnetic stirrer is used in the laboratory for general use and for all those applications that require high regulation for stirring speed and heating liquids.

The Aluminum alloy heating plate ensures:

- Optimum heat distribution and a high specific power thanks to the circular configuration
- Temperature homogeneity
- High resistance to thermal stress and thermal shock

Magnetic stirring is generated by an AlNiCo magnet, driven by brushless motor which offers a virtually unlimited duration.

Size and shape of the magnetic stir bar determines the stirring efficiency at any given speed.

The stirring bar code which satisfies most applications is A00001056 (6x35 mm).

NOTE: The vessel must be made of a suitable material to withstand the foreseen temperature.

NOTE: Using the heating plate at high temperatures may cause discoloring. This does not alter the thermal, mechanical and chemical resistance of the plate in any way.

2. Assembly and installation

- Unpacking
 - Check the integrity of the unit after unpacking.
- The box includes
 - ARE/AREX 6 heating magnetic stirrer
 - Power supply cord
 - Instruction manual
- First installation
 - Place the unit on not-flammable surface
 - Make sure that the values on the rating plate, correspond to those of the power supply
 - Move the main switch to the OFF position
 - Ensure that the socket provided with grounding is accordant to current safety norms and easy to reach. Use only the cable provided with the instrument
 - Insert the mains power cable into the socket

NOTE: the mains cable must remain far away from the hot plate. It can be substituted only by main cables with same features (T=90°C, connector C15).

3. Operating controls

Commissioning	➤ Switch on the instrument through the main switch ➤ The right led between the two potentiometers turns on
Stirring	➤ Adjust motor speed set point by turning the speed control potentiometer ➤ While the instrument is stirring, the led over the speed control potentiometer turns on ➤ A microprocessor ensures constant speed even when the viscosity changes (counter-reaction) ➤ Switch off the stirring by turning the speed control potentiometer to value zero
Heating	➤ Adjust heating plate temperature set point by turning the temperature control potentiometer ➤ Start heating clicking the central button. The button lights on. ➤ While the instrument is heating, the led over the temperature control potentiometer turns on. ➤ Switch off the heating clicking again the central button. ➤ When the heating plate temperature exceeds 50 °C, the left led between the two potentiometers turns on until temperature falls below 50 °C NOTE: this warning is not active if the instrument is not powered.

4. External thermometers

VTF (only for AREX 6)	➤ Move the main switch to off position ➤ Screw the threaded support rod into its seat on the back of the instrument ➤ Fasten the VTF thermoregulator onto the support rod ➤ Place the temperature probe in the flask ➤ Plug the VTF into the dedicated socket on the back of the instrument ➤ Switch on the instrument through the main switch and set the temperature set point on the VTF ➤ Move the temperature potentiometer to the max value* and click the central button
------------------------------	---

* The temperature control function of the heating plate can also be used as a safety thermostat. In this case the maximum temperature of the heating plate will not exceed the temperature setting on the AREX 6 meaning that a longer heating time is required in order to reach the VTF thermoregulator temperature set reducing temperature oscillation at set point value.

5. Setting Mode

The unit can be set in two different starting mode:

Default Value	Range	Description
Stop	Stop - Run	Set starting mode: ➤ Stop: if the instrument is switched on through the main switch, the heating of the instrument doesn't work until the central button is pressed by the operator; ➤ Run: when the instrument is switched on, it restarts to work with the last set points set. The central button lights on and the instrument restarts to heat without any manual operation;

It's possible to switch from Stop-to-run mode following these steps:

- Switch off the instrument through the main switch.
- Keep pressed the central button and switch on the instrument through the main switch.
- When the button starts to light on intermittently move the finger.
- Clicking the central button only once, it switches off.
- If the instrument is switched off through the main switch with the central button switched off, the stop mode is saved.
- If the central button is clicked twice, will cause it to light continuously.
- Switching off the instrument through the main switch with the central button lighted on, the run mode is saved.

If the instrument is switched off when the central button lights on intermittently, the last starting mode remains set.

6. Error messages

The unit is fitted with safety devices which cut-off the power supply to the heating plate in the case of malfunctions.

The leds on the frontal panel indicate the type of malfunction:

Display indicator	Malfunction:
Flashing Heating led and Hot plate led ON	Overtemperature of the working probe (T> 430°C)
Flashing Heating led and Hot plate led OFF	Overtemperature of the safety probe (T>430°C)
Flashing Stirring led	Stirring system doesn't run correctly

Should any of the above occur, please contact your nearest VELP Scientifica service centre.

7. Maintenance and cleaning

Maintenance	➤ No routine or extraordinary maintenance is necessary ; ➤ Repairs must be carried out by authorized Velp personnel only ; ➤ Instrument must be transported in its original packaging any indications present on the original packaging must be followed (e.g. palletized); ➤ It is the responsibility of the user, to properly decontaminate the unit in case of hazardous substances remaining on the surface or interior of the device. If in doubt about the compatibility of a cleaning or decontamination product, contact the manufacturer or distributor.
Cleaning	➤ Disconnect the unit from the power supply and use a cloth dampened with an non-inflammable non-aggressive detergent.

8. Technical data

	Models	F20500462, F20500463	F20510462, F20510463
General features	Power supply	230 V / 50-60 Hz (+/-10%)	115V / 60 Hz (+/-10%)
	Dimensions WxHxD	160x105x280 mm (6.5x4.1x11 in)	160x105x280 mm (6.5x4.1x11 in)
	Weight	2.6 Kg (5.7 lb)	2.6 Kg (5.7 lb)
	Power input	630 W	630 W
	Construction material (body)	Aluminium body – Technopolymer enclosure	Aluminium body –Technopolymer enclosure
	Working in continuous	Admitted	Admitted
	Maximum load on the plate	25kg	25kg
	Noisiness	<< 80 dBA	<< 80 dBA
	Environmental temperature admitted	+5...+40 °C	+5...+40 °C
	Storage temperature admitted	-10...+60 °C	-10...+60 °C
	Max humidity	80%	80%
	Level of electrical protection CEI EN60529	IP 42	IP 42
	Overvoltage category	II	II
	Pollution degree CEI EN61010-1	2	2
	Max altitude	2000 m	2000 m
Heating plate	Heat output	600 W	600 W
	Heating plate dimensions	Ø 135mm	Ø 135mm
	Programmable temperature range	50 - 370 °C	50 - 370 °C
	Type of temperature control	Analog	Analog

	Construction material (plate)	Aluminium	Aluminium
	Safety circuit	Separated with dedicated probe	Separated with dedicated probe
	Hot plate alarm	Over 50°C	Over 50°C
	Overtemperature alarm	Over 430°C	Over 430°C
Stir	Programmable speed range	30 – 1700 rpm	30 – 1700 rpm
	Stirring capacity	20 l di H ₂ O	20 l di H ₂ O
	Motor type	BLDC	BLDC
	Motor rating output	10 W	10 W
VTF	Temperature accuracy	+/-0,5°C*	+/-0,5°C*
	Electrical data	12Vdc – 1.2W (max)	12Vdc – 1.2W (max)

* in following conditions: 800ml water in 1 liter glass beaker (diameter 105mm), stirring bar 6x35mm, 600rpm, 50°C.

9. Accessories / Spare parts

F208B0063	Electronic Vertex VTF -10...+300 °C	A00000326	MultiAluBlockTM, 4 pos. Ø28 x h 24mm
40000781	Spiral electric cable complete of plug *	A00000327	MultiAluBlockTM, 4 pos. Ø21 x h 31 mm
A00001056	Magnetic stir bar, Ø 6x35mm	A00000328	MultiAluBlockTM, 8 pos. Ø17 x h 26mm
A00001057	Magnetic stir bar, Ø6x20 mm	A00000329	MultiAluBlockTM, 11 pos. Ø15 x h 20mm
A00001061	Magnetic stir bar, Ø10x60 mm	A00000337	MultiAluBlockTM, 11 pos. Ø12 x h 14 mm
A00001062	Magnetic stir bar, Ø 3x6mm	A00000338	MonoAluBlock, 17 pos. Ø28 x h 24 mm
A00001063	Magnetic stir bar, Ø 4.5x12mm	A00000339	MonoAluBlock, 17 pos. Ø28 x h 30 mm
A00000336	Magnetic cross shape stir bar, Ø10x5 mm	A00000340	MonoAluBlock, 17 pos. Ø28 x h 43 mm
A00000352	Magnetic cross shape stir bar, Ø20x8 mm	A00000341	MonoAluBlock, 40 pos. Ø12 x h 14 mm
A00000354	Magnetic disc stir bar, Ø10x6 mm	A00000351	Handle for AluBlock removal
A00000355	Magnetic disc stir bar, Ø20x10 mm	A00000342	PTFE Safety cover for bowl 50 ml **
A00000356	Magnetic stir bar, Ø8x40 mm	A00000343	PTFE Safety cover for bowl 100 ml **
A00000357	Magnetic octagonal stir bar, Ø10x13 mm	A00000344	PTFE Safety cover for bowl 250 ml **
A00000330	Hemispheric bowl for 50ml flasks	A00000345	PTFE Safety cover for bowl 500 ml **
A00000331	Hemispheric bowl for 100ml flasks	A00000346	PTFE Safety cover for bowl 1000 ml **
A00000332	Hemispheric bowl for 250ml flasks	A00001069	Support rod
A00000333	Hemispheric bowl for 500ml flasks	A00000335	Protective cover ARE/AREX 6 series
A00000334	Hemispheric bowl for 1000ml flasks		
A00000323	MultiAluBlockTM Base Ø135 mm	10000239	Bumpon 13Dx5H embedded
A00000324	MultiAluBlockTM, 4 pos. Ø28 x h 43mm	10007258	Knob 35D blue with notch
A00000325	MultiAluBlockTM, 4 pos. Ø28 x h 30mm		

* For connection to other contact thermometers other than VTF

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We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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EN 61326-1	Electrical equipment for measurement, control and laboratory use - EMC requirements -
2015/863/EU (RoHS III)	Restriction of the use of certain hazardous substances
2012/19/EU (WEEE)	Waste of electric and electronic equipment
EN 61010-1	Safety requirements for electrical equipment for measurement, control, and laboratory use
EN 61010-2-051	Particular requirements for laboratory equipment for mixing and stirring
EN 61010-2-010	Particular requirements for laboratory equipment for the heating of material

and satisfies the essential requirements of the following directives:

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11. Declaration of conformity ^{UK}CA

We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following regulations:

S.I. 2016/1101 Electrical Equipment (Safety) Regulations 2016

S.I. 2016/1091 Electromagnetic Compatibility Regulations 2016

according to the relevant designated standards:

EN 61010-1 Safety requirements for electrical equipment for measurement, control, and laboratory use
Part 1: General requirements

EN 61010-2-010 Particular requirements for laboratory equipment for the heating of material

EN 61010-2-051 Particular requirements for laboratory equipment for mixing and stirring

EN 61326-1 Electrical equipment for measurement, control and laboratory use - EMC requirements -
Part 1: General requirements

and satisfies the essential requirements of regulations:

S.I. 2008/1597 Supply of Machinery (Safety) Regulations 2008

S.I. 2012/3032 Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic
Equipment Regulations 2012

S.I. 2013/3113 Waste Electrical and Electronic Equipment Regulations 2013

plus modifications.

Thank you for having chosen VELP!

Established in 1983, VELP is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. VELP works according to **ISO 9001**, **ISO14001** and **OHSAS 18001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

Analytical instruments

Elemental Analyzers
Digestion Units
Distillation Units
Solvent Extractors
Fiber Analyzers
Dietary Fiber Analyzers
Oxidation Stability Reactor
Consumables

Laboratory Equipment

Magnetic Stirrers
Heating Magnetic Stirrers
Heating Plates
Overhead stirrers
Vortex mixers
Homogenizers
COD Thermoreactors
BOD and Respirometers
Cooled Incubators
Flocculators
Overhead Shakers
Turbidimeter
Radiation Detector
Open Circulating Baths
Pumps

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VELP Scientifica Srl