**Instruction Manual**  
**Manuale di istruzioni**  
**Manuel d'instructions**  
**Manual de instrucciones**  
**Bedienungsanleitung**

## **AREX 6 CONNECT PRO Heating Magnetic Stirrer**

F20500465, F20510465

## **AREX 6 CONNECT Heating Magnetic Stirrer**

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**This unit must be used for laboratory applications indoor only.** The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions. If the product is used in a not specified way by the manufacturer or with non specified accessories, product's safety may be compromised.

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**This unit has been designed and manufactured in compliance with the following standards:**

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Safety requirements for electrical equipment for measurement, control and for laboratory use

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Electrical equipment for laboratory use

**IEC/EN 61010-1**

**IEC/EN 61010-2-051**

**IEC/EN 61010-2-010**

**UL 61010-1**

General requirement - Canadian electrical code

**CAN/CSA-C22.2 No.61010-1**

VELP reserves the right to modify the characteristics of its products with the aim to constantly improving their quality.

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The heated solution may release toxic, dangerous or poisonous gases. Adequate safety measures must be taken, in accordance with the safety regulations in force, including the presence of hood and personal protective equipment (masks, gloves, goggles, etc.).

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Do not use with explosive or dangerous materials for which the equipment is not designed. The stirrer must not be used in explosive atmospheres, in bain-marie or to stir combustible liquids that have a low combustion temperature. The product is intended for use with very small quantities of flammable liquids or flammable liquids that have a fire point higher than 625°C and a flash point higher than 600°C.

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It is responsibility of the user appropriately decontaminate the instrument in case of dangerous substances fall on or in it accordingly to the safety datasheet of substances used and to the current laboratories safety standards. It is not possible to decontaminate the product under steam.

It is also responsibility of the user to use substances for cleaning or decontaminating which do not react with internal parts of the instrument or with the material contained in it. In case of doubts on the compatibility of a cleaning solution, contact the manufacturer or local distributor.

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
## 1. Introduction

EN

The AREX 6 Connect PRO heating magnetic stirrer is used in the laboratory for general use and for all those applications that require a precise regulation for stirring speed and heating liquids.

The integrated Wi-Fi module transmits and receives information in order to allow the equipment monitoring and control through VELP Ermes platform.

- |   |                                           |   |                               |
|---|-------------------------------------------|---|-------------------------------|
| 1 | Aluminum alloy plate with ceramic coating | A | Probe led                     |
| 2 | Main switch                               | B | Heating ON led                |
| 3 | Speed control knob                        | C | Temperature set point led     |
| 4 | Temperature control knob                  | D | Timer led                     |
| 5 | Display                                   | E | Wi-Fi led                     |
|   |                                           | F | Autoreverse led               |
|   |                                           | G | Timer/autoreverse/lock button |

The Aluminum alloy heating plate with ceramic coating (1)

- Optimum heat distribution and a high specific power thanks to the circular configuration
- Temperature homogeneity
- High resistance to thermal stress and thermal shock
- High resistance to chemical agents and surface abrasion

ensures:

Magnetic stirring is generated by an AlNiCo magnet, driven by brushless motor which offers a virtually unlimited duration. Size and shape of the magnetic stir bar determines the stirring efficiency at any given speed. The stirring bar code which satisfies most applications is Ø8x40mm.

**NOTE:** The vessel must be made of a suitable material to withstand the foreseen temperature.

**NOTE:** Using the heating plate at high temperatures may cause discoloring. This does not alter the thermal, mechanical and chemical resistance of the plate in any way.

## 2. Assembly and installation

- Unpacking
  - Check the integrity of the unit after unpacking.
- The box includes
  - AREX 6 Connect/Connect PRO heating magnetic stirrer
  - Power supply cord
  - Instruction manual
- First installation
  - Place the unit on not-flammable surface
  - Make sure that the values on the rating plate, correspond to those of the power supply
  - Move the main switch (2) to the OFF position
  - Ensure that the socket provided with grounding is accordant to current safety norms and easy to reach. Use only the cable provided with the instrument.
  - Insert the mains power cable into the socket

**NOTE:** the mains cable must remain far away from the hot plate. It can be substituted only by main cables with same features (T=90°C, connector C15).

### 3. VELP Ermes Configuration

VELP Ermes is a revolutionary cloud platform that transforms and improves your laboratory experience by creating an ecosystem of instruments, people and data. The VELP Ermes platform is able to reduce distances and accelerate scientific processes in total safety. In order to access on ERMES, you need to enable your VELP account by selecting "Configure your VELP ERMES account" at <http://www.velp.com/en/login>.

To be able to communicate, the instrument needs to be in the operating range of laboratory Wi-Fi (2.4 GHz) and be configured as follows:

- Switch on the AREX 6 Connect PRO and using the PC/Tablet/Mobile phone, select the AREX6\_SERIAL NUMBER available on the Wi-Fi list, in order to connect directly to the instrument.
- Open a browser on the PC/Tablet/Mobile phone and insert the address 192.168.0.50 to reach the configuration page.
- Set the parameters required to connect to your Wi-Fi (network name, password, security, mac address, etc.) according to your internal procedure and save. As Passkey, insert "velp". If necessary, contact your IT administrator.
- Restart the AREX 6 Connect PRO (turn off and on) and wait until the Wi-Fi LED is ON constantly. If the LED is still flashing after 2 minutes, try the procedure again selecting "STAR AP YES" from the parameter menu (see chapter 6).
- Select OTC - - - from the parameter menu (see chapter 6) and proceed with the product registration from the VELP Ermes platform. For more information see FAQ on VELP website.

**NOTE:** To access to VELP Ermes is necessary to have a VELP account.

### 4. Operating controls

|                      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Commissioning</b> | <ul style="list-style-type: none"> <li>➤ Switch on the instrument through the main switch (2)</li> <li>➤ Display shows software version, the last set point values and OFF on both displays (if "Start mode" is set on Stop, see chapter 5)</li> </ul>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| <b>Stirring</b>      | <ul style="list-style-type: none"> <li>➤ Adjust motor speed set point by turning the speed control knob (3)</li> <li>➤ Click the knob to start stirring</li> <li>➤ The speed increases until set point achievement</li> <li>➤ A microprocessor ensures constant speed even when the viscosity changes (counter-reaction)</li> <li>➤ Switch off the stirring by clicking the knob</li> </ul>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| <b>Heating</b>       | <ul style="list-style-type: none"> <li>➤ Adjust heating plate temperature set point by turning the temperature control knob (4)</li> <li>➤ During all the time of set point adjusting, the temperature set point led (C) is switched on</li> <li>➤ Click the knob to start heating</li> <li>➤ Display shows real plate temperature for 3 seconds and temperature set point for 1 second alternatively</li> <li>➤ When temperature set point is visualized, the relative led (C) turns on</li> <li>➤ While the instrument is heating, the Heating ON led (B) turns on</li> <li>➤ Switch off the heating by clicking the knob</li> <li>➤ When heating function is off and the heating plate temperature exceeds 50 °C, the message "Hot" is shown on the temperature display until temperature falls below 50 °C</li> </ul> <p><b>NOTE:</b> this warning is not active if the instrument is not powered.</p> <p><b>NOTE:</b> in case of black-out, once the power is back the device restarts in the set start mode (see chapter 5 "Start mode")</p> |
| <b>Timer</b>         | <ul style="list-style-type: none"> <li>➤ To access the timer function, push the Timer/autoreverse/lock button (G)</li> <li>➤ On the displays  are visualized</li> <li>➤ Set the timer turning the speed control knob until the desired value and click to confirm</li> <li>➤ Timer led (D) turns on</li> <li>➤ Timer countdown starts when one of the main functions is on or when temperature reaches set point value (see chapter 5, "Time Strt")</li> <li>➤ When countdown is active, led timer flashes</li> <li>➤ Esc timer menu: a) click twice timer/autoreverse button (G)<br/>b) wait 5 seconds without touching any knob</li> <li>➤ When timer countdown ends,  is shown until operator action. It's possible stirring continues if the parameter "time Sped" is ON (see chapter 5)</li> </ul>                                                                                                                                                                                                                                            |
| <b>Autoreverse</b>   | <ul style="list-style-type: none"> <li>➤ To access autoreverse function push twice the Timer/autoreverse/lock button (G)</li> <li>➤ On the displays  are visualized</li> <li>➤ Set the auto reverse turning the speed control knob until the desired value and click to confirm</li> <li>➤ The autoreverse led (F) turns on</li> <li>➤ Autoreverse starts when stirring is on</li> </ul>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
| <b>Lock</b>          | <ul style="list-style-type: none"> <li>➤ Holding the Timer/autoreverse/lock button (G) for 3 seconds, the instrument will lock its settings during operations.</li> <li>➤ Unlock the control panel by holding the Timer/autoreverse/lock button (G) for 3 seconds.</li> <li>➤ If other buttons are clicked while the instrument is locked, the two LEDs around the lock button will blink for few seconds.</li> </ul>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |

## 5. External thermometers

|                                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
|------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>VTF (only for AREX 6 Connect Pro)</b> | <ul style="list-style-type: none"> <li>➤ Move the main switch to off position</li> <li>➤ Screw the threaded support rod into its seat on the back of the instrument</li> <li>➤ Fasten the VTF thermoregulator onto the support rod</li> <li>➤ Place the temperature probe in the flask</li> <li>➤ Plug the VTF into the dedicated socket on the back of the instrument</li> <li>➤ Switching on the instrument through the main switch, the probe led (A) turns on and display shows OFF</li> <li>➤ Select operating temperature on VTF and after clicking temperature control knob, a dashed line (----) appears on the display</li> </ul> <p><b>NOTE:</b> while VTF is inserted, it is not possible to set any temperature turning the knob on the instrument. The plate can heat until the maximum temperature set in the full scale (see chapter 5, "Temp Limt"). It is possible to set the timer only on the VTF.</p> |
| <b>External probe</b>                    | <ul style="list-style-type: none"> <li>➤ Move the main switch to off position</li> <li>➤ Screw the threaded support rod into its seat on the back of the instrument</li> <li>➤ Fasten the clamp onto the support rod</li> <li>➤ Place the temperature probe in the flask</li> <li>➤ Plug the probe into the dedicated socket on the back of the instrument</li> <li>➤ Switching on the instrument through the main switch, the probe led (A) turns on and display shows OFF</li> <li>➤ Select operating temperature using the temperature control knob</li> <li>➤ Click the knob to start heating</li> <li>➤ The display shows the real temperature obtained by the probe. Each 3 seconds the temperature set point is shown on the display and the temperature set point led (C) turns on</li> </ul>                                                                                                                     |

## 6. Setting mode

Press both knobs for 3 seconds to enter into the setting mode when heating and stirring are off.

Clicking the Timer/autoreverse button (G) it is possible to pass from one of the following parameter to the next ones:

| Parameter shown |           | Default value | Range          | Description                                                                                                                                                                                                                                                                                                                                                                                                                                            |
|-----------------|-----------|---------------|----------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Display 1       | Display 2 |               |                |                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| TENP            | LINt      | 370           | OFF - 50 - 370 | Set full scale temperature: <ul style="list-style-type: none"> <li>➤ It limits the maximum value of temperature set point for the heating plate* (from 50 to 370°C step 10°C)</li> <li>➤ If "OFF" is set, it is not possible to perform any action with the temperature control knob</li> </ul>                                                                                                                                                        |
| SPEd            | LINt      | 1700          | 100 - 1700     | Set full scale speed: <ul style="list-style-type: none"> <li>➤ It limits the maximum value of speed set point for the motor (from 100 to 1700rpm, step 100rpm)</li> </ul>                                                                                                                                                                                                                                                                              |
| Strt            | NODE      | StOP          | StOP - run     | Set starting mode: <ul style="list-style-type: none"> <li>➤ Stop: when the instrument is switched on through the main switch, OFF is visualized on both displays</li> <li>➤ Run: when the instrument is switched on, it restarts to work with the last set points set</li> </ul>                                                                                                                                                                       |
| Ctrl            | tyPE      | FlNE          | FlNE - FAST    | Set control temperature with external probe: <ul style="list-style-type: none"> <li>➤ Fine: slow temperature set point reaching. Overshoot optimized</li> <li>➤ Fast: fast temperature set point reaching. High overshoot</li> </ul>                                                                                                                                                                                                                   |
| PrOb            | AlAr      | On            | On - OFF       | Set external probe safety: <ul style="list-style-type: none"> <li>➤ If "on" it gives AL6 alarm when a slow temperature increase is detected by the external probe. This function is active in case of: <ul style="list-style-type: none"> <li>- Probe temperature &lt; 50°C</li> <li>- T set point - T probe &gt; 5°C</li> </ul> </li> <li>➤ If "on" it gives AL7 alarm when a fast temperature decrease is detected by the external probe.</li> </ul> |
| PrOb            | CAL       | 0.0           | -10.0 ÷ 10.0   | Alignment external probe: <ul style="list-style-type: none"> <li>➤ It allows external probe temperature alignment to a reference thermometer</li> </ul>                                                                                                                                                                                                                                                                                                |
| tINE            | Strt      | SnAP          | SnAP - SEtP    | Set start timer countdown: <ul style="list-style-type: none"> <li>➤ Snap: timer countdown starts as soon as one main function (Temperature/Speed) is on</li> <li>➤ SetP: timer countdown starts when temperature set point is reached</li> </ul>                                                                                                                                                                                                       |
| TINE            | SPEd      | On            | On - OFF       | Continue stirring after timer end: <ul style="list-style-type: none"> <li>➤ If "on", stirring continues after timer end</li> </ul>                                                                                                                                                                                                                                                                                                                     |

|       |      |      |              |                                                                                                                                                                                                                                                                                                                                         |
|-------|------|------|--------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| ANb.t | CAL  | 0.0  | -10.0 ÷ 10.0 | Heating plate probe calibration:<br>➤ It allows plate ambient temperature alignment to a reference thermometer.                                                                                                                                                                                                                         |
| COun  | NOtO |      |              | Running time. Days (24H):<br>➤ Motor and plate working time are displayed in hours until 9999 hours (around 416 days).<br>After, they are shown in days starting from 417 and a dot turns on as following, to differentiate days from hours<br><div style="border: 1px solid black; padding: 2px; display: inline-block;">. 4 1 7</div> |
| COun  | HEAt |      |              |                                                                                                                                                                                                                                                                                                                                         |
| rSEt  |      | nO   | YES - nO     | Reset:<br>➤ It allows to restore default setting parameters.                                                                                                                                                                                                                                                                            |
| OtC   | cFG  | ---- |              | Select this when you need to register the instrument using ERMES, in order to receive the unique security code.                                                                                                                                                                                                                         |
| StAr  | AP   | nO   | YES-nO       | Select this if you need to change the connection parameters to your WiFi network (chapter 3).                                                                                                                                                                                                                                           |

Click the stirring control knob to visualize the set value. Turn the knob to set the new value and click to confirm.  
To esc setting mode wait 10 seconds without touching any knob or push both knobs at the same time.

\* The temperature control function of the heating plate can also be used as a safety thermostat. In this case the maximum temperature of the heating plate will not exceed the full scale temperature set on the AREX 6 Digital/Digital PRO meaning that a longer heating time is required in order to reach the temperature set with VTF or external probe inserted, but reducing temperature oscillation at set point value.

## 7. Error messages

When the display shows an error message, the stirring and heating functions stop automatically.  
To remove the error message, disconnect the instrument from the power supply.

| Error code | Cause                                                                      |
|------------|----------------------------------------------------------------------------|
| AL1        | Overtemperature (T > 430 °C)                                               |
| AL2        | Excessive heating time                                                     |
| AL3        | The stirring system doesn't run correctly                                  |
| AL4        | Overtemperature of the external probe (T > 310 °C) – (with probe inserted) |
| AL5        | Overtemperature of the safety probe (T > 430 °C)                           |
| AL6        | Slow temperature increase read by the external probe                       |
| AL7        | Fast temperature decrease read by the external probe                       |

If an error message appears on the display, please contact VELP Scientifica's technical service department.

## 8. Maintenance and cleaning

|                    |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
|--------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Maintenance</b> | <ul style="list-style-type: none"> <li>➤ No routine or extraordinary maintenance is necessary;</li> <li>➤ Repairs must be carried out by authorized Velp personnel only;</li> <li>➤ Instrument must be transported in its original packaging any indications present on the original packaging must be followed (e.g. palletized);</li> <li>➤ It is the responsibility of the user, to properly decontaminate the unit in case of hazardous substances remaining on the surface or interior of the device. If in doubt about the compatibility of a cleaning or decontamination product, contact the manufacturer or distributor.</li> </ul> |
| <b>Cleaning</b>    | <ul style="list-style-type: none"> <li>➤ Disconnect the unit from the power supply and use a cloth dampened with an non-inflammable non-aggressive detergent.</li> </ul>                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |

## 9. Technical data

|                         |                                            |                                          |                                          |
|-------------------------|--------------------------------------------|------------------------------------------|------------------------------------------|
| <b>General features</b> | Models                                     | F20500465, F20500466                     | F20510465, F20510466                     |
|                         | Power supply                               | 230 V / 50-60 Hz (+/-10%)                | 115V / 60 Hz (+/-10%)                    |
|                         | Dimensions (WxHxD)                         | 160x105x280 mm (6.5x4.1x11 in)           | 160x105x280 mm (6.5x4.1x11 in)           |
|                         | Weight                                     | 2.6 Kg (5.7 lb)                          | 2.6 Kg (5.7 lb)                          |
|                         | Power input                                | 630W, 2,7A                               | 630W, 5,5A                               |
|                         | Construction material (body)               | Aluminium body – Technopolymer enclosure | Aluminium body – Technopolymer enclosure |
|                         | Working in continuous                      | Admitted                                 | Admitted                                 |
|                         | Settable restart modality                  | Stop or work                             | Stop or work                             |
|                         | Maximum load on the plate                  | 25kg                                     | 25kg                                     |
|                         | Noisiness                                  | << 80 dBa                                | << 80 dBa                                |
|                         | Environmental temperature admitted         | +5...+40 °C                              | +5...+40 °C                              |
|                         | Storage temperature admitted               | -10...+60 °C                             | -10...+60 °C                             |
|                         | Max humidity                               | 80%                                      | 80%                                      |
|                         | Level of electrical protection CEI EN60529 | IP 42                                    | IP 42                                    |

|                       |                                |                                             |                                             |
|-----------------------|--------------------------------|---------------------------------------------|---------------------------------------------|
|                       | Overvoltage category           | II                                          | II                                          |
|                       | Pollution degree CEI EN61010-1 | 2                                           | 2                                           |
|                       | Max altitude                   | 2000 m                                      | 2000 m                                      |
| <b>Heating plate</b>  | Heat output                    | 600 W                                       | 600 W                                       |
|                       | Heating plate dimensions       | Ø 135mm                                     | Ø 135mm                                     |
|                       | Programmable temperature range | 0 - 370 °C                                  | 0 - 370 °C                                  |
|                       | Type of temperature control    | Digital                                     | Digital                                     |
|                       | Temperature selection          | 1 °C                                        | 1 °C                                        |
|                       | Reading temperature resolution | 1 °C                                        | 1 °C                                        |
|                       | Hot plate alarm                | Over 50°C                                   | Over 50°C                                   |
|                       | Overtemperature alarm          | Over 430°C                                  | Over 430°C                                  |
|                       | Construction material (plate)  | Aluminium                                   | Aluminium                                   |
|                       | Full scale temperature range   | From 50°C to 370°C                          | From 50°C to 370°C                          |
|                       | Safety circuit                 | Separated with dedicated probe              | Separated with dedicated probe              |
|                       |                                |                                             |                                             |
| <b>External Probe</b> | Type                           | Pt100 Class A – Ø 3mm                       | Pt100 Class A – Ø 3mm                       |
|                       | Temperature control type       | 2 modalities: Fine, Fast                    | 2 modalities: Fine, Fast                    |
|                       | Programmable temperature range | 0 – 300°C                                   | 0 – 300°C                                   |
|                       | Temperature selection          | 1°C                                         | 1°C                                         |
|                       | Reading temperature resolution | 1°C                                         | 1°C                                         |
|                       | Accuracy                       | +/- 1°C*                                    | +/- 1°C*                                    |
|                       | Calibration                    | Da -10.0 a +10.0°C                          | Da -10.0 a +10.0°C                          |
|                       | Position alarm                 | Not inserted in the sample                  | Not inserted in the sample                  |
|                       | Overtemperature alarm          | 310°C                                       | 310°C                                       |
|                       | Electrical data                | 5 Vdc – 1W (max)                            | 5 Vdc – 1W (max)                            |
| <b>Stir</b>           | Stirring capacity              | 20 l di H <sub>2</sub> O                    | 20 l di H <sub>2</sub> O                    |
|                       | Programmable speed range       | 30 – 1700 rpm                               | 30 – 1700 rpm                               |
|                       | Motor type                     | BLDC                                        | BLDC                                        |
|                       | Speed selection                | 5 rpm step                                  | 5 rpm step                                  |
|                       | Autoreverse                    | From 5s to 99m:59s                          | From 5s to 99m:59s                          |
|                       | Autoreverse selection          | 1 s                                         | 1 s                                         |
|                       | Full scale speed range         | From 100 to 1700 rpm                        | From 100 to 1700 rpm                        |
|                       | Stirring alarm                 | Motor fault                                 | Motor fault                                 |
|                       | Motor rating output            | 10 W                                        | 10 W                                        |
| <b>Timer</b>          | Timer                          | From 1min to 99h:59min                      | From 1min to 99h:59min                      |
|                       | Timer selection                | 1 min                                       | 1 min                                       |
|                       | 2 modalities                   | Immediate start or at temperature set point | Immediate start or at temperature set point |
| <b>Counters</b>       | Motor counter                  | Hours of work                               | Hours of work                               |
|                       | Heating plate counter          | Hours of work                               | Hours of work                               |
| <b>VTF</b>            | Temperature accuracy           | +/-0,5°C*                                   | +/-0,5°C*                                   |
|                       | Electrical data                | 12Vdc-1.2W                                  | 12Vdc-1.2W                                  |
| <b>Wi-Fi</b>          | Frequency                      | 2.4GHz                                      | 2.4GHz                                      |
|                       | Protocols                      | 802.11 b/g/n                                | 802.11 b/g/n                                |
|                       | Output Power                   | 15.42 dBm EIRP                              | 15.42 dBm EIRP                              |

\* in the following conditions: 800ml water in 1 liter glass beaker (diameter 105mm), stirring bar 8x40mm, 600rpm, 50°C.

## 10. Accessories / Spare parts

|           |                                         |           |                                        |
|-----------|-----------------------------------------|-----------|----------------------------------------|
| A00001056 | Magnetic stir bar Ø 6x35mm              | A00000325 | MultiAluBlockTM, 4 pos. Ø28 x h 30mm   |
| A00001057 | Magnetic stir bar, Ø6x20 mm             | A00000326 | MultiAluBlockTM, 4 pos. Ø28 x h 24mm   |
| A00001061 | Magnetic stir bar, Ø10x60 mm            | A00000327 | MultiAluBlockTM, 4 pos. Ø21 x h 31 mm  |
| A00001062 | Magnetic stir bar Ø 3x6mm               | A00000328 | MultiAluBlockTM, 8 pos. Ø17 x h 26mm   |
| A00001063 | Magnetic stir bar Ø 4.5x12mm            | A00000329 | MultiAluBlockTM, 11 pos. Ø15 x h 20mm  |
| A00000336 | Magnetic cross shape stir bar, Ø10x5 mm | A00000337 | MultiAluBlockTM, 11 pos. Ø12 x h 14 mm |
| A00000352 | Magnetic cross shape stir bar, Ø20x8 mm | A00000338 | MonoAluBlock, 17 pos. Ø28 x h 24 mm    |
| A00000354 | Magnetic disc stir bar, Ø10x6 mm        | A00000339 | MonoAluBlock, 17 pos. Ø28 x h 30 mm    |
| A00000355 | Magnetic disc stir bar, Ø20x10 mm       | A00000340 | MonoAluBlock, 17 pos. Ø28 x h 43 mm    |
| A00000356 | Magnetic stir bar, Ø8x40 mm             | A00000341 | MonoAluBlock, 40 pos. Ø12 x h 14 mm    |
| A00000357 | Magnetic octagonal stir bar, Ø10x13 mm  | A00000351 | Handle for AluBlock removal            |
| A00000330 | Hemispheric bowl for 50ml flasks        | A00000342 | PTFE Safety cover for bowl 50 ml *     |
| A00000331 | Hemispheric bowl for 100ml flasks       | A00000343 | PTFE Safety cover for bowl 100 ml *    |
| A00000332 | Hemispheric bowl for 250ml flasks       | A00000344 | PTFE Safety cover for bowl 250 ml *    |
| A00000333 | Hemispheric bowl for 500ml flasks       | A00000345 | PTFE Safety cover for bowl 500 ml *    |
| A00000334 | Hemispheric bowl for 1000ml flasks      | A00000346 | PTFE Safety cover for bowl 1000 ml *   |
| A00000280 | Clamp with probe support                | A00001069 | Support rod                            |
| A00000268 | Temp. Probe AISI 316 Ti Ø3mm            | A00000335 | Protective cover ARE/AREX 6 series     |
| A00000349 | Temp. probe glass-coated Ø5mm           |           |                                        |
| A00000323 | MultiAluBlockTM Base Ø135 mm            | 10000239  | Bumpon 13Dx5H embedded                 |
| A00000324 | MultiAluBlockTM, 4 pos. Ø28 x h 43mm    | 10007175  | Knob 35D blue                          |

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We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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|                           |                                                                                                                                                        |
|---------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------|
| EN 61010-1                | Safety requirements for electrical equipment for measurement, control, and laboratory use<br>Part 1: General requirements                              |
| EN 61010-2-010            | Particular requirements for laboratory equipment for the heating of material                                                                           |
| EN 61010-2-051            | Particular requirements for laboratory equipment for mixing and stirring                                                                               |
| ETSI EN 301 489-1 V2.2.0  | Electromagnetic compatibility (EMC) standard for radio equipment and services;<br>Part 1: Common technical requirements                                |
| ETSI EN 301 489-17 V3.2.0 | Electromagnetic Compatibility (EMC) standard for radio equipment and services;<br>Part 17: Specific conditions for Broadband Data Transmission Systems |
| EN 62311                  | Assessment of electronic and electrical equipment related to human exposure restrictions for<br>electromagnetic fields (0 Hz to 300 GHz)               |
| ETSI EN 300 328 V1.9.1    | Wideband transmission systems; Data transmission equipment operating in the 2,4 GHz<br>band                                                            |

and satisfies the essential requirements of the following directives:

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|                        |                                                        |
|------------------------|--------------------------------------------------------|
| 2006/42/EC             | Machinery directive                                    |
| 2014/53/EU             | Radio Equipment directive (RED)                        |
| 2015/863/EU (RoHS III) | Restriction of the use of certain hazardous substances |
| 2012/19/EU (WEEE)      | Waste of electric and electronic equipment             |

plus modifications / più modifiche / plus modifications / más sucesivas modificaciones / in der jeweils gültigen Fassung.

## 12. Declaration of conformity <sup>UK</sup><sub>CA</sub>

We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following regulations:

S.I. 2017/1206                      Radio Equipment Regulations 2017

according to the relevant designated standards:

EN 61010-1                      Safety requirements for electrical equipment for measurement, control, and laboratory use  
Part 1: General requirements

EN 61010-2-010                  Particular requirements for laboratory equipment for the heating of material

EN 61010-2-051                  Particular requirements for laboratory equipment for mixing and stirring

ETSI EN 301 489-1 V2.2.0      Electromagnetic compatibility (EMC) standard for radio equipment and services;  
Part 1: Common technical requirements

ETSI EN 301 489-17 V3.2.0    ElectroMagnetic Compatibility (EMC) standard for radio equipment and services;  
Part 17: Specific conditions for Broadband Data Transmission Systems

EN 62311                         Assessment of electronic and electrical equipment related to human exposure restrictions for  
electromagnetic fields (0 Hz to 300 GHz)

ETSI EN 300 328 V1.9.1        Wideband transmission systems; Data transmission equipment operating in the 2,4 GHz  
band

and satisfies the essential requirements of regulations:

S.I. 2008/1597                      Supply of Machinery (Safety) Regulations 2008

S.I. 2012/3032                      Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic  
Equipment Regulations 2012

S.I. 2013/3113                      Waste Electrical and Electronic Equipment Regulations 2013

plus modifications.

## Thank you for having chosen VELP!

Established in 1983, VELP is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. VELP works according to **ISO 9001**, **ISO14001** and **OHSAS 18001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

### Analytical instruments

Elemental Analyzers  
Digestion Units  
Distillation Units  
Solvent Extractors  
Fiber Analyzers  
Dietary Fiber Analyzers  
Oxidation Stability Reactor  
Consumables

### Laboratory Equipment

Magnetic Stirrers  
Heating Magnetic Stirrers  
Heating Plates  
Overhead stirrers  
Vortex mixers  
Homogenizers  
COD Thermoreactors  
BOD and Respirometers  
Cooled Incubators  
Flocculators  
Overhead Shakers  
Turbidimeter  
Radiation Detector  
Open Circulating Baths  
Pumps

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VELP Scientifica Srl