Instruction Manual
Manuale di istruzioni
Manuel d'instructions
Manual de instrucciones
Bedienungsanleitung
操作指南

HSC Heating Magnetic Stirrer

F20500552, F20510552

General Information / Informazioni Generali / Informations Générales / Información General / Allgemeine Hinweise / 基本信息

Before using the unit, please read the following instruction manual carefully.

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Attention, surface chaude!

Prudencia, superficie caliente!

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Do not dispose of this equipment as urban waste, in accordance with EEC directive 2012/19/UE.

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The product can be used with flammable liquids.

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This unit must be used for laboratory applications indoor only.

The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions. If the product is used in a not specified way by the manufacturer or with non specified accessories, product's safety may be compromised.

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This unit has been designed and manufactured in compliance with the following standards:

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Safety requirements for electrical equipment for measurement, control and for laboratory use

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IEC/EN 61010-2-051

IEC/EN 61010-2-010

Electrical equipment for laboratory use

UL 61010-1

General requirement - Canadian electrical code

CAN/CSA-C22.2 No.61010-1

VELP reserves the right to modify the characteristics of its products with the aim to constantly improving their quality.

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The plug disconnects the instrument. Therefore, place the instrument where it can be quickly disconnected.

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Use only the power supply cable provided with the instrument.

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The heated solution may release toxic, dangerous or poisonous gases. Adequate safety measures must be taken, in accordance with the safety regulations in force, including the presence of hood and personal protective equipment (masks, gloves, goggles, etc.).

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The vessel must be made of a suitable material to withstand the foreseen temperature.

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Magnetic field doesn't interfere with the function of cardiac pacemakers or data media if they are farther than 20 cm to the instrument.

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Position the instrument on a flat surface, with a distance from the wall of 30 cm (at least).

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Do not use explosive and dangerous materials for which the equipment is not designed for. The stirrer must not be used in explosive atmospheres, in bain-marie and to stir combustible liquids that have a low combustion temperature. The product is intended for use with very small quantities of flammable liquids or flammable liquids that have a fire point higher than 935 °C and a flash point higher than 910 °C.

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The unit is fitted with two fuses (2xT5 A L 250 V (for 230V), 2xT8 A 250 V (for 115 V)), found in the socket on the back. To replace one or more disconnect the mains cable and, using a screwdriver, lift up the small cover on the fuse box.

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It is responsibility of the user appropriately decontaminate the instrument in case of dangerous substances fall on or in it accordingly to the safety datasheet of substances used and to the current laboratories safety standards. It is not possible to decontaminate the product under steam.

It is also responsibility of the user to use substances for cleaning or decontaminating which do not react with internal parts of the instrument or with the material contained in it. In case of doubts on the compatibility of a cleaning solution, contact the manufacturer or local distributor.

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1. Introduction

The HSC heating magnetic stirrer with ceramic plate is used to heat and mix liquids inside a suitable container placed on the ceramic plate. The stirring is possible thanks to a magnetic drive stir bar placed into the vessel.

The instrument is used in the laboratory for general use and for all those applications that require a regulation of the stirring speed and of the heating plate temperature.

The ceramic heating plate allows to reach higher temperature than a traditional aluminium heating plate; furthermore, the ceramic is an inert and very hard material, so it offers great resistance to almost any type of chemical or mechanical aggression keeping the surface characteristics unaltered over time.

The heating plate works up to 400 °C.

The motor is able to stir up to 1300 rpm.

Note: it's important to choose the most suitable magnetic stirring bar related to the quantity and type of the liquid to be stirred as well as to the type of the beaker. Size and shape of the magnetic stir bar determines the stirring efficiency at any given speed. The stirring bar which satisfies most applications is code A00000356 (8 x 40 mm).

2. Installation

- Unpacking
 - Check the integrity of the unit after unpacking
- The box includes
 - HSC, Heating Magnetic Stirrer
 - Instruction manual
 - Power supply cable
- Installation
 - Place the unit on a non-flammable surface
 - Make sure that the technical specification of the unit and the power supply line are the same
 - Make sure the electrical network is grounded
 - Make sure that the main switch is on position "0" (OFF)
 - Connect the instrument to an easily accessible socket (compliant with the current safety norms), using only the provided power supply cable
 - Note:** keep the power supply cable far away from the hot plate.
 - Place the container with the liquid and the stir bar on the stirring plate

3. Operation

Power-up

- Turn the unit on using the main switch

Stirring

- Rotate the right knob to set the stirring speed
- The dedicated LED is on when the stirring function is active

Heating

- Rotate the left knob to set the heating plate temperature
- The dedicated LED is on when the heating element is powered

4. Maintenance

No routine or extraordinary maintenance is necessary except periodically cleaning.

Cleaning

Disconnect the unit from the power supply and use a cloth dampened with a non-flammable non-aggressive detergent.

Repair

Repairs must be carried out by authorized Velp personnel only.

The instrument must be transported using the original shockproof packaging when freight forwarders, couriers or other are taking care of the unit transport. Follow the indications on the original packaging (e.g. place the instrument on a pallet).

It is the responsibility of the user, to properly decontaminate the unit in case of hazardous substances remaining on the surface or interior of the device. If in doubt about the compatibility of a cleaning or decontamination product, contact the manufacturer or distributor.

5. Technical data

General features	Model	F20500552	F20510552
	Voltage	230 V – 50/60 Hz	115 V – 60 Hz
	Power input	800 W	
	Current consumption	3.5 A	7 A
	Dimensions (WxHxD)	203 x 94 x 344 mm (8 x 3.7 x 13.5 in)	
	Weight	3.2 kg (7 lb)	
	Construction material	Technopolymer	
	Working in continuous	Admitted	
	Maximum load on the plate	25 kg	
	Noisiness	<< 80 dBA	
	Environmental temperature admitted	+5...+40 °C	
	Storage temperature admitted	-10...+60 °C	
	Max humidity	80%	
	Overvoltage category	II	
Heating	Pollution degree CEI EN61010-1	2	
	Max altitude	4000 m	
	Heating plate power output	770 W	
	Heating plate dimensions	180 x 180 mm (7 x 7 in)	
Stir	Temperature range	0 ÷ 400 °C	
	Heating plate construction material	Ceramic	
	Stirring capacity	15 L H ₂ O	
	Speed range	0 ÷ 1300 rpm	
	Motor rating output	10 W	

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
We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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EN 61010-1	Safety requirements for electrical equipment for measurement, control, and laboratory use Part 1: General requirements
EN 61010-2-010	Particular requirements for laboratory equipment for the heating of material
EN 61010-2-051	Particular requirements for laboratory equipment for mixing and stirring
EN 61326-1	Electrical equipment for measurement, control and laboratory use - EMC requirements - Part 1: General requirements

and satisfies the essential requirements of the following directives:

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2006/42/EC	Machinery directive
2014/30/EU	Electromagnetic compatibility directive
2015/863/EU (RoHS III)	Restriction of the use of certain hazardous substances
2012/19/EU (WEEE)	Waste of electric and electronic equipment

plus modifications / più modifiche / plus modifications / más sucesivas modificaciones / in der jeweils gültigen Fassung /
加修改.

9. Declaration of conformity UK CA

We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following regulations:

S.I. 2016/1101 Electrical Equipment (Safety) Regulations 2016

S.I. 2016/1091 Electromagnetic Compatibility Regulations 2016

according to the relevant designated standards:

EN 61010-1 Safety requirements for electrical equipment for measurement, control, and laboratory use
Part 1: General requirements

EN 61010-2-010 Particular requirements for laboratory equipment for the heating of material

EN 61010-2-051 Particular requirements for laboratory equipment for mixing and stirring

EN 61326-1 Electrical equipment for measurement, control and laboratory use - EMC requirements -
Part 1: General requirements

and satisfies the essential requirements of regulations:

S.I. 2008/1597 Supply of Machinery (Safety) Regulations 2008

S.I. 2012/3032 Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic
Equipment Regulations 2012

S.I. 2013/3113 Waste Electrical and Electronic Equipment Regulations 2013

plus modifications.

Thank you for having chosen VELP!

Established in 1983, VELP is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. VELP works according to **ISO 9001**, **ISO 14001** and **OHSAS 18001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

Analytical instruments

Elemental Analyzers
Digestion Units
Distillation Units
Solvent Extractors
Fiber Analyzers
Dietary Fiber Analyzers
Oxidation Stability Reactor
Consumables

Laboratory Equipment

Magnetic Stirrers
Heating Magnetic Stirrers
Heating Plates
Overhead stirrers
Vortex mixers
Homogenizers
COD Thermoreactors
BOD and Respirometers
Cooled Incubators
Flocculators
Overhead Shakers
Turbidimeter
Radiation Detector
Open Circulating Baths
Pumps

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VELP Scientifica Srl