Instruction Manual
Manuale di istruzioni
Manuel d'instructions
Manual de instrucciones
Bedienungsanleitung
操作指南

AREC 7 Connect Heating Magnetic Stirrer

F20500550, F20510550

General Information / Informazioni Generali / Informations Générales / Información General / Allgemeine Hinweise / 基本信息

Before using the unit, please read the following instruction manual carefully.
Prima dell'utilizzo dello strumento si raccomanda di leggere attentamente il seguente manuale operativo.
Avant d'utiliser l'instrument, il est recommandé de lire attentivement le présent manuel d'instructions.
Antes de utilizar el instrumento, le recomendamos que lea con atención el siguiente manual de instrucciones.
Bitte lesen Sie vor Inbetriebnahme des Geräts diese Bedienungsanleitung sorgfältig durch.
使用前请仔细阅读以下使用说明书。

Caution, hot surface!
Attenzione, superficie calda!
Attention, surface chaude!
Prudencia, superficie caliente!
Vorsicht, heiße Oberfläche!
注意：加热面高温

Do not dispose of this equipment as urban waste, in accordance with EEC directive 2012/19/UE.
Non smaltire l'apparecchiatura come rifiuto urbano, secondo quanto previsto dalla Direttiva 2012/19/UE.
Ne pas recycler l'appareil comme déchet solide urbain, conformément à la Directive 2012/19/UE.
No tirar el equipo en los residuos urbanos, como exige la Directiva 2012/19/UE.
Dieses Gerät unterliegt der Richtlinie 2012/19/UE und darf nicht mit dem normalen Hausmüll entsorgt werden.
按照EEC指令2002/96/CE，不要将该设备作为城市垃圾处理。

The product can be used with flammable liquids.
Il prodotto può essere utilizzato con liquidi infiammabili.
Le produit peut être utilisé avec des liquides inflammables.
El producto puede utilizarse con líquidos inflamables.
Das Produkt kann mit brennbaren Flüssigkeiten verwendet werden.
本产品可与易燃液体混合使用。

This unit must be used for laboratory applications indoor only.
The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions. If the product is used in a not specified way by the manufacturer or with non specified accessories, product's safety may be compromised.

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IEC/EN 61010-1
IEC/EN 61010-2-051
IEC/EN 61010-2-010

Electrical equipment for laboratory use

UL 61010-1

General requirement - Canadian electrical code

CAN/CSA-C22.2 No.61010-1

VELP reserves the right to modify the characteristics of its products with the aim to constantly improving their quality.
Nell'impegno di migliorare costantemente la qualità dei prodotti, VELP si riserva la facoltà di variarne le caratteristiche.
Dans le but d'améliorer constamment la qualité de ses produits, VELP se réserve le droit d'apporter des modifications aux caractéristiques de ceux-ci.
VELP se reserva el derecho de modificar las características de productos con el fin de mejorar constantemente su calidad.
VELP behält sich zum Zwecke der ständigen Verbesserung der Produktqualität das Recht auf Änderung der Geräteeigenschaften vor.
VELP保留修改其产品特性的权利，以不断提高其质量。

Safety Regulations / Norme di Sicurezza / Consignes de Sécurité / Advertencias de Seguridad / Sicherheitshinweise / 安全法规

The plug disconnects the instrument. Therefore, place the instrument where it can be quickly disconnected.

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Use only the power supply cable provided with the instrument.

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The heated solution may release toxic, dangerous or poisonous gases. Adequate safety measures must be taken, in accordance with the safety regulations in force, including the presence of hood and personal protective equipment (masks, gloves, goggles, etc.).

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The vessel must be made of a suitable material to withstand the foreseen temperature.

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Magnetic field doesn't interfere with the function of cardiac pacemakers or data media if they are farther than 20 cm to the instrument.

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Position the instrument on a flat surface, with a distance from the wall of 30 cm (at least).

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Do not use explosive and dangerous materials for which the equipment is not designed for. The stirrer must not be used in explosive atmospheres, in bain-marie and to stir combustible liquids that have a low combustion temperature. The product is intended for use with very small quantities of flammable liquids or flammable liquids that have a fire point higher than 935 °C and a flash point higher than 910 °C.

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The unit is fitted with two fuses (2xT5 A L 250 V (for 230V), 2xT8 A 250 V (for 115 V)), found in the socket on the back. To replace one or more disconnect the mains cable and, using a screwdriver, lift up the small cover on the fuse box.

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It is responsibility of the user appropriately decontaminate the instrument in case of dangerous substances fall on or in it accordingly to the safety datasheet of substances used and to the current laboratories safety standards. It is not possible to decontaminate the product under steam.

It is also responsibility of the user to use substances for cleaning or decontaminating which do not react with internal parts of the instrument or with the material contained in it. In case of doubts on the compatibility of a cleaning solution, contact the manufacturer or local distributor.

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The AREC 7 Connect heating magnetic stirrer with ceramic plate is used to heat and mix liquids inside a suitable container placed on the ceramic plate. The stirring is possible thanks to a magnetic drive stir bar placed into the vessel.

The instrument is used in the laboratory for general use and for all those applications that require a precise regulation of the stirring speed and of the heating plate temperature.

The ceramic heating plate allows to reach higher temperature than a traditional aluminium heating plate; furthermore, the ceramic is an inert and very hard material, so it offers great resistance to almost any type of chemical or mechanical aggression keeping the surface characteristics unaltered over time.

The heating plate works up to 550 °C.

The powerful brushless motor is able to stir from 30 up to 1700 rpm. SpeedServo technology maintains constant speed as viscosity changes.


The instrument is equipped with advanced functions as auto-reverse of the stirring direction, intermittent mode programs and integrated timer to meet users' workflow and safety lock feature to avoid possible alteration of the working conditions.

The AREC 7 Connect can work in combination with the Pt100 external probe for precise and accurate control of the sample temperature.

The integrated Wi-Fi module transmits and receives information in order to monitor and control the unit through the VELP Ermes cloud platform.

- 1 Ceramic heating plate
- 2 Main switch
- 3 Speed control knob
- 4 Temperature control knob
- 5 Display

- A Probe led
- B Wi-Fi led
- C Temperature setpoint led
- D Timer led
- E Autoreverse led
- F Intermittent mode led
- G Menu and lock button

Note: it's important to choose the most suitable magnetic stirring bar related to the type and the quantity of the liquid to be stirred as well as to the type of the beaker. Size and shape of the magnetic stir bar determines the stirring efficiency at any given speed. The stirring bar which satisfies most applications is code A00000356 (8 x 40 mm).

2. Installation

- Unpacking
 - Check the integrity of the unit after unpacking
- The box includes
 - AREC 7 Connect, Heating Magnetic Stirrer
 - Instruction manual
 - Magnetic stir bar
 - Power supply cable
- Installation
 - Place the unit on a non-flammable surface
 - Make sure that the technical specification of the unit and the power supply line are the same
 - Make sure that the electrical network is grounded
 - Make sure that the main switch is on position "0" (OFF)
 - Connect the instrument to an easily accessible socket (compliant with the current safety norms), using only the provided power supply cable
 - Note:** *keep the power supply cable far away from the hot plate.*
 - Place the container with the liquid and the stir bar on the stirring plate

3. Operation

Power-up

- Turn the unit on using the main switch
- The display shows the software version and the last setpoint values

If the set starting mode is *StOP* (see chap. 4), the heating and the stirring functions are not active and both displays show OFF. Otherwise, the unit starts working at the last set temperature and speed setpoints.

Stirring

- Rotate the right knob to set the stirring speed
- Push the knob to start stirring
- The display shows the current motor speed
- Push the knob to stop stirring

Heating

- Rotate the left knob to set the heating plate temperature
- Push the knob to start heating
- The display shows the current heating plate temperature for 3 seconds and the setpoint for 1 second

Note: *when the display shows the temperature setpoint, the dedicated led (C) is on.*

- Push the knob to stop stirring

Note: *if the heating function is not active and the heating plate temperature is more than 50 °C, the display shows the blinking message "Hot". This warning is not active if the instrument is not powered.*

Timer, autoreverse and intermittent mode functions

Press the central button repeatedly to access in sequence the following function menus:

- Timer
 - Rotate the right knob to set the timer (hours/minutes)
 - Press it to confirm

Note: when the timer is set, the dedicated led (D) is on. When the countdown starts, the led is blinking.

Note: select HHnn to disable the function.

After the timer has elapsed, the heating will be stopped (if *Time Heat* is OFF, see chap. 4), the stirring will be stopped (if *Time Sped* is OFF, see chap. 4) and the display will show *End*.

- Autoreverse
 - Rotate the right knob to set the autoreverse time (minutes/seconds)
 - Press it to confirm

Note: when the autoreverse is set, the dedicated led (E) is on.

Note: select nnSS to disable the function.

The rotation direction will change after each set interval time.

- Intermittent mode
 - Rotate the right knob to set the run time of the intermittent mode (minutes/seconds)
 - Press it to confirm
 - Rotate the knob to set the pause time of the intermittent mode (minutes/seconds)
 - Press it to confirm

Note: when the intermittent mode is set, the dedicated led (F) is on during the run time and it is blinking during the pause time.

Note: select nnSS to disable the function.

The stirring function will be interrupted cyclically.

To exit the function menu, push the central button again or do not press any buttons/knobs for 5 seconds.

Note: each function can be set individually or in combination with the others.

Lock function

- While the unit is working, push the central button for 3 seconds to enable the lock function.

If the lock function is enabled, it is not possible to change the work conditions.

The led around the central button blinks if the knobs or the central button are pressed/rotated.

- Push the central button for 3 seconds to disable the function.

External probe connection

- Turn off the instrument using the main switch
- Screw the threaded support rod into its seat on the back of the instrument (optional)
- Fasten the clamp onto the support rod (optional)
- Place the external temperature probe into the clamp (optional); then place it into the sample contained inside the flask
- Plug the probe into the dedicated socket on the back of the instrument
- Turn the unit on using the main switch

When the external probe is connected, the dedicated led (A) is on.

The heating function works as described above but the settable and showed temperature will be that read by the external probe.

4. Setting mode

Press both knobs for 3 seconds when heating and stirring are OFF to access the setting mode.

Use the “menu” buttons to browse through the available menus, listed below.

Press the speed control knob to access the selected menu. Then, turn the knob to modify the value of the parameter (when it is possible). Press the knob to confirm.

To exit the setting mode, do not press any buttons/knobs for 10 seconds or press both knobs at the same time.

Text shown		Default value	Range	Menu's name and description
Display 1	Display 2			
TENP	LINt	550 °C	OFF, 50 ÷ 550 °C 10 °C steps	<u>Heating plate temperature limit</u> It limits the maximum value of temperature setpoint for the heating plate. If OFF is selected, the heating function is “disabled”.
SPEd	LINt	1700	100 ÷ 1700 rpm 100 rpm steps	<u>Stirring speed limit</u> It limits the maximum value of speed setpoint for the motor.
Strt	NOdE	StOP	StOP – run	<u>Starting mode</u> ➤ StOP: when the instrument is powered, the heating and the stirring functions are not active and both displays show OFF ➤ run: when the instrument is powered, it starts working at the last temperature and speed setpoints
Ctrl	tyPE	FINE	FinE - FAST	<u>Thermoregulation mode with external probe</u> ➤ FINE: good and optimized temperature control, minimized overshoot and oscillations, slow rise in temperature ➤ FAST: fast rise in temperature, increased overshoot
PrOb	ALar	On	On - OFF	<u>External probe safety alarm</u> It allows to enable / disable AL 6 and AL 7 error messages (see chap. 5).
PrOb	CAL	0.0 °C	-10 ÷ +10 °C 0,1 °C steps	<u>Probe calibration</u> It allows to align the external temperature probe reading to a reference thermometer.
tINE	Strt	SNAP	SNAP - SEtP	<u>Starting mode of timer countdown</u> ➤ SNAP: the timer countdown starts as soon as one function (heating, stirring) is active ➤ SEtP: the timer countdown starts when temperature setpoint is reached
tINE	SPEd	On	On - OFF	<u>Keep stirring at timer end</u> If On is selected, the stirring function remains active when the timer countdown ends.
tINE	HEAt	OFF	On - OFF	<u>Keep heating at timer end</u> If On is selected, the heating function remains active when the timer countdown ends.
ANBt	CAL	0.0 °C	-10 ÷ +10 °C 0,1 °C steps	<u>Heating plate probe calibration</u> It allows to align the heating plate probe reading to a reference thermometer.
COun	NOtO	---		<u>Motor and heating element operating time</u> It shows motor and heating element operating times. Operating times are showed in hours up to 9999 hours (around 416 days). Then they are showed in days (a digits' decimal point is turned on).
COun	HEAt			
UIFI		NO	YES – NO	<u>Wi-Fi connection</u> It allows to enable / disable the Wi-Fi connection.
rSEt		NO	YES – NO	<u>Reset</u> It allows to reset all the system settings to the default values.
OtC	cFG	---		<u>One time code</u> It shows the One Time Code during the registration of the device on the VELP Ermes cloud platform.
StAr	AP	NO	YES – NO	<u>Access point</u> It allows to enable / disable the Access Point mode of the Wi-Fi connection.

5. Error and warning messages

When the display shows an error message, the instrument's functions are stopped automatically.

AL1	Heating plate overtemperature	
AL2	Excessive heating time	
AL3	The stirring system doesn't run correctly	
AL4	External probe overtemperature	Only if the external probe is connected
AL6	Slow temperature increase read by the external probe	Only if the external probe is connected and if <i>PrOb ALAr</i> is On
AL7	Fast temperature decrease read by the external probe	

To remove the error message, restart the instrument.

If the alarm persists on the display, please contact VELP Scientifica's technical service department.

6. VELP Hermes cloud platform configuration

VELP Hermes is a revolutionary cloud platform that transforms and improves your laboratory experience.

VELP Hermes creates a connected ecosystem of devices, people and data that transmit information between each other cutting down distances and expanding your scientific potential.

VELP Hermes collects and stores your data with the maximum level of encryption following the highest cyber-security standards.

- Activate VELP Hermes account:
 - Go to www.velp.com through your internet browser
 - Create your VELP account or login if you already have it
 - Access to your profile page and select "Configure Your VELP Hermes Account".
Activate your VELP Hermes account after having accepted Terms and Conditions.
 - You will receive an email containing your ID, Password and a link to proceed to VELP Hermes authentication at the email address that you entered during the registration
 - At the first login on VELP Hermes platform, a new password will be requested
- Connect the heating magnetic stirrer to your Wi-Fi network:
 - Turn on the AREC 7 Connect
 - Enable the Wi-Fi connection and the access point mode as described above
 - Using a PC, tablet or smartphone select the *VELP_AREC_(Product Code)_(Serial Number)* network available on the Wi-Fi list
 - Open a web browser on your PC, tablet or smartphone; insert 192.168.4.1 into the address bar and press "Enter"
Enter "admin" as User ID and "velp_ermes" as Password
 - Configure the laboratory Wi-Fi network parameters. Then press "Push" to confirm.
NOTE: if necessary, contact your IT specialist to get all your laboratory Wi-Fi connection details
 - Wait until the Wi-Fi led stops blinking (it has to be constantly ON)
NOTE: if the led is still blinking after 2 minutes, please try to repeat the procedure and check that the network parameters are correct
- Register the heating magnetic stirrer on VELP Hermes platform:
 - On the AREC 7 Connect, access the One Time Code (OTC) setting menu as described above
 - On the PC, tablet or smartphone, login on VELP Hermes with your account and select "Device Registration"
 - Enter the AREC 7 Connect product code, serial number and purchase date.
Press "Next"
 - Enter the One Time Code showed on the display of the AREC.
Press "Next" and wait for the connection
NOTE: if the OTC is not showed on the display of the AREC, please try to repeat the procedure and check that the network parameters are correct
 - Give a name to the AREC 7 Connect and enter the address where the instrument is located
Press "Next"
 - Assign the AREC 7 Connect to an existing virtual laboratory or create a new one through "Laboratory Management"

At the end of the procedure, you will be able to monitor and control the connected heating magnetic stirrer through "Instruments" menu.

NOTE: the instrument can be connected only to a 2.4 GHz Wi-Fi network.

NOTE: for further information, please visit VELP website.

7. Maintenance

No routine or extraordinary maintenance is necessary except periodically cleaning.

Cleaning

Disconnect the unit from the power supply and use a cloth dampened with a non-flammable non-aggressive detergent.

Repair

Repairs must be carried out by authorized Velp personnel only.

The instrument must be transported using the original shockproof packaging when freight forwarders, couriers or other are taking care of the unit transport. Follow the indications on the original packaging (e.g. place the instrument on a pallet).

It is the responsibility of the user, to properly decontaminate the unit in case of hazardous substances remaining on the surface or interior of the device. If in doubt about the compatibility of a cleaning or decontamination product, contact the manufacturer or distributor.

8. Technical data

General features	Model	F20500550	F20510550
	Voltage	230 V – 50/60 Hz	115 V – 60 Hz
	Power input	800 W	
	Current consumption	3.5 A	7 A
	Dimensions (WxHxD)	203 x 94 x 344 mm (8 x 3.7 x 13.5 in)	
	Weight	3.3 kg (7.3 lb)	
	Construction material	Technopolymer	
	Working in continuous	Admitted	
	Maximum load on the plate	25 kg	
	Noisiness	<< 80 dBa	
	Environmental temperature admitted	+5...+40 °C	
	Storage temperature admitted	-10...+60 °C	
	Max humidity	80%	
	Overvoltage category	II	
	Pollution degree CEI EN61010-1	2	
	Max altitude	2000 m	
Heating	Heating plate power output	770 W	
	Heating plate dimensions	180 x 180 mm (7 x 7 in)	
	Temperature range	0 ÷ 550 °C (1 °C steps)	
	Temperature resolution	1 °C	
	Heating plate construction material	Ceramic	
External probe	Type	Pt100 Class A – Ø 3mm	
	Temperature range	0 ÷ 300 °C (1 °C steps)	
	Temperature resolution	1 °C	
	Accuracy	± 0.5 °C *	
	Electrical data	3.3 VDC – 1 W (max)	
Stir	Stirring capacity	20 L H ₂ O	
	Speed range	30 ÷ 1700 rpm (5 rpm steps)	
	Speed resolution	1 rpm	
	Motor type	BLDC	
	Autoreverse / Intermittent mode	5 s ÷ 99 min 59 s (1 s steps)	
	Motor rating output	10 W	
Timer	Countdown range	1 min ÷ 99 h 59 min (1 min steps)	
Counters	Motor	Operating hours	
	Heating element		
Wi-Fi	Frequency	2.4 GHz	
	Protocol	802.11 b/g/n	
	Power output	17.5 ÷ 20 dBm @ 11b mode	

* in the following conditions: 800 ml water in 1 liter glass beaker (diameter 105 mm), stirring bar 8 x 40 mm, 600 rpm, 50 °C.

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
We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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EN 61010-1	Safety requirements for electrical equipment for measurement, control, and laboratory use Part 1: General requirements
EN 61010-2-010	Particular requirements for laboratory equipment for the heating of material
EN 61010-2-051	Particular requirements for laboratory equipment for mixing and stirring
EN 62479	Assessment of the compliance of low-power electronic and electrical equipment with the basic restrictions related to human exposure to electromagnetic fields (10 MHz to 300 GHz)
ETSI EN 301 489-1 V2.2.3	Electromagnetic compatibility (EMC) standard for radio equipment and services; Part 1: Common technical requirements
ETSI EN 301 489-17 V3.2.4	ElectroMagnetic Compatibility (EMC) standard for radio equipment and services; Part 17: Specific conditions for Broadband Data Transmission Systems
EN 61326-1	Electrical equipment for measurement, control and laboratory use - EMC requirements - Part 1: General requirements
ETSI EN 300 328 V2.2.2	Wideband transmission systems; Data transmission equipment operating in the 2,4 GHz band

and satisfies the essential requirements of the following directives:

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2006/42/EC	Machinery directive
2014/53/EU	Radio Equipment directive (RED)
2015/863/EU (RoHS III)	Restriction of the use of certain hazardous substances
2012/19/EU (WEEE)	Waste of electric and electronic equipment

plus modifications / più modifiche / plus modifications / más sucesivas modificaciones / in der jeweils gültigen Fassung /
加修改.

Notified body: Nemko S.p.A. (NB 2051)

EU Type Certificate: 2051-RED-211502

12. Declaration of conformity ^{UK} _{CA}

We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following regulations:

S.I. 2017/1206 Radio Equipment Regulations 2017

according to the relevant designated standards:

EN 61010-1	Safety requirements for electrical equipment for measurement, control, and laboratory use Part 1: General requirements
EN 61010-2-010	Particular requirements for laboratory equipment for the heating of material
EN 61010-2-051	Particular requirements for laboratory equipment for mixing and stirring
EN 62479	Assessment of the compliance of low-power electronic and electrical equipment with the basic restrictions related to human exposure to electromagnetic fields (10 MHz to 300 GHz)
ETSI EN 301 489-1 V2.2.3	Electromagnetic compatibility (EMC) standard for radio equipment and services; Part 1: Common technical requirements
ETSI EN 301 489-17 V3.2.4	ElectroMagnetic Compatibility (EMC) standard for radio equipment and services; Part 17: Specific conditions for Broadband Data Transmission Systems
EN 61326-1	Electrical equipment for measurement, control and laboratory use - EMC requirements - Part 1: General requirements
ETSI EN 300 328 V2.2.2	Wideband transmission systems; Data transmission equipment operating in the 2,4 GHz band

and satisfies the essential requirements of regulations:

S.I. 2008/1597 Supply of Machinery (Safety) Regulations 2008

S.I. 2012/3032 Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012

S.I. 2013/3113 Waste Electrical and Electronic Equipment Regulations 2013

plus modifications.

Notified body: Nemko North America Inc. (NB 1622)

UK Type Certificate: 1622-UK-432035

Thank you for having chosen VELP!

Established in 1983, VELP is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. VELP works according to **ISO 9001**, **ISO 14001** and **ISO 45001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

Analytical instruments

Elemental Analyzers
Digestion Units
Distillation Units
Solvent Extractors
Fiber Analyzers
Dietary Fiber Analyzers
Oxidation Stability Reactor
Consumables

Laboratory Equipment

Magnetic Stirrers
Heating Magnetic Stirrers
Heating Plates
Overhead stirrers
Vortex mixers
Dispersers
COD Thermoreactors
BOD and Respirometers
Cooled Incubators
Flocculators
Overhead Shakers
Turbidimeter
Open Circulating Baths
Pumps

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