Instruction Manual
Manuale di istruzioni
Manuel d'instructions
Manual de instrucciones
Bedienungsanleitung
操作指南

ECO 6 Thermoreactor

F10100120, F10110120

ECO 16 Thermoreactor

F10100126, F10110126

General Information / Informazioni Generali / Informations Générales / Información General / Allgemeine Hinweise / 基本信息

Before using the unit, please read the following instruction manual carefully.

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Attention, surface chaude!

Prudencia, superficie caliente!

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Do not dispose of this equipment as urban waste, in accordance with EEC directive 2012/19/UE.

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This unit must be used for laboratory applications indoor only.
The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions. If the product is used in a not specified way by the manufacturer or with non specified accessories, product's safety may be compromised.

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The plug disconnects the instrument. Therefore, place the instrument where it can be quickly disconnected.

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Use only the power supply cable provided with the instrument.

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The values indicated on the rating plate of the instrument must correspond to those of the power supply.

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The heated solution may release toxic, dangerous or poisonous gases. Adequate safety measures must be taken, in accordance with the safety regulations in force, including the presence of hood and personal protective equipment (masks, gloves, goggles, etc.).

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The vessel must be made of a suitable material to withstand the foreseen temperature.

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Position the instrument on a flat surface, with a distance from the wall of 30 cm (at least).

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It is responsibility of the user appropriately decontaminate the instrument in case of dangerous substances fall on or in it accordingly to the safety datasheet of substances used and to the current laboratories safety standards. It is not possible to decontaminate the product under steam.

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The ECO 6/16 thermoreactor is used for Chemical Oxygen Demand (COD) analysis.

The instrument allows simple and easy programming of time and temperature.

The working time can be set up to 999 minutes.

The working temperature can be set up to 200 °C for ECO 6 and up to 165 °C for ECO 16.

Through the digital display, it is possible to monitor the test at any time, with the analysis time constantly shown.

The unit automatically shuts down when the analysis time expires, whilst safety is of primary importance.

The ECO 6 offers the possibility to process 6 test tubes Ø 42 mm.

Dedicated accessories allow using 18 test tubes Ø 16 mm or 6 test tubes Ø 22 mm.

The ECO 16 offers the possibility to process 14 tubes Ø 16 mm and 2 test tubes Ø 22 mm.

2. Installation

- Unpacking
 - Check the integrity of the unit after unpacking
- The box includes
 - ECO 6/16 thermoreactor
 - Instruction manual
 - Power supply cable
- Installation
 - Place the unit on a non-flammable surface
 - Make sure that the technical specification of the unit and the power supply line is the same
 - Make sure the electrical network is grounded
 - Make sure that the main switch is on position "0" (OFF)
 - Connect the instrument to an easily accessible socket (compliant with the current safety norms), using only the provided power supply cable

Note: keep the power supply cable far away from the heating block.


3. Operation

Power-up

- Turn the unit on using the main switch.
- The display shows the instrument model and the software version.

- Then, it shows the main screen with the current heating block temperature on the left and the set program on the right.

Analysis

- Place the test tubes in the heating block.
- Set and select the desired program, as explained in the chapter 4.
- Press the ENTER button to start the analysis.
- The instrument carries out automatically the analysis:
 - The block is heated up to the set temperature
 - The temperature is maintained stable for the set timeDuring the analysis, the display constantly shows the current and set working parameters

Note: once the set temperature is reached, by pressing ENTER it is possible to activate the continuous mode.

- After the analysis time has elapsed, the heating is stopped and the display shows END
Press any buttons to return to the main screen

Note: the analysis can be manually interrupted by pressing ESC at any time; the instrument will ask for confirmation

4. Menu

From the main screen, press the DOWN button to access the main menu.
Use the UP and DOWN buttons to browse through the available menus, listed below.
Press the ENTER button to access the selected menu.
Then, use the UP and DOWN buttons to modify the value of the parameter (when it is possible).
Press the ENTER button twice to confirm.

To exit the menu, press the ESC button.

Note: please see the dedicated paragraph for the program setting.

Menu	Default value	Range	Description
PROGRAMS	---		It allows setting the working parameters. Please see the dedicated paragraph.
LANGUAGE	English	English – Italiano	It allows setting the software language.
TEMPERATURE UNIT	Celsius	Celsius – Fahrenheit	It allows setting the temperature unit.
SOUNDS	Yes	Yes – No	It allows enabling/disabling acoustic signals.
BLACKOUT	Yes	Yes – No	It allows enabling / disabling the blackout function. If Yes is selected and a blackout occurs, when the power supply is restored, the working cycle restarts automatically.
COUNTERS	---		It shows the number of analyses performed and the heating element operating time.
SERVICE	---		Only for authorized technicians.

Programs setting

From the PROGRAMS menu, use the UP and DOWN buttons to browse through the 20 available programs.
Press the ENTER button to access the selected program.
Press the ENTER button to select the program for the analysis.
Press the DOWN button to modify it:

- Press the ENTER button
- Use the UP and DOWN buttons to modify the working temperature; then press the ENTER button to confirm.
- Use the UP and DOWN buttons to modify the analysis time; then press the ENTER button to confirm.

5. Error messages

When the display shows an error message, the instrument's functions are stopped automatically.

AL2	Temperature probe broken
AL4	Overtemperature

To remove the error message, restart the instrument.
If the alarm persists on the display, please contact VELD Scientifica technical service department.

6. Maintenance

No routine or extraordinary maintenance is necessary except periodically cleaning.

Cleaning

Disconnect the unit from the power supply and use a cloth dampened with a non-flammable non-aggressive detergent.

Repair

Repairs must be carried out by authorized Velp personnel only.

The instrument must be transported using the original shockproof packaging when freight forwarders, couriers or other are taking care of the unit transport. Follow the indications on the original packaging (e.g. place the instrument on a pallet).

It is the responsibility of the user, to properly decontaminate the unit in case of hazardous substances remaining on the surface or interior of the device. If in doubt about the compatibility of a cleaning or decontamination product, contact the manufacturer or distributor.

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
We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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EN 61010-1	Safety requirements for electrical equipment for measurement, control, and laboratory use Part 1: General requirements
EN 61010-2-010	Particular requirements for laboratory equipment for the heating of material
EN 61326-1	Electrical equipment for measurement, control and laboratory use - EMC requirements - Part 1: General requirements

and satisfies the essential requirements of the following directives:

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2006/42/EC	Machinery directive
2015/863/EU (RoHS III)	Restriction of the use of certain hazardous substances
2012/19/EU (WEEE)	Waste of electric and electronic equipment

plus modifications / più modifiche / plus modifications / más sucesivas modificaciones / in der jeweils gültigen Fassung./
加修改.

11. Declaration of conformity

We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following regulations:

S.I. 2016/1101 Electrical Equipment (Safety) Regulations 2016

S.I. 2016/1091 Electromagnetic Compatibility Regulations 2016

according to the relevant designated standards:

EN 61010-1 Safety requirements for electrical equipment for measurement, control, and laboratory use
Part 1: General requirements

EN 61010-2-010 Particular requirements for laboratory equipment for the heating of material

EN 61326-1 Electrical equipment for measurement, control and laboratory use - EMC requirements -
Part 1: General requirements

and satisfies the essential requirements of regulations:

S.I. 2008/1597 Supply of Machinery (Safety) Regulations 2008

S.I. 2012/3032 Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment
Regulations 2012

S.I. 2013/3113 Waste Electrical and Electronic Equipment Regulations 2013

plus modifications.

Thank you for having chosen VELP!

Established in 1983, VELP is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. VELP works according to **ISO 9001**, **ISO 14001** and **OHSAS 18001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

Analytical instruments

Elemental Analyzers
Digestion Units
Distillation Units
Solvent Extractors
Fiber Analyzers
Dietary Fiber Analyzers
Oxidation Stability Reactor
Consumables

Laboratory Equipment

Magnetic Stirrers
Heating Magnetic Stirrers
Heating Plates
Overhead stirrers
Vortex mixers
Homogenizers
COD Thermoreactors
BOD and Respirometers
Cooled Incubators
Flocculators
Overhead Shakers
Turbidimeter
Radiation Detector
Open Circulating Baths
Pumps

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