Instruction Manual
Manuale di istruzioni
Manuel d'instructions
Manual de instrucciones
Bedienungsanleitung
操作指南

AREC 4

F20500670, F20510670

AREC 7

F20500680, F20510680, F20580680

General Information / Informazioni Generali / Informations Générales / Información General / Allgemeine Hinweise / 基本信息

Before using the unit, please read the following instruction manual carefully.
Prima dell'utilizzo dello strumento si raccomanda di leggere attentamente il seguente manuale operativo.
Avant d'utiliser l'instrument, il est recommandé de lire attentivement le présent manuel d'instructions.
Antes de utilizar el instrumento, le recomendamos que lea con atención el siguiente manual de instrucciones.
Bitte lesen Sie vor Inbetriebnahme des Geräts diese Bedienungsanleitung sorgfältig durch.
使用前请仔细阅读以下使用说明书。

Caution, hot surface!
Attenzione, superficie calda!
Attention, surface chaude!
Prudencia, superficie caliente!
Vorsicht, heiße Oberfläche!
注意：加热面高温

Do not dispose of this equipment as urban waste, in accordance with EEC directive 2012/19/UE.
Non smaltire l'apparecchiatura come rifiuto urbano, secondo quanto previsto dalla Direttiva 2012/19/UE.
Ne pas recycler l'appareil comme déchet solide urbain, conformément à la Directive 2012/19/UE.
No tirar el equipo en los residuos urbanos, como exige la Directiva 2012/19/UE.
Dieses Gerät unterliegt der Richtlinie 2012/19/UE und darf nicht mit dem normalen Hausmüll entsorgt werden.
按照 EEC 指令 2002/96/CE，不要将该设备作为城市垃圾处理。

The product can be used with flammable liquids.
Il prodotto può essere utilizzato con liquidi infiammabili.
Le produit peut être utilisé avec des liquides inflammables.
El producto puede utilizarse con líquidos inflamables.
Das Produkt kann mit brennbaren Flüssigkeiten verwendet werden.
本产品可与易燃液体混合使用。

This unit must be used for indoor laboratory applications only.
The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions. If the product is used in a way not specified by the manufacturer or with non specified accessories, product's safety may be compromised.

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IEC/EN 61010-1
IEC/EN 61010-2-051
IEC/EN 61010-2-010

Electrical equipment for laboratory use

UL 61010-1

General requirement - Canadian electrical code

CAN/CSA-C22.2 No.61010-1

VELP reserves the right to modify the characteristics of its products with the aim to constantly improve their quality.
Nell'impegno di migliorare costantemente la qualità dei prodotti, VELP si riserva la facoltà di variarne le caratteristiche.
Dans le but d'améliorer constamment la qualité de ses produits, VELP se réserve le droit d'apporter des modifications aux caractéristiques de ceux-ci.
VELP se reserva el derecho de modificar las características de productos con el fin de mejorar constantemente su calidad.
VELP behält sich zum Zwecke der ständigen Verbesserung der Produktqualität das Recht auf Änderung der Geräteeigenschaften vor.
VELP 保留修改其产品特性的权利，以不断提高其质量。

Safety Regulations / Norme di Sicurezza / Consignes de Sécurité / Advertencias de Seguridad / Sicherheitshinweise / 安全法规

The plug disconnects the instrument. Therefore, place the instrument where it can be quickly disconnected.

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Use only the power cable provided with the instrument.

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The heated solution may release toxic, dangerous or poisonous gases. Adequate safety measures must be taken, in accordance with the safety regulations in force, including the presence of hood and personal protective equipment (masks, gloves, goggles, etc.).

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The vessel must be made of a suitable material to withstand the foreseen temperature.

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Magnetic field doesn't interfere with the function of cardiac pacemakers or data media if they are farther than 20 cm to the instrument.

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Position the instrument on a flat surface, with a distance from the wall of 30 cm (at least) and between instruments of 10 cm (at least).

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Do not use with explosive or dangerous materials for which the equipment is not designed. The stirrer must not be used in explosive atmospheres, in bain-marie or to stir combustible liquids that have a low combustion temperature. The product is intended for use with very small quantities of flammable liquids or flammable liquids that have a fire point higher than 625°C and a flash point higher than 600°C.

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It is responsibility of the user to appropriately decontaminate the instrument in case of dangerous substances fall on or in it accordingly to the safety datasheet of substances used and to the current laboratories safety standards. It is not possible to decontaminate the product under steam.

It is also responsibility of the user to use substances for cleaning or decontaminating which do not react with internal parts of the instrument or with the material contained in it. In case of doubts on the compatibility of a cleaning solution, contact the manufacturer or local distributor.

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The AREC 4 / AREC 7 heating magnetic stirrer with a ceramic plate is designed for general laboratory use and for applications requiring precise control of the stirring speed and heating plate temperature.

The ceramic heating plate means that higher temperatures can be obtained compared to the traditional aluminium heating plate and it takes less time to heat the sample. Ceramic is an inert and tough material resistant to almost any chemical or mechanical aggression; its surface characteristics remain unaltered over time.

The AREC 4 / AREC 7 is used to stir liquids using a magnetically driven stir bar placed inside the vessel. The powerful brushless motor can stir from 100 to 1500 rpm.

SpeedServo™ technology maintains constant speed as viscosity changes.

Note: *it's important to choose the most suitable magnetic stirring bar related to the type and the quantity of the liquid to be stirred as well as to the type of the beaker. Size and shape of the magnetic stir bar determines the stirring efficiency at any given speed. The stirring bar which satisfies most applications is code A00000356 (8 x 40 mm).*

2. Installation

- Unpacking
 - Check the integrity of the unit after unpacking
- The box includes
 - AREC 4 / AREC 7 heating magnetic stirrer
 - Instruction manual
 - Power cable
 - Magnetic stir bar Ø 8 x 40 mm (A00000356)
- Installation
 - Place the unit on a non-flammable surface
 - Make sure that the technical specification of the unit and of the power supply line are the same
 - Make sure that the electrical network is grounded
 - Make sure that the main switch is on position "0" (OFF)
 - Connect the instrument to an easily accessible socket (compliant with the current safety norms), using only the provided power cable

Note: *keep the power cable far away from the hot plate.*

- Place the container with the liquid and the stir bar on the stirring plate

3. Operation

Power-up

- Turn the unit on using the main switch
- The right led between the two knobs lights up

Stirring

- Rotate the right knob to set the stirring speed, the led over the speed knob lights up
- The real speed increases until setpoint achievement
- A microprocessor ensures constant speed even if the viscosity changes (SpeedServo™)
- Turn the knob to the 0 positions to stop stirring

Heating

- Rotate the left knob to set the heating plate temperature
- If the set temperature is higher than the actual one, the led over the temperature knob lights up
- When the heating plate temperature exceeds 50 °C, the left led between the two knobs lights up until the temperature drops below 50 °C

Note: *This warning is active even if the main switch is off but not if the instrument is not powered.*

- Turn the knob to the 0 positions to stop heating

4. Error and warning messages

When the control panel shows one of the following error messages, the instrument functions automatically stop.

Display indicator	Malfunction:
Flashing Heating led and Hot plate led ON	Heating plate overtemperature ($T > 570\text{ }^{\circ}\text{C}$)
Flashing Stirring led	The stirring system doesn't run correctly

To remove the error message, restart the instrument from the main switch.

If the alarm persists on the control panel, please contact VELP Scientifica's technical assistance service.

5. Maintenance

No routine or extraordinary maintenance is required except periodic cleaning.

Cleaning

Disconnect the unit from the power supply and use a cloth dampened with a non-flammable non-aggressive detergent.

Repair

Repairs must be carried out by authorized Velp personnel only.

The transport of the instrument by freight forwarders, couriers or others must be carried out using the original shockproof packaging. Follow any instructions on the original packaging (e.g. palletizing).

It is the responsibility of the user to properly decontaminate the unit in case of hazardous substances remaining on the surface or interior of the device. If in doubt about the compatibility of a cleaning or decontamination product, contact the manufacturer or distributor of cleaning solution

6. Technical data

AREC 4			
General features	Model	F20500670	F20510670
	Voltage	230 V 50/60 Hz	115 V 60 Hz
	Power input	280 W	
	Current consumption	1,2 A	2,4 A
	Dimensions (WxHxD)	160x100x270 mm (6,3x3,9x10,6 in)	
	Weight	1,9 kg (4,2 lb)	
	Construction material	Aluminum body – Technopolymer enclosure	
	Working in continuous	Admitted	
	Maximum load on the plate	25 kg	
	Noisiness	<< 80 dBa	
	Environmental temperature admitted	+5...+40 °C	
	Storage temperature admitted	-10...+60 °C	
	Max humidity	80%	
	Overvoltage category	II	
	Pollution degree CEI EN61010-1	2	
	Max altitude	4000 m	
Heating	Heating plate power output	250 W	
	Heating plate dimensions	100 x 100 mm (4"x4")	
	Temperature range	50 ÷ 550 °C	
	Heat control type	Analogic	
	Heating plate construction material	Ceramic	
	Overtemperature alarm	Over 570°C	
Stir	Stirring capacity	5 l H ₂ O	
	Speed range	100 ÷ 1500 rpm	
	Stir control type	Analogic	
	Motor type	BLDC	
	Motor rating output	10 W	

AREC 7				
General features	Model	F20500680	F20510680	F20580680
	Voltage	230 V – 50/60 Hz	115 V – 60 Hz	100 V – 50/60 Hz
	Power input	800 W		
	Current consumption	3,5 A	7 A	8 A
	Dimensions (WxHxD)	200x100x340 mm (7,87 x 3,94 x 13,39 in)		
	Weight	3,5 kg (7,7 lb)		
	Construction material	Aluminium body – Technopolymer enclosure		
	Working in continuous	Admitted		
	Maximum load on the plate	30 kg (66 lb)		
	Noisiness	<< 80 dBa		
	Environmental temperature admitted	+5...+40 °C		
	Storage temperature admitted	-10...+60 °C		
	Max humidity	80%		
	Overvoltage category	II		
	Pollution degree CEI EN61010-1	2		
	Max altitude	4000 m		
Heating	Heating plate power output	770 W		
	Heating plate dimensions	180 x 180 mm (7"x7")		
	Temperature range	0 ÷ 550 °C (1 °C step)		
	Temperature resolution	1 °C		
	Heating plate construction material	Ceramic		
	Hot plate alarm	Over 50°C		
	Overtemperature alarm	Over 570°C		
Stir	Stirring capacity	15 l H ₂ O		
	Speed range	100 ÷ 1500 rpm		
	Stir control type	Analogic		
	Motor type	BLDC		
	Motor rating output	10 W		

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
We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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EN 61010-1	Safety requirements for electrical equipment for measurement, control, and laboratory use Part 1: General requirements
EN 61010-2-010	Particular requirements for laboratory equipment for the heating of material
EN 61010-2-051	Particular requirements for laboratory equipment for mixing and stirring
EN 61326-1	Electrical equipment for measurement, control and laboratory use - EMC requirements - Part 1: General requirements

and satisfies the essential requirements of the following directives:

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2006/42/EC	Machinery directive
2014/30/EU	Electromagnetic compatibility directive
2015/863/EU (RoHS III)	Restriction of the use of certain hazardous substances
2012/19/EU (WEEE)	Waste of electric and electronic equipment

plus modifications / più modifiche / plus modifications / más sucesivas modificaciones / in der jeweils gültigen Fassung /
加修改.

11. Declaration of conformity ^{UK}CA

We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following regulations:

S.I. 2016/1101 Electrical Equipment (Safety) Regulations 2016

S.I. 2016/1091 Electromagnetic Compatibility Regulations 2016

according to the relevant designated standards:

EN 61010-1 Safety requirements for electrical equipment for measurement, control, and laboratory use
Part 1: General requirements

EN 61010-2-010 Particular requirements for laboratory equipment for the heating of material

EN 61010-2-051 Particular requirements for laboratory equipment for mixing and stirring

EN 61326-1 Electrical equipment for measurement, control and laboratory use - EMC requirements -
Part 1: General requirements

and satisfies the essential requirements of regulations:

S.I. 2008/1597 Supply of Machinery (Safety) Regulations 2008

S.I. 2012/3032 Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic
Equipment Regulations 2012

S.I. 2013/3113 Waste Electrical and Electronic Equipment Regulations 2013

plus modifications.

Thank you for having chosen VELP!

Established in 1983, VELP is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. VELP works according to **ISO 9001**, **ISO 14001** and **ISO 45001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

Analytical instruments

Elemental Analyzers
Digestion Units
Distillation Units
Solvent Extractors
Fiber Analyzers
Dietary Fiber Analyzers
Oxidation Stability Reactor
Consumables

Laboratory Equipment

Magnetic Stirrers
Heating Magnetic Stirrers
Heating Plates
Overhead stirrers
Vortex and Shakers
Dispersers
COD Thermoreactors
BOD and Respirometers
Cooled Incubators
Flocculators
Turbidimeter
Open Circulating Baths
Pumps

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