Instruction Manual
Manuale di istruzioni
Manuel d'instructions
Manual de instrucciones
Bedienungsanleitung
指导手册

JLT 4 - JLT 6 Flocculation Tester

F105A0108, F105A0109

General Information / Informazioni Generali / Informations Générales / Información General / Allgemeine Hinweise / 基本信息

Before using the unit, please read the following instruction manual carefully.
Prima dell'utilizzo dello strumento si raccomanda di leggere attentamente il seguente manuale operativo.
Avant d'utiliser l'instrument, il est recommandé de lire attentivement le présent manuel d'instructions.
Antes de utilizar el instrumento, le recomendamos que lea con atención el siguiente manual de funcionamiento.
Bitte lesen Sie vor Inbetriebnahme des Geräts diese Bedienungsanleitung sorgfältig durch.
在使用本装置之前, 请仔细阅读以下使用说明书。

Do not dispose of this equipment as urban waste, in accordance with EEC directive 2002/96/CE.
Non smaltire l'apparecchiatura come rifiuto urbano, secondo quanto previsto dalla Direttiva 2002/96/CE.
Ne pas recycler l'appareil comme déchet solide urbain, conformément à la Directive 2002/96/CE.
No tirar el aparato en los desechos urbanos, como exige la Directiva 2002/96/CE.
Dieses Gerät unterliegt der Richtlinie 2002/96/EG und darf nicht mit dem normalen Hausmüll entsorgt werden.
根据 EEC 指令 2002/96/CE, 请不要将本设备作为城市垃圾处理。

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IEC/EN 61010-1
IEC/EN 61010-2-051

Electrical equipment for laboratory use
General requirement - Canadian electrical code

UL 61010-1
CAN/CSA-C22.2 No.61010-1

VELP reserves the right to modify the characteristics of its products with the aim to constantly improve their quality.
Nell'impegno di migliorare costantemente la qualità dei prodotti, VELP si riserva la facoltà di variane le caratteristiche.
Dans le but d'améliorer constamment la qualité de ses produits, VELP se réserve le droit d'apporter des modifications aux caractéristiques de ceux-ci.
VELP se reserva el derecho de modificar las características de productos con el fin de mejorar constantemente su calidad.
VELP behält sich zum Zwecke der ständigen Verbesserung der Produktqualität das Recht auf Änderung der Geräteeigenschaften vor.
VELP 保留修改其产品特性的权利，以不断提高其质量。

Safety Regulations / Norme di Sicurezza / Consignes de Sécurité / Advertencias de Seguridad / Sicherheitshinweise / 安全法规

Use only the power supply provided with the instrument.
Utilizzare solo l'alimentatore fornito con lo strumento.
Utilisez uniquement le bloc d'alimentation fourni avec l'instrument.
Utilizar únicamente el dispositivo de alimentación suministrado con el instrumento.
Verwenden Sie nur das mit dem Instrument gelieferte Netzteil.
只能使用仪器随附的电源。

Position the instrument on a flat surface, with a distance from the wall of 30 cm (at least).
Posizionare lo strumento su superfici piane, ad una distanza dalle pareti di almeno 30 cm.
Positionner l'appareil sur une surface plat, avec une distance de la paroi de 30 cm (au moins).
Coloque la unidad sobre una superficie plana, con una distancia de la pared de 30 cm (por lo menos).
Stellen Sie das Gerät auf einer ebenen Fläche mit einem Abstand zur Wand von 30 cm (mindestens).
将仪器放置在一个水平平台上，与墙壁的距离至少为 30 厘米。

It is responsibility of the user to appropriately decontaminate the instrument in case of dangerous substances fall on or in it accordingly to the safety datasheet of substances used and to the current laboratories' safety standards. It is not possible to decontaminate the product under steam.
It is also responsibility of the user to use substances for cleaning or decontaminating which do not react with internal parts of the instrument or with the material contained in it. In case of doubts about the compatibility of a cleaning solution, contact the manufacturer or local distributor.

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
11. DECLARATION OF CONFORMITY **UK**
CA 38

During chemical-physical treatments of waters aimed at coagulating and precipitating pollutants, chemical coagulants are evaluated and their dosing is optimized by jar tests. Because wastewater has a variable composition, the chemical dosing and the test frequency are to be adapted to local conditions. JLT Series is characterized by a settable rotation speed, from 10 to 300 rpm, with the real speed shown on the relative display. The epoxy-painted steel structure is studied and tested to give the instrument excellent resistance against the attack of chemical agents and corrosion in general. The stainless-steel stirring rods are adjustable in height by a self-blocking chuck. The examined samples can be lighted by a back LED lamp controlled by a switch located on the front panel.

	Description	Function
1	Main switch	To turn the instrument ON
2	Light switch	To turn the back LED lamp ON
A	Speed setting buttons	To adjust the speed setpoint
B	Speed LED display	To show the speed value
C	Timer setting buttons	To adjust the timer value
D	Timer LED display	To show the timer countdown
E	SET button	To set working mode: continuous, timer (minutes or hours)

2. Assembly and installation

- Unpacking
 - Check the integrity of the unit after unpacking
- The box includes
 - JLT 4 / JLT 6 Flocculation Tester
 - Instruction manual
 - Power cable for power supply
 - Switching power supply 100-240V 50/60Hz / 12V
- Installation
 - Place the unit on a non-flammable surface
 - Make sure that the values on the rating plate correspond to those of the switching power supply
 - Make sure that the technical specifications of the switching power supply and power supply line are the same
 - Make sure that the *main switch* (1) is on position "0" (OFF)

- Power the unit using the provided power supply. Make sure that the socket is compliant with the current safety norms
- Lift the stirring rods to place the beakers containing the sample. Once the beakers are positioned, adjust the stirring rods to the correct height

3. Operation

- Power-up
 - Turn the unit on using the *main switch* (1)
 - The display shows the software version
 - The *speed LED display* (B) shows the speed setpoint. The instrument starts to work and reaches the setpoint
- Stirring
 - Using the *speed setting buttons* (A), it is possible to set the desired speed in the range OFF / 10÷300rpm. The progression speed is increased by holding down one of the two buttons, to reach the desired value more quickly
 - The real speed increases/decreases until setpoint achievement
 - A microprocessor ensures constant speed even if the viscosity changes (SpeedServo)
- Timer
 - Click the *SET button* (E). “tIM” appears on the *timer LED display* (D) and the motor stops rotating. The speed setpoint on the *speed LED display* (B) starts blinking. To exit timer mode, press the *SET button* (E) again. When the timer is not active, the *timer LED display* (D) shows “...” to indicate that the instrument is operating in continuous mode
 - Using the *timer setting buttons* (C), it is possible to set the desired time value. The progression speed is increased by holding down one of the two buttons, to reach the desired value more quickly
 - Press the *SET button* (E) to start the timer. The instrument reaches the speed setpoint and the *timer LED display* (D) shows the countdown. To exit timer mode during the countdown, keep pressing the *SET button* (E) for 3 seconds. The instrument disables the timer and works in continuous mode
 - During the countdown, it is always possible to change the speed and timer values using the appropriate buttons
 - When the countdown finishes, the motor stops rotating and an acoustic signal indicates the end of the stirring process
- Timer modality
 - The timer can operate with two different scales: minutes and hours
 - Click the *SET button* (E). “tIM” appears on the display (minutes time base). Using the *timer setting buttons* (C), it is possible to set the time value in minutes of work.
 - By holding down the *SET button* (E) for 2 seconds, “t.IM” appears on the display indicating that the timer base is now hours. Using the *timer setting button* (C), it is now possible to set the time value in hours of work. By holding down the *SET button* (E) for 2 seconds again, the timer scale returns to minutes.
 - The timer range in the two scales is:
 - “tIM”: from 1 to 999 minutes
 - “t.IM”: from 1 to 99 hours

4. Error and warning messages

When the display shows an error message, the instrument functions automatically stop.

AL1 | The motor doesn't start or not reach the set speed

To remove the error message, restart the instrument from the main switch.

If the alarm persists on the display, please contact VELP Scientifica's technical service department.

5. Maintenance

No routine or extraordinary maintenance is necessary for the unit.

Cleaning

Disconnect the unit from the power supply and use a cloth dampened with a non-flammable non-aggressive detergent. During cleaning operations, prevent water and humidity residues from entering the unit.

Repair

Repairs must be carried out by authorized Velp personnel only.

The transport of the instrument by freight forwarders, couriers or others must be carried out using the original shockproof packaging. Follow any instructions on the original packaging (e.g., palletizing).
It is the responsibility of the user to properly decontaminate the unit in case hazardous substances remain on the surface or interior of the device. If in doubt about the compatibility of a cleaning or decontamination product, contact the manufacturer or distributor.

6. Technical data

	Models	F105A0108	F105A0109
General	Power supply	100-240 V / 50-60 Hz	
	Dimensions (WXHxD)	655x404x296 mm (25.8x15.9x11.7 in)	935x404x296 mm (36.8x15.9x11.7 in)
	Weight	13 Kg 29 lb	17 Kg 38 lb
	Power consumption	19 W	
	Construction material	Epoxy painted metal structure	
	Noisiness	<< 80 dBa	
	Environmental temperature admitted	+5...+40 °C	
	Storage temperature admitted	-10...+60 °C	
	Max humidity	80%	
	Overvoltage category	II	
	Pollution degree CEI EN61010-1	2	
	Max altitude	2.000 m	
Stirring	Programmable speed range	10 ÷ 300 rpm	
	Motor type	BLDC	
	Speed selection	1 rpm step	
Safety	Alarm if motor can't start or not reach the set speed	Yes	
Timer	Working time	Minutes 1 - 999	
		Hours 1 - 99	

7. Accessories / Spare parts

00001000	Plastic beaker 1000 ml	A00001001	Glass beaker 1000 ml
A00001002	Transparent plastic Imhoff cone	A00001003	Glass graduated Imhoff cone
A00001004	Stand for two Imhoff cones	10000235	Foot

8. Analytical methods

Coagulation and flocculation of wastewaters

Chemical coagulants are added to wastewater for the primary purpose of removing suspended solids and also phosphorus. The addition of chemicals reduces also heavy metal concentrations and improves disinfection efficiency. Chemical coagulation can be obtained by adding to wastewater defined amounts of lime (calcium hydroxide), alum (aluminium sulphate) or iron salts (ferric or ferrous). The coagulation-sedimentation process typically involves:

- Injection and mixing of the coagulant that neutralizes the electric charge, prevalently negative, on suspended particles. Good and rapid mixing of the coagulant and wastewater is important to ensure efficient use of the chemical. Typical detention times in plant basins for turbulent mixing are 15 to 120 seconds
- Agglomeration of the coagulated particles into large settleable flocs. This agglomeration is accomplished by stirring the water slowly to allow the solids formed by the addition of chemicals to grow in size so that they can be removed by gravity settling. The stirring is usually performed by slowly rotating paddles for a period of 10 to 30 minutes.
- Sedimentation of the flocculated material by gravity separation. The solids are allowed to settle by gravity from the chemically treated standing wastewater. The effluent from the chemical settling tank is very clear and does not usually require further filtration.

The quantity of chemical coagulant required to perform the test varies with time and from wastewater to wastewater. Typical coagulant doses are 75 to 250 g/m³ for alum, 45 to 90 g/m³ for ferric chloride, and 200 to 400 g/m³ for lime.

Reagents

The choice and dosing of the chemical coagulant to be adopted for the removal of suspended solids from wastewater derive from the results of laboratory evaluations by the so-called Jar test. Multiple stirrers with reproducible stirring speeds allow to adopt standard conditions for the test, which are the basic requirement for reproducible results.

The standard conditions are related to:

- dimensions and shape of the vessel and stirring paddle
- time and speed of stirring during turbulent mixing
- evaluation criteria for the obtained results
- volume of wastewater
- time and speed of stirring during flocculation

The most diffused adopted conditions are:

- 1000 ml glass beakers (jars), tall form, Ø 105 mm
- stirring paddles 25 mm high, 75 mm wide, 1 mm thickness
- turbulent stirring after chemical addition: 120 rpm, 120 seconds
- first evaluation of results after 5 minutes of sedimentation
- 600 ml wastewater samples
- stirring height of paddles, the middle height of the sample
- slow speed flocculation: 30 rpm, 25 minutes

The results can be evaluated based on different criteria:

- A) Floc dimension evaluation with numerical degree (0=no flocs, 2=very small, 4=small, 6=medium, 8=large, 10=very large).
- B) Time from the addition of chemicals to the first appearance of flocs.
- C) Evaluation of residual turbidity of supernatant, after a determined sedimentation time, by turbidimeter.
- D) Measurement of the electro-kinetic potential of suspended particles on a sample taken immediately after the addition and mixing of chemicals. Coagulation should be easier when the zeta potential of particles is very low.
- E) Evaluation of filterability of clarified water by standardized membrane filters under pressure. The reduction of water flow with time is related to the degree of clogging of filters due to residual unsettled suspended matter.

Other measurements are related to the pH of water after chemical treatment (lime raises the value while iron and aluminium salts lower it) and the temperature of water at which the experiments are performed.

References: American Society for Testing Materials Norm ASTM D-2035-G4 T; Passino R. and M.Beccari (1970) "Standardization of jar test in coagulation-flocculation processes for the removal of inorganic turbidities (in Italian). *Acqua e Aria*, Milano, 1-10, October 1970"; Degremont (1978) "Etude de la coagulation et de la flocculation des eaux. Methode n.703, 948-950. *Memento technique de l'eau*. 1200 pages, Paris"; U.S. Association of Environmental Engineering Professors (1972) "Environmental Engineering Unit Operations and Unit Processes Laboratory Manual. J.T. O'Connor Ed. 350 ppg".

Leaching test

Leaching tests are used to simulate the behaviour of waste in a dump that is submitted to the leaching activity of rainwater. There are two different methods that use leaching solutions of different types:

- F) Acetic acid, proposed for dumps receiving both organic and inorganic wastes;
- G) Carbon dioxide-saturated water, for dumps receiving only inorganic wastes.

Waste sample to be examined: if a liquid phase is present, it must be removed by filtration or centrifugation and preserved in a refrigerator at 2-4 °C. The granulometry of the solid phase must be lower than 9.5 mm (standard sieve).

Waste extraction: 20-25 g of the solid phase are weighed with a precision of ± 0.1 g.

Method a) A volume of distilled water corresponding to 16 times the weight of the sample is added.

The pH of the mixture is maintained at 5 ± 0.2 using 0.5 N acetic acid. The stirring is performed for 24 hours, controlling the pH every 15 minutes during the first hour and then every hour. The pH is corrected to 5 ± 2 using again 0.5 N acetic acid. If after 24 hours of stirring the pH is higher than 5.2, a new correction is made and the stirring is prolonged for another 4 hours.

Method b) The leaching solution is given by a volume corresponding to 20 times the sample weight of distilled water saturated by carbon dioxide for not less than 15 minutes. The stirring is performed for 6 hours without any pH correction.

Stirring speed: must be able to maintain continuously mixed the solid phase with the leaching solution.

Analysis of extract: the leaching solution and the liquid phase if present in the original sample are filtered by a membrane filter of 0.45 μm porosity and submitted to analysis for the single metals by the methods commonly used for waters.

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Riferimenti: American Society for Testing Materials Norm ASTM D-2035-G4 T; Passino R. and M.Beccari (1970) "Standardization of jar test in coagulation-flocculation processes for the removal of inorganic turbidities (in Italian). Acqua e Aria, Milano, 1-10, October 1970"; Degremont (1978) "Etude de la coagulation et de la floculation des eaux. Methode n.703, 948-950. Memento technique de l'eau. 1200 pages, Paris"; U.S. Association of Environmental Engineering Professors (1972) "Environmental Engineering Unit Operations and Unit Processes Laboratory Manual. J.T. O'Connor Ed. 350 ppg".

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参考文献： American Society for Testing Materials Norm ASTM D-2035-G4 T; Passino R. and M.Beccari (1970) "Standardization of jar test in coagulation-flocculation processes for the removal of inorganic turbidities (in Italian). Acqua e Aria, Milano, 1-10, October 1970"; Degremont (1978) "Etude de la coagulation et de la floculation des eaux. Methode n.703, 948-950. Memento technique de l'eau. 1200 pages, Paris"; U.S. Association of Environmental Engineering Professors (1972) "Environmental Engineering Unit Operations and Unit Processes Laboratory Manual. J.T. O'Connor Ed. 350 pp".

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
We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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EN 61010-1	Safety requirements for electrical equipment for measurement, control, and laboratory use Part 1: General requirements
EN 61010-2-051	Particular requirements for laboratory equipment for mixing and stirring
EN 61326-1	Electrical equipment for measurement, control and laboratory use - EMC requirements - Part 1: General requirements

and satisfies the essential requirements of the following directives:

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2006/42/EC	Machinery directive
2014/30/EU	Electromagnetic compatibility directive
2015/863/EU (RoHS III)	Restriction of the use of certain hazardous substances
2012/19/EU (WEEE)	Waste of electric and electronic equipment

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11. Declaration of conformity ^{UK}CA

We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following regulations:

S.I. 2016/1101 Electrical Equipment (Safety) Regulations 2016

S.I. 2016/1091 Electromagnetic Compatibility Regulations 2016

according to the relevant designated standards:

EN 61010-1 Safety requirements for electrical equipment for measurement, control, and laboratory use
Part 1: General requirements

EN 61010-2-051 Particular requirements for laboratory equipment for mixing and stirring

EN 61326-1 Electrical equipment for measurement, control and laboratory use - EMC requirements -
Part 1: General requirements

and satisfies the essential requirements of regulations:

S.I. 2008/1597 Supply of Machinery (Safety) Regulations 2008

S.I. 2012/3032 Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment
Regulations 2012

S.I. 2013/3113 Waste Electrical and Electronic Equipment Regulations 2013

plus modifications.

Thank you for having chosen VELP!

Established in 1983, VELP is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. VELP works according to **ISO 9001**, **ISO14001** and **ISO 45001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

Analytical instruments

Elemental Analyzers
Digestion Units
Distillation Units
Solvent Extractors
Fiber Analyzers
Dietary Fiber Analyzers
Oxidation Stability Reactor
Consumables

Laboratory Equipment

Magnetic Stirrers
Heating Magnetic Stirrers
Heating Plates
Overhead stirrers
Vortexers and Shakers
Dispersers
COD Thermoreactors
BOD and Respirometers
Cooled Incubators
Flocculators
Turbidimeter
Open Circulating Baths
Pumps

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