Instruction Manual
Manuale di istruzioni
Manuel d'instructions
Manual de instrucciones
Bedienungsanleitung
指导手册

AM4 Heating Magnetic Stirrer, multiple positions

F20500420

AM4X Heating Magnetic Stirrer, multiple positions & VTF

F20500421

General Information / Informazioni Generali / Informations Générales / Información General / Allgemeine Hinweise / 一般信息

Before using the unit, please read the following instruction manual carefully.

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Do not dispose of this equipment as urban waste, in accordance with EEC directive 2002/96/CE.

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This unit must be used for laboratory applications only.

The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions.

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This unit has been designed and manufactured in compliance with the following standards:

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Safety requirements for electrical equipment for measurement, control and for laboratory use

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IEC/EN 61010-1

IEC/EN 61010-2-051

Electrical equipment for laboratory use

UL 61010-1

General requirement - Canadian electrical code

CAN/CSA-C22.2 No.61010-1

VELP reserves the right to modify the characteristics of its products with the aim to constantly improving their quality.

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The plug disconnects the instrument. Therefore, place the instrument where it can be quickly disconnected.

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The heated solution may release toxic, dangerous or poisonous gases. Adequate safety measures must be taken, in accordance with the safety regulations in force, including the presence of hood and personal protective equipment (masks, gloves, goggles, etc.).

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Beware of the effect of the magnetic field on cardiac pacemakers and data media.

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Position the instrument on a flat surface, with a distance from the wall of 30 cm (at least).

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Do not use with explosive and dangerous materials for which the equipment is not designed. The stirrer must not be used in explosive atmospheres, in bain-marie and to stir combustible liquids that have a low combustion temperature. The minimum fire point of flammable solution is 750 °C. Only small amounts (< 50 ml) of flammable liquid can be used with the device.

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It is responsibility of the user appropriately decontaminate the instrument in case of dangerous substances fall on or in it.

It is also responsibility of the user to use safety substances for cleaning or decontaminating, which do not react with internal parts of the instrument or with the material contained in it. In case of doubts on the compatibility of a cleaning solution, contact the manufacturer or local distributor.

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The application of new technology has led to the creation of these modern magnetic stirrers for laboratory use whose basic features are safety, reproducibility of results, high performance and energy saving. The structure of the unit is made out of aluminum treated with epoxy resins offering a high resistance to the many chemical aggressions typical of the laboratory environment. A wide central recess on both sides offers a good grip and facilitates bench-top handling whilst the handle moulded directly into the rear of the die-cast structure can be used to transport the instrument once it has been disconnected from the power supply. The structure is designed so that accidental liquid spills cannot reach the internal parts of the unit. The electrical connections are recessed into the rear of the instrument and offer optimum electrical protection of the external electrics in the case of liquid spills. The engineering of the unit facilitates internal inspection if technical assistance is required.

The heating plate is made out of aluminum alloy with special coating and is specially treated to ensure:

- Optimum heat distribution and a high specific power thanks to the circular configuration
- Optimum temperature homogeneity across the plate
- High resistance to thermal stress and thermal shock, chemical aggressions, scratches and surface abrasions.

Magnetic stirring is generated by the VELP Scientifica Patented Composite Magnet driven by an asynchronous mono-phase brushless motor in alternating current which offers a virtually unlimited duration.

NOTE: using the heating plate at high temperatures may cause discoloring. This does not alter the thermal, mechanical and chemical resistance of the plate in any way.

2. Assembly and installation

Check the integrity of the unit after unpacking. The box includes:

- Heating Magnetic Stirrer
- Instruction manual

2.1 Electrical connections

After having unpacked the instrument, place the unit on the laboratory bench.

Before connecting the instrument to the power supply, make sure that the values on the rating plate correspond to those of the power supply. The equipment shall be connected to a power socket with protective earth connection, using only the power supply cord provided with the instrument.

Ensure that the socket and the relative cut-off device conform to current safety norms and are easy to reach.

NOTE: the mains cable must remain far away from the hot plate.

2.2 Start-up

Rotate the speed (right) and temperature (left) knobs completely to the left. Place the flask containing the sample and a suitable magnetic stirring bar on the stirring plate. Then, set the speed and temperature by turning the dedicated knobs.

3. Operating controls

REGULATION KNOBS

The right knob on the front of the unit can be used for quick precision regulation of mixing speed from 100 to 1500 rpm. The left one is for temperature, up to 370 °C.

ON-OFF SWITCH

The on-off switch turns the unit on and off. If the switch is in the "OFF" position the unit is off; if the switch is in the "ON" position the unit is on.

Always turn the unit off after use.

3.1 Using the thermoregulator VTF (Available only for AM4X model)

Screw the threaded support rod (optional accessory, A00001069) into its seat on the back of the instrument and fasten the VTF thermoregulator onto the support rod. Place the temperature probe in the receptacle making sure that it is completely immersed in the sample. Connect the two instruments (Heating Magnetic Stirrer and VTF) by plugging the VTF into the dedicated socket on the back of the instruments. Select the operating temperature required on the VTF thermoregulator. Turn the temperature control knob on the front panel of the Heating Magnetic Stirrer to maximum. The heating magnetic stirrer always has primary control of the heating plate temperature. When using the VTF thermoregulator always select the max temperature on the Heating magnetic Stirrer. The temperature control function of the heating plate can also be used as a safety thermostat. In this case the maximum temperature of the heating plate will not exceed the temperature setting on the Heating magnetic Stirrer meaning that a longer heating time is required in order to reach the VTF thermoregulator temperature setting.

3.2 Error messages

The unit is fitted with safety devices which cut-off the power supply to the heating plate in the case of malfunctions.

Display

Flashing Stirrer LED (once every second)
Flashing Heating LED (once every second)
Flashing Heating LED (2 times every second)
Flashing Heating LED (once every 3 seconds)

Cause

Fault in the stirring system
Overheated heating plate or thermocouple open
Fault in the temperature reading circuit
Thermocouple wires reversed or white wire of potentiometer board not connected

Should any of the above occur, please contact your nearest VELP Scientifica service centre.

4. Maintenance

No routine or extraordinary maintenance is necessary apart from periodically cleaning the unit as described in this manual. In compliance with the product guarantee law, replace the fuses and repairs to our units must be carried out in our factory, unless previously agreed otherwise with local distributors and using Velp spare part only.

The instrument must be transported in its original packaging and any indications present on the original packaging must be followed (e.g. palletized).

It is the responsibility of the user, to properly decontaminate the unit in case of hazardous substances remaining on the surface or interior of the device. If in doubt about the compatibility of a cleaning or decontamination product, contact the manufacturer or distributor.

4.1 Cleaning

Disconnect the unit from the power supply and use a cloth dampened with an non-inflammable non-aggressive detergent.

5. Technical data

Power supply	230V / 50-60Hz (+/-10%)
Dimensions WxHxD	715x115x246 (28.1x4.5x9.7)
Weight	8.6 Kg (AM4) ; 9.0 Kg (AM4X)
Power input	630W x 4 (2520W)
Heat output	600W x 4 (2400W)
Diameter of the heating plate	135mm x 4
Programmable temperature range	50 – 370 °C
Type of temperature control	Analog
Overtemperature protection	Yes
Motor rating output	10W x 4
Stirring capacity	20 liters of H ₂ O x 4
Programmable speed range	100 – 1500 rpm
Counter-reaction	Constant speed even when the viscosity changes
Type of motor control	Electronic speed control
Temperature range	+5...+40 °C
Storage temperature range	-10...+60 °C
Max humidity	80%
Pollution degree CEI EN61010-1	2
Overvoltage category	II
Noise level	<< 80 dBa
Max altitude	4000 m

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
We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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2015/863/EU (RoHS III)

2012/19/EU (WEEE)

and satisfies the essential requirements of the following directives:

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9. Declaration of conformity

We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following regulations:

S.I. 2016/1101 Electrical Equipment (Safety) Regulations 2016

S.I. 2016/1091 Electromagnetic Compatibility Regulations 2016

according to the relevant designated standards:

EN 61010-1 Safety requirements for electrical equipment for measurement, control, and laboratory use
Part 1: General requirements

EN 61010-2-010 Particular requirements for laboratory equipment for the heating of material

EN 61010-2-051 Particular requirements for laboratory equipment for mixing and stirring

EN 61326-1 Electrical equipment for measurement, control and laboratory use - EMC requirements -
Part 1: General requirements

and satisfies the essential requirements of regulations:

S.I. 2008/1597 Supply of Machinery (Safety) Regulations 2008

S.I. 2012/3032 Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic
Equipment Regulations 2012

S.I. 2013/3113 Waste Electrical and Electronic Equipment Regulations 2013

plus modifications.

Thank you for having chosen VELP!

Established in 1983, VELP is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. VELP works according to **ISO 9001**, **ISO14001** and **OHSAS 18001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

Analytical instruments

Elemental Analyzers
Digestion Units
Distillation Units
Solvent Extractors
Fiber Analyzers
Dietary Fiber Analyzers
Oxidation Stability Reactor
Consumables

Laboratory Equipment

Magnetic Stirrers
Heating Magnetic Stirrers
Heating Plates
Overhead stirrers
Vortex mixers
Homogenizers
COD Thermoreactors
BOD and Respirometers
Cooled Incubators
Flocculators
Overhead Shakers
Turbidimeter
Radiation Detector
Open Circulating Baths
Pumps

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