Instruction Manual
Manuale di istruzioni
Manuel d'instructions
Manual de instrucciones
Bedienungsanleitung
操作指南

AM4 Digital PRO Heating Magnetic Stirrer, multi-positions
AM4 Digital PRO Agitatore Magnetico Riscaldante Multiposizione
AM4 Digital PRO Agitateur magnétique chauffant à positions multiples
AM4 Digital PRO Agitador Magnético Calefactor, posiciones múltiples
AM4 Digital PRO Magnetrührer mit Heizplatte, mehrere Positionen
AM4 数字式专业加热磁力搅拌器，多位置

F20500425
SA20500425

**General Information / Informazioni Generali / Informations Générales / Información General /
Allgemeine Hinweise / 基本信息**

Before using the unit, please read the following instruction manual carefully.
Prima dell'utilizzo dello strumento si raccomanda di leggere attentamente il seguente manuale operativo.
Avant d'utiliser l'instrument, il est recommandé de lire attentivement le présent manuel d'instructions.
Antes de utilizar el instrumento, le recomendamos que lea con atención el siguiente manual de funcionamiento.
Bitte lesen Sie vor Inbetriebnahme des Geräts diese Bedienungsanleitung sorgfältig durch
使用前请仔细阅读以下使用说明书。

Caution, hot surface! / Attenzione, superficie calda! / Attention, surface chaude! / Prudencia, superficie caliente! /
Vorsicht, heiße Oberfläche! / 注意：加热面高温

Do not dispose of this equipment as urban waste, in accordance with EEC directive 2002/96/CE.
Non smaltire l'apparecchiatura come rifiuto urbano, secondo quanto previsto dalla Direttiva 2002/96/CE.
Ne pas recycler l'appareil comme déchet solide urbain, conformément à la Directive 2002/96/CE.
No tirar el aparato en los desechos urbanos, como exige la Directiva 2002/96/CE.
Dieses Gerät unterliegt der Richtlinie 2002/96/EG und darf nicht mit dem normalen Hausmüll entsorgt werden.
按照 EEC 指令 2002/96/CE，不要将该设备作为城市垃圾处理。

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This unit has been designed and manufactured in compliance with the following standards:

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Safety requirements for electrical equipment for measurement, control and for laboratory use

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Electrical equipment for laboratory use

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General requirement - Canadian electrical code

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Allgemeine Anforderung - Canadian electrical code

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VELP reserves the right to modify the characteristics of its products with the aim to constantly improving their quality.

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The plug disconnects the instrument. Therefore, place the instrument where it can be quickly disconnected.

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The heated solution may release toxic, dangerous or poisonous gases. Adequate safety measures must be taken, in accordance with the safety regulations in force, including the presence of hood and personal protective equipment (masks, gloves, goggles, etc.).

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Beware of the effect of the magnetic field on cardiac pacemakers and data media.

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Position the instrument on a flat surface, with a distance from the wall of 30 cm (at least).

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Do not use with explosive and dangerous materials for which the equipment is not designed. The stirrer must not be used in explosive atmospheres, in bain-marie and to stir combustible liquids that have a low combustion temperature. The minimum fire point of flammable solution is 750 °C. Only small amounts (< 50 ml) of flammable liquid can be used with the device.

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It is responsibility of the user appropriately decontaminate the instrument in case of dangerous substances fall on or in it.

It is also responsibility of the user to use safety substances for cleaning or decontaminating, which do not react with internal parts of the instrument or with the material contained in it. In case of doubts on the compatibility of a cleaning solution, contact the manufacturer or local distributor.

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AM4 Digital PRO heating magnetic stirrer with four aluminium alloy plate is used in laboratory for general use and for all those applications that require precise stirring speed regulation and temperature of the heating plate temperature. Aluminum structure is treated with epoxy resins offering high resistance to the most chemical aggressions typical of laboratory environment. The structure is designed so that accidental liquid spills cannot reach the internal parts of the unit. The electrical connections are recessed into the rear of the instrument and offer optimum electrical protection of the external electrics in the case of liquid spills. The engineering of the unit facilitates internal inspection if technical assistance is required.

The heating plate is made out of aluminum alloy with special coating and it is specially treated to ensure:

- Optimum heat distribution and a high specific power thanks to the circular configuration
- Optimum temperature homogeneity across the plate
- High resistance to thermal stress and thermal shock, chemical aggressions, scratches and surface abrasions.

Magnetic stirring is generated by the VELP Scientifica Patented Composite Magnet driven by an asynchronous mono-phase brushless motor in alternating current which offers a virtually unlimited duration.

NOTE: If the heating function is used the vessel must be made of a suitable material to withstand the foreseen temperature.

NOTE: It is therefore advisable to choose the most suitable magnetic stir bar in relation to the quantity and characteristics of the liquid and the type of vessel used.

NOTE: using the heating plate at high temperatures may cause discoloring. This does not alter the thermal, mechanical and chemical resistance of the plate in any way.

2. Assembly and installation

Check the integrity of the unit after unpacking. The box includes:

- Heating Magnetic Stirrer
- Instruction manual

2.1 Electrical connections

After having unpacked the instrument, place the unit on the laboratory bench.

Before connecting the instrument to the power supply, make sure that the values on the rating plate correspond to those of the power supply. The equipment shall be connected to a power socket with protective earth connection, using only the power supply cord provided with the instrument.

Ensure that the socket and the relative cut-off device conform to current safety norms and are easy to reach.

NOTE: the mains cable must remain far away from the hot plate.

2.2 Start-up

Rotate the speed (right) and temperature (left) knobs to the OFF position. Place the flask containing the sample and a suitable magnetic stirring bar on the stirring plate. Then, set the speed and temperature by turning the dedicated knobs.

3. Operating controls

3.1 Turning on and regulating

Connect the unit to mains and turn it on using the on-off switch. In the next 5 seconds, displays shows the software version and the restarting mode (A or b, chapter 3.4 "Setting mode").

Start the stirring function by turning the stirrer knob "Stirrer rpm" on the front panel. Speeds from 50 to 1500 rpm (10 rpm resolution) can be selected. The speed shown flashes until it reaches the set-point. A microprocessor ensures constant speed even when the viscosity changes (counter-reaction).

To start the heating function turn the heating knob "Heating °C" on the front panel. Temperatures from 5 to 370 °C can be selected at 5 °C intervals. The display shows the real temperature but every 5 seconds temperature set point is shown for 1 second, indicated by the central led on the display "Set" turned on.

NOTE: in case of black-out, once the power is back the device will restart in the restarting mode set (Opt tYpE, chapter 3.4). When heating is turned Off (Heating knob OFF) and the heating plate temperature exceeds 50 °C, the message "Hot" flashes on the display until the temperature of the heating plate falls to below 50 °C (or until the instrument is turned off using the ON/OFF switch).

NOTE: this warning is not active if the instrument is not powered.

3.2 Using the thermoregulator VTF

Screw the threaded support rod (optional accessory, A00001069) into its seat on the back of the instrument and fasten the VTF thermoregulator onto the support rod. Place the temperature probe in the receptacle making sure that it is completely inserted in the sample. Connect the VTF into the dedicated socket on the back of the AM4 Digital PRO. Select the operating temperature required on the VTF thermoregulator. Turn the temperature control knob on the frontal panel to maximum value. The heating magnetic stirrer always has primary control on the heating plate temperature. With VTF thermoregulator always select the max temperature on the Heating magnetic Stirrer. The temperature control function of the heating plate can also be used as a safety thermostat. In this case the maximum temperature of the heating plate will not exceed the temperature set on the Heating magnetic Stirrer meaning that a longer heating time is required in order to reach the VTF thermoregulator temperature setting.

3.3 Using the external probe

Screw the threaded support rod (optional accessory, code A00001069) into its seat on the back of the instrument and fasten the clamp (optional accessory, code A00000004) onto the support rod. Fit the temperature probe A00000268 into the clamp making sure that it is completely immersed in the sample.

Plug the temperature probe A00000268 into the dedicated socket on the back of AM4 Digital PRO with the instrument off. A led on the left-hand side of the display indicates that the probe is connected, the display indicates the probe temperature reading. Select the operating temperature using the heating knob on the AM4 Digital PRO. A led at the center of the display indicates that the set-point temperature is shown; after 5 seconds the led goes out and the display indicates the probe temperature.

3.4 Setting mode

Turn the left knob twice from the minimum to the maximum during the visualization of the software version in order to enter into the setting mode.

Once there, by turning the left knob, it is possible to pass from one of the following parameter to the next ones:

Parameter shown		Default value	Range	Description
CAL	P100	0,0°C	-10,0 ÷ +10,0°C	External probe calibration
CAL	tA	0,0°C	-10,0 ÷ +10,0°C	Probe calibration at 20°C
FS	tCJ	370°C	0 ÷ 370°C	Maximum temperature settable
OPt	tYpE	A	A – b	(*) Restarting mode
OPt	P100	FASt	FInE - FASt	(**) Type of control with external probe
ALL	P100	POn	POn-POFF	(***) External probe safety

Turn the right knob from the minimum to the maximum and then back to the minimum in order to pass to the parameter you want to modify. The display will show the value saved for 3 seconds, and then the new value according to the position of the right knob.

The new value is automatically saved after 4 seconds, if during this time the right knob is not moved (no change of value). "dOnE" is shown and then the system goes back to the main menu.

Once the new values are set, switch off and restart the instrument.

(*) When **A** is selected: in case of black-out, once the power is back the device will remain off

When **b** is selected: in case of black-out, once the power is back the device will automatically restart with the latest settings

(**) Select **FInE** to minimize overshoot and fluctuations of temperature.

Select **FASt** to save time.

(***) When **POn** is selected: if the instrument don't measure an increase of temperature, a signal of alarm is shown.

3.5 Error messages

If an error message appears on the display, please contact VELP Scientifica's technical service department.

When the display shows an error message, the stirring and heating functions stop automatically.

NOTE: To remove the error message disconnect the instrument from the power supply.

AL1 Overtemperature (temperature is higher than 430 °C)

AL2 Excessive heating time

AL3 The stirring system doesn't run correctly

AL4 Overtemperature of the external probe (temperature is higher than 260 °C) - only when A00000268 is connected

AL5 Overtemperature of the safety probe (temperature is higher than 430 °C)

AL6 Probe not inserted inside the sample

4. Maintenance

No routine or extraordinary maintenance is necessary apart from periodically cleaning the unit as described in this manual. In compliance with the product guarantee law, replace the fuses and repairs to our units must be carried out in our factory, unless previously agreed otherwise with local distributors and using Velp spare part only.

The instrument must be transported in its original packaging and any indications present on the original packaging must be followed (e.g. palletized).

It is the responsibility of the user, to properly decontaminate the unit in case of hazardous substances remaining on the surface or interior of the device. If in doubt about the compatibility of a cleaning or decontamination product, contact the manufacturer or distributor.

4.1 Cleaning

Disconnect the unit from the power supply and use a cloth dampened with an non-inflammable non-aggressive detergent.

5. Technical data

Power supply	230V / 50-60Hz (+/-10%)
Dimensions WxHxD	715x115x246 (28.1x4.5x9.7)
Weight	9.0 Kg
Power input	630W x 4 (2520W)
Heat output	600W x 4 (2400W)
Motor rating output	10 W x 4
Continuous operating mode	Admitted
Diameter of the heating plate	135mm x 4
Programmable temperature range	Room temp. – 370 °C
Type of temperature control	Digital
Temperature selection	5 °C steps
Reading temperature resolution	1 °C steps
Overtemperature protection	Yes
Stirring capacity	20 liters of H ₂ O x 4
Programmable speed range	50 – 1500 rpm
Type of motor control	Digital
Speed selection	10 rpm steps
External connection	VTF
Electric rating VTF	5 Vdc (max)
Accuracy of temperature with VTF	+/-0,5°C*
Accuracy of temperature with external probe	+/-1°C*
Noise level	<< 80 dBa
Temperature range	+5...+40 °C
Storage temperature range	-10...+60 °C
Max humidity	80%
Overvoltage category	II
Pollution degree CEI EN61010-1	2
Max altitude	4000 m

* in following conditions: 800ml water in 1liter glass beaker (diameter 105mm), stirring bar 8x40mm, 600rpm, 50°C.

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
We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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EN 61010-2-051

EN 61010-2-010

and satisfies the essential requirements of the following directives:

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10. Declaration of conformity

We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following regulations:

S.I. 2016/1101 Electrical Equipment (Safety) Regulations 2016

S.I. 2016/1091 Electromagnetic Compatibility Regulations 2016

according to the relevant designated standards:

EN 61010-1 Safety requirements for electrical equipment for measurement, control, and laboratory use
Part 1: General requirements

EN 61010-2-010 Particular requirements for laboratory equipment for the heating of material

EN 61010-2-051 Particular requirements for laboratory equipment for mixing and stirring

EN 61326-1 Electrical equipment for measurement, control and laboratory use - EMC requirements -
Part 1: General requirements

and satisfies the essential requirements of regulations:

S.I. 2008/1597 Supply of Machinery (Safety) Regulations 2008

S.I. 2012/3032 Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic
Equipment Regulations 2012

S.I. 2013/3113 Waste Electrical and Electronic Equipment Regulations 2013

plus modifications.

Thank you for having chosen VELP!

Established in 1983, VELP is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. VELP works according to **ISO 9001**, **ISO14001** and **OHSAS 18001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

Analytical instruments

Elemental Analyzers
Digestion Units
Distillation Units
Solvent Extractors
Fiber Analyzers
Dietary Fiber Analyzers
Oxidation Stability Reactor
Consumables

Laboratory Equipment

Magnetic Stirrers
Heating Magnetic Stirrers
Heating Plates
Overhead stirrers
Vortex mixers
Homogenizers
COD Thermoreactors
BOD and Respirometers
Cooled Incubators
Flocculators
Overhead Shakers
Turbidimeter
Radiation Detector
Open Circulating Baths
Pumps

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