**Instruction Manual**  
**Manuale di istruzioni**  
**Manuel d'instructions**  
**Manual de instrucciones**  
**Bedienungsanleitung**

**ARE Heating Magnetic Stirrer**

F20500162, F20510162

**AREX Heating Magnetic Stirrer**

F20500413, F20510413

**General Information / Informazioni Generali / Informations Générales / Información General /  
Allgemeine Hinweise**

Before using the unit, please read the following instruction manual carefully.

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Do not dispose of this equipment as urban waste, in accordance with EEC directive 2002/96/CE.

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**This unit must be used for laboratory applications only.**

The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions.

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**This unit has been designed and manufactured in compliance with the following standards:**

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**IEC/EN 61010-1  
IEC/EN 61010-2-051**

Electrical equipment for laboratory use

**UL 61010-1**

General requirement - Canadian electrical code

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Do not use with explosive and dangerous materials for which the equipment is not designed. The stirrer must not be used in explosive atmospheres, in bain-marie and to stir combustible liquids that have a low combustion temperature. The minimum fire point of flammable solution is 750 °C. Only small amounts (< 50 ml) of flammable liquid can be used with the device.

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The unit is fitted with two fuses (2xT5 A L 250 V (for 230V), 2xF8 A 250 V (for 115 V)), found in the socket on the back. To replace one or more disconnect the mains cable and, using a screwdriver, lift up the small cover on the fuse box.

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It is responsibility of the user appropriately decontaminate the instrument in case of dangerous substances fall on or in it. It is also responsibility of the user to use safety substances for cleaning or decontaminating, which do not react with internal parts of the instrument or with the material contained in it. In case of doubts on the compatibility of a cleaning solution, contact the manufacturer or local distributor.

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## Contents / Indice / Index / Índice / Inhalt

1.	INTRODUCTION.....	5
2.	ASSEMBLY AND INSTALLATION.....	5
	2.1 ELECTRICAL CONNECTIONS.....	5
	2.2 START-UP.....	5
3.	OPERATING CONTROLS .....	5
	3.1 USING THE THERMOREGULATOR VTF.....	5
4.	MAINTENANCE .....	6
	4.1 CLEANING.....	6
5.	TECHNICAL DATA .....	6
6.	ACCESSORIES / SPARE PARTS .....	7

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The application of new technology has led to the creation of these modern magnetic stirrers for laboratory use whose basic features are safety, reproducibility of results, high performance and energy saving.

The structure of the unit is made out of pressure die-cast aluminum treated with epoxy resins offering a high resistance to the many chemical aggressions typical of the laboratory environment. The meticulously studied design of the pressure die-cast casing combines aesthetic fluidity with ergonomics. A wide central recess on both sides offers a good grip and facilitates bench-top handling whilst the handle moulded directly into the rear of the die-cast structure can be used to transport the instrument once it has been disconnected from the power supply.

The pressure die-cast structure is designed so that accidental liquid spills cannot reach the internal parts of the unit. The electrical connections are recessed into the rear of the instrument and offer optimum electrical protection of the external electrics in the case of liquid spills.

The engineering of the unit facilitates internal inspection if technical assistance is required.

The heating plate is made out of aluminum alloy with special coating and is specially treated to ensure:

- Optimum heat distribution and a high specific power thanks to the circular configuration
- Optimum temperature homogeneity across the plate
- High resistance to thermal stress and thermal shock, chemical aggressions, scratches and surface abrasions.

Magnetic stirring is generated by the VELP Scientifica Patented Composite Magnet driven by an asynchronous mono-phase brushless motor in alternating current which offers a virtually unlimited duration.

**NOTE:** using the heating plate at high temperatures may cause discoloring. This does not alter the thermal, mechanical and chemical resistance of the plate in any way.

## 2. Assembly and installation

Check the integrity of the unit after unpacking. The box includes:

- Heating Magnetic Stirrer
- Power supply cord
- Instruction manual

### 2.1 Electrical connections

After having unpacked the instrument, place the unit on the laboratory bench.

Before connecting the instrument to the power supply, make sure that the values on the rating plate correspond to those of the power supply. Connect the unit to the power supply using the cord supplied.

Ensure that the socket and the relative cut-off device conform to current safety norms and are easy to reach.

### 2.2 Start-up

Rotate the speed (right) and temperature (left) knobs completely to the left. Place the flask containing the sample and a suitable magnetic stirring bar on the stirring plate. Then, set the speed and temperature by turning the dedicated knobs.

## 3. Operating controls

### REGULATION KNOBS

The right knob on the front of the unit can be used for quick precision regulation of mixing speed from 100 to 1500 rpm. The left one is for temperature, up to 370 °C.

### ON-OFF SWITCH

The on-off switch turns the unit on and off. If the switch is in the "OFF" position the unit is off; if the switch is in the "ON" position the unit is on.

Always turn the unit off after use.

### 3.1 Using the thermoregulator VTF (Only for models AREX)

Screw the threaded support rod (optional accessory, A00001069) into its seat on the back of the instrument and fasten the VTF thermoregulator onto the support rod. Place the temperature probe in the receptacle making sure that it is completely immersed in the sample. Connect the two instruments (AREX and VTF) by plugging the VTF into the dedicated socket on the back of the AREX.

Select the operating temperature required on the VTF thermoregulator. Turn the temperature control knob on the front panel of the AREX to maximum.

The heating magnetic stirrer always has primary control of the heating plate temperature.

When using the VTF thermoregulator always select the max temperature on the AREX. The temperature control function of the heating plate can also be used as a safety thermostat. In this case the maximum temperature of the heating plate will not exceed the temperature setting on the AREX meaning that a longer heating time is required in order to reach the VTF thermoregulator temperature setting.

### 3.2 Error messages

The unit is fitted with safety devices which cut-off the power supply to the heating plate in the case of malfunctions. The led display on the front panel indicates the type of malfunction:

#### Display

Flashing Stirrer LED (once every second)  
Flashing Heating LED (once every second)  
Flashing Heating LED (2 times every second)  
Flashing Heating LED (once every 3 seconds)

#### Cause

Fault in the stirring system  
Overheated heating plate or thermocouple open  
Fault in the temperature reading circuit  
Thermocouple wires reversed or white wire of potentiometer board not connected

Should any of the above occur, please contact your nearest VELP Scientifica service centre.

## 4. Maintenance

No routine or extraordinary maintenance is necessary apart from periodically cleaning the unit as described in this manual. In compliance with the product guarantee law, repairs to our units must be carried out in our factory, unless previously agreed otherwise with local distributors.

The instrument must be transported in its original packaging and any indications present on the original packaging must be followed (e.g. palletized).

It is the responsibility of the user, to properly decontaminate the unit in case of hazardous substances remaining on the surface or interior of the device. If in doubt about the compatibility of a cleaning or decontamination product, contact the manufacturer or distributor.

### 4.1 Cleaning

Disconnect the unit from the power supply and use a cloth dampened with a non-flammable non-aggressive detergent.

## 5. Technical data

Power supply	230 V / 50-60 Hz o 115V / 60 Hz (+/-10%)
Dimensions WxHxD	165x115x280 mm
Weight	2.5 kg
Overall power	630 W
Power of the heating plate	600 W
Diameter of the heating plate	135mm
Programmable temperature range	50 – 370 °C
Type of temperature control	Analog
Overtemperature protection	Yes
Motor rating output	10 W
Stirring capacity	15 liters of H <sub>2</sub> O (ARE) - 20 liters of H <sub>2</sub> O (AREX)
Programmable speed range	100 – 1500 rpm
Counter-reaction	Constant speed even when the viscosity changes
Type of motor control	Electronic
Temperature range	+5...+40 °C
Storage temperature range	-10...+60 °C
Max humidity	80%
Level of electrical protection CEI EN60529	IP 42
Pollution degree CEI EN61010-1	2
Overvoltage category	II
Noise level	<< 80 dBa
Max altitude	4000 m

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
We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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EN61010-2-051(2015)                      EN 61010-2-010 (2014)

and satisfies the essential requirements of the following directives:

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- Machinery directive 2006/42/EC / Macchine 2006/42/CE / Machines 2006/42/CE / Máquinas 2006/42/CE / Maschinen 2006/42/EG

- Low voltage directive 2014/35/EU / Bassa tensione 2014/35/EU / Basse tension 2014/35/EU / Baja tensión 2014/35/EU / Niederspannung 2014/35/EU

- Electromagnetic compatibility directive 2014/30/EU / Compatibilità elettromagnetica 2014/30/EU / Compatibilité électromagnétique 2014/30/EU / Compatibilidad electromagnética 2014/30/EU / Elektromagnetische Verträglichkeit 2014/30/EU

- plus modifications / più modifiche / plus modifications / más sucesivas modificaciones / in der jeweils gültigen Fassung.

## Thank you for having chosen VELP!

Established in 1983, VELP is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. VELP works according to **ISO 9001**, **ISO14001** and **OHSAS 18001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

### Analytical instruments

Elemental Analyzers  
Digestion Units  
Distillation Units  
Solvent Extractors  
Fiber Analyzers  
Dietary Fiber Analyzers  
Oxidation Stability Reactor  
Consumables

### Laboratory Equipment

Magnetic Stirrers  
Heating Magnetic Stirrers  
Heating Plates  
Overhead stirrers  
Vortex mixers  
Homogenizers  
COD Thermoreactors  
BOD and Respirometers  
Cooled Incubators  
Flocculators  
Overhead Shakers  
Turbidimeter  
Radiation Detector  
Open Circulating Baths  
Pumps

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VELP Scientifica Srl