**Instruction Manual**  
**Manuale di istruzioni**  
**Manuel d'instructions**  
**Manual de instrucciones**  
**Bedienungsanleitung**

## AREX DIGITAL PRO Heating Magnetic Stirrer

F20500410, F20510410

### General Information / Informazioni Generali / Informations Générales / Información General / Allgemeine Hinweise

Before using the unit, please read the following instruction manual carefully.

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Do not dispose of this equipment as urban waste, in accordance with EEC directive 2002/96/CE.

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**This unit must be used for laboratory applications only.** The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions. If the product is used in a not specified way by the manufacturer or with non specified accessories, product's safety may be compromised.

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**This unit has been designed and manufactured in compliance with the following standards:**

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**IEC/EN 61010-1**  
**IEC/EN 61010-2-051**

Electrical equipment for laboratory use

**UL 61010-1**

General requirement - Canadian electrical code

**CAN/CSA-C22.2 No.61010-1**

VELP reserves the right to modify the characteristics of its products with the aim to constantly improving their quality.  
Nell'impegno di migliorare costantemente la qualità dei prodotti, VELP si riserva la facoltà di variarne le caratteristiche.  
Dans le but d'améliorer constamment la qualité de ses produits, VELP se réserve le droit d'apporter des modifications aux caractéristiques de ceux-ci.  
VELP se reserva el derecho de modificar las características de productos con el fin de mejorar constantemente su calidad.  
VELP behält sich zum Zwecke der ständigen Verbesserung der Produktqualität das Recht auf Änderung der Geräteeigenschaften vor.

### **Safety Regulations / Norme di Sicurezza / Consignes de Sécurité / Advertencias de Seguridad / Sicherheitshinweise**

The plug disconnects the instrument. Therefore, place the instrument where it can be quickly disconnected. / La spina è il mezzo di disconnessione dell'apparecchio. Pertanto, non posizionare l'apparecchio in modo che sia difficile azionare il mezzo di disconnessione. / Le bouchon est le moyen de déconnexion de l'appareil. Par conséquent, placer l'appareil où il peut être rapidement débranché. / El tapón es el medio de desconexión del dispositivo. No coloque el dispositivo en una forma que es difícil de desconectar. / Der Stecker trennt das Gerät. Daher Stellen Sie das Instrument, wo es schnell getrennt werden kann.

Hotplate temperature: up to 370 °C. / Temperatura piastra riscaldante: fino a 370 °C. / Température de la plaque chauffante: jusqu'à 370 °C. / Temperatura de la placa calefactora: hasta 370 °C. / Temperaturbereich Heizplatte: bis zu 370 °C.

The heated solution may release toxic, dangerous or poisonous gases. Adequate safety measures must be taken, in accordance with the safety regulations in force, including the presence of hood and personal protective equipment (masks, gloves, goggles, etc.).

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Do not use with explosive and dangerous materials for which the equipment is not designed. The stirrer must not be used in explosive atmospheres, in bain-marie and to stir combustible liquids that have a low combustion temperature. The minimum fire point of flammable solution is 750 °C. Only small amounts (< 50 ml) of flammable liquid can be used with the device.

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The unit is fitted with two fuses (2xT5 A L 250 V (for 230V), 2xF8 A 250 V (for 115 V)), found in the socket on the back. To replace one or more disconnect the mains cable and, using a screwdriver, lift up the small cover on the fuse box.

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It is responsibility of the user appropriately decontaminate the instrument in case of dangerous substances fall on or in it.

It is also responsibility of the user to use safety substances for cleaning or decontaminating, which do not react with internal parts of the instrument or with the material contained in it. In case of doubts on the compatibility of a cleaning solution, contact the manufacturer or local distributor.

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The AREX Digital PRO heating magnetic stirrer with aluminium alloy plate, is used in the laboratory for general use and for all those applications that require a precise regulation of the stirring speed and of the heating plate temperature.

The heating plate (1) is made out of aluminum alloy with ceramic coating and is specially treated to ensure:

- Optimum heat distribution and a high specific power thanks to the circular configuration
- Optimum temperature homogeneity across the plate
- High resistance to thermal stress and thermal shock, chemical aggressions, scratches and surface abrasions.

The AREX Digital PRO is used to mix liquids using a magnetically driven stir bar placed inside the vessel. The size and shape of the magnetic stir bar determines the stirring efficiency at any given speed. Magnetic stirring is generated by the VELP Scientifica Patented Composite Magnet driven by an asynchronous mono-phase brushless motor in alternating current which offers a virtually unlimited duration.

**NOTE:** It is therefore advisable to choose the most suitable magnetic stir bar in relation to the quantity and characteristics of the liquid and the type of vessel used.

The stirring bar which satisfies most applications is code A00000356 (8x40 mm).

The instrument is turned on using the ON/OFF switch (2). The stirring speed (3) and the temperature (4) of the heating plate are adjusted using the relative knobs on the front panel, the values are indicated on the digital display (5).

**NOTE:** If the heating function is used the vessel must be made of a suitable material to withstand the foreseen temperature.

**NOTE:** using the heating plate at high temperatures may cause discoloring. This does not alter the thermal, mechanical and chemical resistance of the plate in any way

Fig. 1

## 2. Assembly and installation

Check the integrity of the unit after unpacking. The box includes:

- AREX Digital PRO Hot Plate Stirrer
- Power supply cord
- Instruction manual

### 2.1 Electrical connections

After having unpacked the instrument, place the unit on the laboratory bench.

Before connecting the instrument to the power supply, make sure that the values on the rating plate correspond to those of the power supply. The equipment shall be connected to a power socket with protective earth connection, using only the power supply cord provided with the instrument.

Ensure that the socket and the relative cut-off device conform to current safety norms and are easy to reach.

**NOTE:** the mains cable must remain far away from the hot plate.

### 2.2 Start-up

Rotate the speed (right) and temperature (left) knobs to the OFF position. Place the flask containing the sample and a suitable magnetic stirring bar on the stirring plate. Then, set the speed and temperature by turning the dedicated knobs.

## 3. Operating controls

### 3.1 Turning on and regulating

Connect the unit to mains and turn it on using the on-off switch. The green led indicates that the instrument is On. In the next 5 seconds, displays shows the software version and the restarting mode (A or b, chapter 3.4 "Setting mode").

Start the stirring function by turning the stirrer knob "Stirrer rpm" on the front panel. Speeds of from 50 to 1500 rpm (10 rpm resolution) can be selected. The speed shown flashes until it reaches the set-point. A microprocessor ensures constant speed even when the viscosity changes (counter-reaction).

To start the heating function turn the heating knob "Heating °C" on the front panel. Temperatures from 5 to 370 °C can be selected at 5 °C intervals. The display shows the real temperature but every 5 seconds the temperature set is shown for 1 second, indicated by the central led on the display "Set" turned on.

**NOTE:** in case of black-out, once the power is back the device will restart in the restarting mode set (Opt tYpE, chapter 3.4).

When heating is turned Off (Heating knob OFF) and the heating plate temperature exceeds 50 °C, the message “Hot” flashes on the display until the temperature of the heating plate falls to below 50 °C (or until the instrument is turned off using the ON/OFF switch).

**NOTE:** this warning is not active if the instrument is not powered.

### 3.2 Using the thermoregulator VTF

Screw the threaded support rod into its seat on the back of the instrument and fasten the VTF thermoregulator onto the support rod. Place the temperature probe in the flask making sure that it is completely immersed in the sample. Connect the two instruments (AREX Digital PRO and VTF) by plugging the VTF into the dedicated socket on the back of the AREX Digital PRO.

Select the operating temperature required on the VTF thermoregulator.

Turn the temperature control knob on the front panel of the AREX Digital PRO to the desired temperature.

**NOTE:** when the two instruments are connected, the display of the AREX DIGITAL PRO does not show the temperature but a dashed line (---). The operating temperature is shown on the display of the VTF thermoregulator.

The AREX DIGITAL PRO always has primary control of the heating plate temperature. When using the VTF thermoregulator always select the max temperature on the AREX Digital PRO. The temperature control function of the heating plate can also be used as a safety thermostat. In this case the maximum temperature of the heating plate will not exceed the temperature setting on the AREX Digital PRO meaning that a longer heating time is required in order to reach the VTF thermoregulator temperature setting. (chap. 3.4)

### 3.3 Using the external probe

Screw the threaded support rod (optional accessory, code A00001069) into its seat on the back of the instrument and fasten the clamp (optional accessory, code A00000004) onto the support rod. Fit the temperature probe A00000268 into the clamp making sure that it is completely immersed in the sample.

Plug the temperature probe A00000268 into the dedicated socket on the back of the AREX Digital PRO with the instrument off. A led on the left-hand side of the display indicates that the probe is connected, the display indicates the probe temperature reading. Select the operating temperature using the heating knob on the AREX Digital PRO. A led at the center of the display indicates that the set-point temperature is shown; after 5 seconds the led goes out and the display indicates the probe temperature.

### 3.4 Setting mode

Turn the left knob twice from the minimum to the maximum during the visualization of the software version in order to enter into the setting mode.

Once there, by turning the left knob, it is possible to pass from one of the following parameter to the next ones:

Parameter shown	Default value	Range	Description
CAL P100	0,0°C	-10,0 ÷ +10,0°C	External probe calibration
CAL tA	0,0°C	-10,0 ÷ +10,0°C	Probe calibration at 20°C
FS tCJ	370°C	0 ÷ 370°C	Maximum temperature settable
OPt tYpE	A	A – b	(*) Restarting mode
OPt P100	FASt	FlnE - FASt	(**) Type of control with external probe
ALL P100	POn	POn-POFF	(***) External probe safety

Turn the right knob from the minimum to the maximum and then back to the minimum in order to pass to the parameter you want to modify. The display will show the value saved for 3 seconds, and then the new value according to the position of the right knob.

The new value is automatically saved after 4 seconds, if during this time the right knob is not moved (no change of value). “dOnE” is shown and then the system goes back to the main menu.

Once the new values are set, switch off and restart the instrument.

(\*) When **A** is selected: in case of black-out, once the power is back the device will remain off

When **b** is selected: in case of black-out, once the power is back the device will automatically restart with the latest settings

(\*\*) Select **FlnE** to minimize overshoot and fluctuations of temperature.

Select **FASt** to save time.

(\*\*\*) When **POn** is selected: if the instrument don't measure an increase of temperature, a signal of alarm is shown.

### 3.5 Error messages

If an error message appears on the display, please contact VELP Scientifica's technical service department. When the display shows an error message, the stirring and heating functions stop automatically.

**NOTE:** To remove the error message disconnect the instrument from the power supply.

**AL1** Overtemperature (temperature is higher than 430 °C)

**AL2** Excessive heating time

**AL3** The stirring system doesn't run correctly

**AL4** Overtemperature of the external probe (temperature is higher than 260 °C) - only when A00000268 is connected

**AL5** Overtemperature of the safety probe (temperature is higher than 430 °C)

**AL6** Probe not inserted inside the sample

## 4. Maintenance

No routine or extraordinary maintenance is necessary apart from periodically cleaning the unit as described in this manual. In compliance with the product guarantee law, repairs to our units must be carried out in our factory, unless previously agreed otherwise with local distributors.

The instrument must be transported in its original packaging and any indications present on the original packaging must be followed (e.g. palletized).

It is the responsibility of the user, to properly decontaminate the unit in case of hazardous substances remaining on the surface or interior of the device. If in doubt about the compatibility of a cleaning or decontamination product, contact the manufacturer or distributor.

### 4.1 Cleaning

Disconnect the unit from the power supply and use a cloth dampened with an non-inflammable non-aggressive detergent.

## 5. Technical data

Power supply	230 V / 50-60 Hz or 115V / 60 Hz (+/-10%)
Dimensions WxHxD	165x115x280 mm (6.5x4.5x11 in)
Weight	2.6 Kg (5.7 lb)
Power input	630 W
Heat output	600 W
Motor rating output	10 W
Continuous operating mode	Admitted
Dimension of the heating plate	Ø 135 mm
Programmable temperature range	Room temp. - 370 °C
Type of temperature control	Digital
Temperature selection	5 °C steps
Reading temperature resolution	1 °C steps
Overtemperature protection	Yes
Stirring capacity	20 liters of H <sub>2</sub> O
Programmable speed range	50 – 1500 rpm
Type of motor control	Digital
Speed selection	10 rpm steps
Connection	VTF - External probe Pt100
Electric rating for external probe	5 Vdc – 1W (max)
Accuracy of temperature with VTF	+/-0,5°C*
Accuracy of temperature with external probe	+/-1°C*
Noise level	<< 80 dBa
Temperature range	+5...+40 °C
Storage temperature range	-10...+60 °C
Max humidity	80%
Level of electrical protection CEI EN60529	IP 42
Overvoltage category	II
Pollution degree CEI EN61010-1	2
Max altitude	4000 m

\* in following conditions: 800ml water in 1 liter glass beaker (diameter 105mm), stirring bar 8x40mm, 600rpm, 50°C.

## 6. Accessories / Spare parts

A00001056	Magnetic stir bar Ø. 6x35mm	A00000340	MonoAluBlock, 17 pos. Ø28 x h 43 mm
A00001057	Magnetic stir bar Ø. 6x20mm	A00000339	MonoAluBlock, 17 pos. Ø28 x h 30 mm
A00001061	Magnetic stir bar Ø 10x60mm	A00000338	MonoAluBlock, 17 pos. Ø28 x h 24 mm
A00001062	Magnetic stir bar, Ø3x6 mm	A00000337	MultiAluBlockTM, 11 pos. Ø12 x h 14 mm
A00001063	Magnetic stir bar, Ø4,5X12 mm	A00000336	Magnetic cross shape stir bar, Ø10x5 mm
A00001069	Support rod	A00000334	Hemispheric bowl 1000ml flasks, plate 135
A00000352	Magnetic cross shape stir bar, Ø20x8 mm	A00000333	Hemispheric bowl 500ml flasks, plate 135
A00000354	Magnetic disc stir bar, Ø10x6 mm	A00000332	Hemispheric bowl 250ml flasks, plate 135
A00000356	Magnetic stir bar, Ø8x40 mm	A00000331	Hemispheric bowl 100ml flasks, plate 135
A00000351	Handle for AluBlock removal	A00000330	Hemispheric bowl 50ml flasks, plate 135
A00000382	Extension for support rod ARE – ARES	A00000327	MultiAluBlockTM, 4 pos. Ø21 x h31 mm
A00000375	Hemispheric bowl 5L flasks, plate 135	A00000325	MultiAluBlockTM, 4 pos. Ø28 x h 30 mm
A00000374	Hemispheric bowl 3L flasks, plate 135	A00000324	MultiAluBlockTM, 4 pos. Ø28 x h 43 mm
A00000373	Hemispheric bowl 25ml flasks, plate 135	A00000284	Protective cover
A00000346	PTFE Safety cover for bowl 1000 ml	A00000280	Clamp with probe support
A00000345	PTFE Safety cover for bowl 500 ml	A00000268	External probe
A00000344	PTFE Safety cover for bowl 250 ml		
A00000343	PTFE Safety cover for bowl 100 ml	10005213	Knob 24D blue
A00000342	PTFE Safety cover for bowl 50 ml	10000239	Bumpon 13Dx5H embedded
A00000341	MonoAluBlock, 40 pos. Ø12 x h 14 mm		

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
We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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EN61010-2-051(2015)

EN 61010-2-010 (2014)

and satisfies the essential requirements of the following directives:

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- Machinery directive 2006/42/EC / Macchine 2006/42/CE / Machines 2006/42/CE / Máquinas 2006/42/CE / Maschinen 2006/42/EG

- Low voltage directive 2014/35/EU / Bassa tensione 2014/35/EU / Basse tension 2014/35/EU / Baja tensión 2014/35/EU / Niederspannung 2014/35/EU

- Electromagnetic compatibility directive 2014/30/EU / Compatibilità elettromagnetica 2014/30/EU / Compatibilité électromagnétique 2014/30/EU / Compatibilidad electromagnética 2014/30/EU / Elektromagnetische Verträglichkeit 2014/30/EU

- plus modifications / più modifiche / plus modifications / más sucesivas modificaciones / in der jeweils gültigen Fassung.

## Thank you for having chosen VELP!

Established in 1983, VELP is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. VELP works according to **ISO 9001**, **ISO14001** and **OHSAS 18001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

### Analytical instruments

Elemental Analyzers  
Digestion Units  
Distillation Units  
Solvent Extractors  
Fiber Analyzers  
Dietary Fiber Analyzers  
Oxidation Stability Reactor  
Consumables

### Laboratory Equipment

Magnetic Stirrers  
Heating Magnetic Stirrers  
Heating Plates  
Overhead stirrers  
Vortex mixers  
Homogenizers  
COD Thermoreactors  
BOD and Respirometers  
Cooled Incubators  
Flocculators  
Overhead Shakers  
Turbidimeter  
Radiation Detector  
Open Circulating Baths  
Pumps

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VELP Scientifica Srl