Instruction Manual
Manuale di istruzioni
Manuel d'instructions
Manual de instrucciones
Bedienungsanleitung
指导手册

BOD Sensor

F102B0133

General Information / Informazioni Generali / Informations Générales / Información General / Allgemeine Hinweise / 一般信息

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This unit has been designed and manufactured in compliance with the following standards:

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Safety requirements for electrical equipment for measurement, control and for laboratory use

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We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

Noi, casa costruttrice VELP SCIENTIFICA, dichiariamo sotto la ns. responsabilità che il prodotto è conforme alle seguenti norme:
Nous, VELP Scientifica, déclarons sous notre responsabilité que le produit est conforme aux normes suivantes:
Nosotros casa fabricante, VELP Scientifica, declaramos bajo nuestra responsabilidad que el producto es conforme con las siguientes normas:
Der Hersteller, VELP Scientifica, erklärt unter eigener Verantwortung, dass das Gerät mit folgenden Normen übereinstimmt:
我们，制造商 VELP Scientifica，根据我们的责任声明，该产品的生产符合以下标准。

EN 61010-1 EN 61326-1 2015/863/EU (RoHS III) 2012/19/EU (WEEE)

and satisfies the essential requirements of the following directives:

e soddisfa i requisiti essenziali delle direttive:
et qu'il satisfait les exigences essentielles des directives:
y cumple con los requisitos esenciales de las directivas:
und den Anforderungen folgender Richtlinien entspricht:
并满足以下指令的基本要求。

- Machinery directive 2006/42/EC / Macchine 2006/42/CE / Machines 2006/42/CE / Máquinas 2006/42/CE / Maschinen 2006/42/EG / 机械指令 2006/42/EC

- Electromagnetic compatibility directive 2004/108/EC / Compatibilità elettromagnetica 2004/108/CE / Compatibilité électromagnétique 2004/108/CE / Compatibilidad electromagnética 2004/108/CE / Elektromagnetische Verträglichkeit 2004/108/EG / 电磁兼容性指令 2004/108/EC

- plus modifications / più modifiche / plus modifications / más sucesivas modificaciones / in der jeweils gültigen Fassung / 加上变化

Declaration of conformity

We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following regulations:

S.I. 2016/1101 Electrical Equipment (Safety) Regulations 2016
S.I. 2016/1091 Electromagnetic Compatibility Regulations 2016

according to the relevant designated standards:

EN 61010-1 Safety requirements for electrical equipment for measurement, control, and laboratory use
Part 1: General requirements
EN 61326-1 Electrical equipment for measurement, control and laboratory use - EMC requirements -
Part 1: General requirements

and satisfies the essential requirements of regulations:

S.I. 2008/1597 Supply of Machinery (Safety) Regulations 2008
S.I. 2012/3032 Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment Regulations 2012
S.I. 2013/3113 Waste Electrical and Electronic Equipment Regulations 2013

plus modifications.

Analysis requirements

BOD Sensor or RESPIROMETRIC Sensor
Potassium hydroxide KOH
Stirring bar
Incubator (20 °C ± 0.5)

Alkali container
500 ml dark glass bottle
Stirring base

Analysis procedures

1. Pour into bottles the volume of sample according to adopted scale, measured by a graduated cylinder.
2. Put into each bottle a magnetic stirring bar.
3. Fill the alkali holder, small container located under bottle cap, with an amount of carbon dioxide absorber, enough to avoid losses through the holes. If some alkali falls into a bottle, wash thoroughly before pouring again a volume of sample under examination.
4. Place the bottles into position in the stirring equipment.
5. Introduce the stirring equipment into the refrigerated thermostat set to the temperature chosen for incubation. Connect magnetic stirrers to mains introducing the plug into the socket. Put into operation the refrigerated thermostat.
6. After 30-40 minutes the apparatus and samples are usually in thermal equilibrium at the chosen temperature.
7. The apparatus is ready to start BOD measurement. Place a BOD Sensor on each bottle and screw tight. Reset each BOD Sensor cancelling any memorized value, choose the most appropriate scale and start the measuring cycle.

Analytical method

Biochemical Oxygen Demand is expressed as weight of oxygen consumed per unit volume of water during a defined period of time at a defined temperature, e.g. mg O₂/l (or ppm = parts per million) during 5 days at 20 °C.

For carbon dioxide absorption:

Potassium hydroxide (KOH) scales, commercial grade **or** not deliquescent soda lime, 1.0-1.7 mm granules

If samples are to be neutralized before incubation:

Sulfuric acid (H₂SO₄) 1.0 N solution

Sodium hydroxide (NaOH) 1.0 N solution

If samples are to be supplemented with nutrients:

Distilled water, high quality, free of copper and organic matter

Calcium chloride, anhydrous (CaCl₂)

Potassium ortophosphate, biacid, (KH₂PO₄)

Sodium ortophosphate, monoacid, (Na₂HPO₄)

Iron chloride, ferric, exahydrate (FeCl₃ • 6 H₂O)

Magnesium sulfate, heptahydrate (MgSO₄ • 7 H₂O)

Potassium ortophosphate, monoacid, (K₂HPO₄)

Ammonium chloride (NH₄Cl)

If samples contain chlorine:

Sodium sulfite (Na₂SO₃)

If nitrification in samples is to be inhibited:

Allyl thiourea (Thiosinamine) C₄H₈N₂S **or** 2-chloro-6 (thrichloromethyl) pyridine, C₆H₃Cl₄N (CTCMP, N-Serve)

Additive for mineral fertilizers aimed to lower nitrogen losses in soils.

Contents / Indice / Index / Índice / Inhalt

DECLARATION OF CONFORMITY / DICHIARAZIONE DI CONFORMITÀ / DECLARATION DE CONFORMITE / DECLARACIÓN DE CONFORMIDAD / KONFORMITÄTSERKLÄRUNG / 符合性声明 	2
DECLARATION OF CONFORMITY 	2
ANALYSIS REQUIREMENTS	3
ANALYSIS PROCEDURES	3
ANALYTICAL METHOD	3
1. INTRODUCTION	5
2. ASSEMBLY AND INSTALLATION	5
3. OPERATING CONTROLS	5
3.1 DESCRIPTION OF KEYS	5
3.2 SCALE SELECTION	5
3.3 ALARMS	5
4. MAINTENANCE	5
5. TECHNICAL DATA	5
6. ACCESSORIES / SPARE PARTS	5
1. INTRODUZIONE	6
2. MONTAGGIO ED INSTALLAZIONE	6
3. CONTROLLI DI FUNZIONAMENTO	6
3.1 DESCRIZIONE DEI TASTI	6
3.2 IMPOSTAZIONE DELLA SCALA	6
3.3 ALLARMI	6
4. MANUTENZIONE	6
5. CARATTERISTICHE TECNICHE	6
6. ACCESSORI / PARTI DI RICAMBIO	6
1. INTRODUCTION	7
2. MONTAGE ET INSTALLATION	7
3. CONTRÔLES DES OPÉRATIONS	7
3.1 DESCRIPTION DES TOUCHES	7
3.2 DEFINITION DE L'ÉCHELLE	7
3.3 ALARMES	7
4. ENTRETIEN	7
5. CARACTÉRISTIQUES TECHNIQUES	7
6. ACCESSOIRES / PIÈCES DE RECHANGE	7
1. INTRODUCCIÓN	8
2. MONTAJE E INSTALACIÓN	8
3. CONTROLES DE FUNCIONAMIENTO	8
3.1 DESCRIPCIÓN DE LOS BOTONES	8
3.2 ESCALAS	8
3.3 ALARMAS	8
4. MANTENIMIENTO	8
5. CARACTERÍSTICAS TÉCNICAS	8
6. ACCESORIOS / REFACCIONES	8
1. EINFÜHRUNG	9
2. MONTAGE UND INSTALLATION	9
3. BEDIENUNGSELEMENTE	9
3.1 BESCHREIBUNG DER TASTEN	9
3.2 SKALA AUSWAHL	9
3.3 ALARME	9
4. WARTUNG	9
5. TECHNISCHE MERKMALE	9
6. ZUBEHÖR / ERSATZTEILE	9
1. 简介	10
2. 装配和安装	10
3. 操作控制	10
3.1 密钥的描述	10
3.2 比例选择	10
3.3 警报	10
4. 维护	10
5. 技术数据	11
6. ACCESSORIES / SPARE PARTS	11

1. Introduction

BOD measurement with BOD Sensor is based on the respirometric method that provides a direct measurement of the oxygen consumed by microorganisms from an air or oxygen-enriched environment in a closed atmosphere, at a constant temperature ($20\text{ °C} \pm 0.5$) and with a soft sample stirring.

2. Assembly and installation

Upon receipt and after having removed the packaging, please check the integrity of the instrument. The box includes:

- BOD Sensor
- Alkali holder
- Instruction manual

BOD Sensor is delivered with 2 lithium batteries installed: remove the insulating strip from the rear without opening it.

3. Operating controls

The display of the electronic measuring unit is normally turned off reducing the consumption of internal batteries. The touch of a key activates the display which automatically is turned off after three seconds. The key activates different functions if the unit is performing a measurement or if it is in memory code.

3.1 Description of keys

Out of a measuring cycle (before starting)

Key SET with display off	The scale is shown
Key SET with display on	The scale is changed
Key START with display off	The scale is shown
Key START with display on	The cycle of measurement is started

During a measuring cycle (after starting)

Key SET	The scale is shown
Key START	The current measure is shown (also after the 5-day period if Reset is not performed)
Key START kept pressed	The memory mode is entered (also after the 5-day period if Reset is not performed)
<u>Memory mode (it is active during a measuring cycle and also after if Reset is not performed)</u>	

Key SET	The stored value is shown
Key START	The order number of stored value is increased
If no values are stored, the display shows three points after a prolonged pressing of key START .	

Reset

Keys **SET+START** kept pressed Measurement cycle stopped, stored values deleted and the scale is shown

Before starting a new measuring cycle, it is requested to reset the stored values.

3.2 Scale selection

Scale	Volume of sample
0 ÷ 90 mg O ₂ /l	400 ml
0 ÷ 250 mg O ₂ /l	250 ml
0 ÷ 600 mg O ₂ /l	150 ml
0 ÷ 999 mg O ₂ /l	100 ml
0 ÷ 4000 mg O ₂ /l	25 ml

The scale is related to the expected BOD value.
The possibility of measuring diluted samples extends the range.

3.3 Alarms

--- The measured value overcomes the maximum value of the chosen scale.

bAt Battery replacement. The charge of batteries does not assure in this case the finishing of a complete cycle.

4. Maintenance

BOD Sensor warns about the required change of batteries when the unit is resetted and the charge is lower than a minimum threshold. Using a coin lever in the rear part of BOD Sensor and open it. Batteries are removed using the fingers. Introduce two fresh batteries at the same time paying attention to the polarity. Battery replacement deletes any stored value.

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Thank you for having chosen VELP!

Established in 1983, VELP is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. VELP works according to **ISO 9001**, **ISO14001** and **OHSAS 18001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

Analytical instruments

Elemental Analyzers
Digestion Units
Distillation Units
Solvent Extractors
Fiber Analyzers
Dietary Fiber Analyzers
Oxidation Stability Reactor
Consumables

Laboratory Equipment

Magnetic Stirrers
Heating Magnetic Stirrers
Heating Plates
Overhead stirrers
Vortex mixers
Homogenizers
COD Thermoreactors
BOD and Respirometers
Cooled Incubators
Flocculators
Overhead Shakers
Turbidimeter
Radiation Detector
Open Circulating Baths
Pumps

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VELP Scientifica Srl