EN - Instruction Manual
IT - Manuale di istruzioni
FR- Manuel d'instructions
ES - Manual de instrucciones
De - Bedienungsanleitung

CONTROLLER Advance

F201A0500

General Information / Informazioni Generali / Informations Générales / Información General / Allgemeine Hinweise

Before using the unit, please read the following instruction manual carefully.
Prima dell'utilizzo dello strumento si raccomanda di leggere attentamente il seguente manuale operativo.
Avant d'utiliser l'instrument, il est recommandé de lire attentivement le présent manuel d'instructions.
Antes de utilizar el instrumento, le recomendamos que lea con atención el siguiente manual de funcionamiento.
Bitte lesen Sie vor Inbetriebnahme des Geräts diese Bedienungsanleitung sorgfältig durch

Do not dispose of this equipment as urban waste, in accordance with EEC directive 2002/96/CE.
Non smaltire l'apparecchiatura come rifiuto urbano, secondo quanto previsto dalla Direttiva 2002/96/CE.
Ne pas recycler l'appareil comme déchet solide urbain, conformément à la Directive 2002/96/CE.
No tirar el aparato en los desechos urbanos, como exige la Directiva 2002/96/CE.
Dieses Gerät unterliegt der Richtlinie 2002/96/EG und darf nicht mit dem normalen Hausmüll entsorgt werden.

This unit must be used for laboratory applications indoor only. The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions. If the product is used in a not specified way by the manufacturer or with not specified accessories, product's safety may be compromised.

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This unit has been designed and manufactured in compliance with the following standards:

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IEC/EN 61010-1

Electrical equipment for laboratory use

UL 61010-1

General requirement - Canadian electrical code

CAN/CSA-C22.2 No.61010-1

VELP reserves the right to modify the characteristics of its products with the aim to constantly improving their quality. Nell'impegno di migliorare costantemente la qualità dei prodotti, VELP si riserva la facoltà di variane le caratteristiche. Dans le but d'améliorer constamment la qualité de ses produits, VELP se réserve le droit d'apporter des modifications aux caractéristiques de ceux-ci. VELP se reserva el derecho de modificar las características de sus productos con el objetivo de mejorar constantemente su calidad. VELP behält sich zum Zwecke der ständigen Verbesserung der Produktqualität das Recht auf Änderung der Geräteeigenschaften vor.

Safety Regulations / Norme di Sicurezza / Consignes de Sécurité / Advertencias de Seguridad / Sicherheitshinweise

The plug disconnects the instrument. Therefore, place the instrument where it can be quickly disconnected. / La spina è il mezzo di disconnessione dell'apparecchio. Pertanto, non posizionare l'apparecchio in modo che sia difficile azionare il mezzo di disconnessione. / Le bouchon est le moyen de déconnexion de l'appareil. Par conséquent, placer l'appareil où il peut être rapidement débranché. / El tapón es el medio de desconexión del dispositivo. No coloque el dispositivo en una forma que es difícil de desconectar. / Der Stecker trennt das Gerät. Daher Stellen Sie das Instrument, wo es schnell getrennt werden kann.

The values indicated on the rating plate of the instrument must correspond to those of the a.c. mains supply.

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Fasten the unit to the support rod using the fixing block. Position the instrument with a distance from the wall of 30 cm (at least).

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The unit is fitted with two fuses (2xT5 A L 250 V), found in the socket on the back. To replace one or more, disconnect the mains cable and, using a screwdriver, lift up the small cover on the fuse box.

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The working speed set on the instrument must be such as to avoid wobbling and/or splashes.

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Do not use with explosive or dangerous materials for which the equipment is not designed. The stirrer must not be used in explosive atmospheres, in bain-marie or to stir harmful liquids prior using protective measures according to the safety standards of the processed products and/or in force in the laboratories including personal protective equipment and the presence of an extraction hood which ensures at least 10-fold air change in accordance with the standards EN 14175 and DIN 12924.

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It is dangerous to run the motor free in air. Always fix the motor before turning the unit on.

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The solution may release toxic, dangerous or poisonous gases. Adequate safety measures must be taken, in accordance with the safety regulations in force, including the presence of hood and personal protective equipment (masks, gloves, goggles, etc.).

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The instrument contains a battery. Its substitution must be carried out by authorized Velp personnel only.

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CONTROLLER Advance is an innovative solution for control and data logging in process reactors. It is able to monitor the temperature in the reactor and the shaft rotation speed of brushless gear-motors for liquid volumes until 150L, even with viscosity changes. The motor is reversible.

All data logging features can be managed through VELP Ermes platform or USB by a PC.

Figure 1. CONTROLLER Advance

- A Display LCD
- B Display keys
- C Speed control knob
- D Lock key
- E Main switch
- F Motor circular connector 6P *
- G Motor circular connector 4P **

- H USB socket
- I Appliance inlet with fuse holder
- J Pt100 socket
- K RS232 socket

*Circular connectors Model MIL-C-505 MS3102A-14S-6

**Circular connectors Model MIL-C-505 MS3102A -14S-2

2 EN - Installation

- Unpacking
 - Check the integrity of the unit after unpacking.
- The box includes
 - Controller advance
 - Power supply cable
 - Instruction manual
- First installation
 - Place the unit on non-flammable surface
 - Fix the instrument to the support through the handle
 - Connect the instrument to one of the following brushless motors 4 poles, 24 V:
1/11HP (max 4.4 A) or 1/4HP (max 8,3 A) or 3/8HP (max 15 A)


Figure 2. Connector side view

- Make sure that the rating value of the instrument corresponds to the one of the power supply
- Ensure that the socket with grounding is compliant with the current safety norms and easy to reach. Use only the cable provided with the instrument.
- Insert the power cable into the socket and switch the instrument on.
- Select the motor type in the SET UP/Motor type menu.

NOTE: the cable can be substituted only by main cables with same features (T=70°C, connector C14).


Figure 3. Display

- | | |
|-----------|--|
| 1 | Wi-Fi symbol |
| 2 | Start Mode symbol |
| 3 | Vibration Sensor symbol |
| 4 | Hour |
| 5 | Current Speed |
| 6 | Set Speed |
| 7 | Timer or Time counter |
| 8 | Set Timer |
| 9 | Current Torque |
| 10 | Current temperature (only when Pt100 probe is connected) |
| 11 | Reverse button |
| 12 | Timer button |
| 13 | Menu button |
| 14 | Current rotation direction |
| 15 | Current Method indication |
| 16 | Motor type |

Commissioning	 ➤ Switch on the instrument using the main switch (E) ➤ Display (A) shows Welcome page and the main screen Figure 4. Main screen
Stirring	➤ Adjust speed set point by turning the speed control knob (C). As soon as the knob is moved, set rpm (6) becomes blue. ➤ Click the speed control knob (C) to start stirring. ➤ Speed increases until set point is achieved. ➤ A microprocessor ensures constant speed even when the viscosity changes (counter-reaction). ➤ Switch off the stirring by clicking the knob (C).
Timer	➤ Click Timer button (12) to select the timer. ➤ Set timer (8) - time becomes blue. Adjust timer by turning the speed control knob (C). ➤ Click the knob (C) to confirm. ➤ If the instrument is already working, timer countdown starts immediately, otherwise timer (7) is fixed as set timer (8) until stirring begins. ➤ If timer is not set, set timer (8) shows hh:mm:ss and timer (7) works as a counter.
CW / CCW	 ➤ Set the rotating direction by clicking CW/CCW button (11) and rotating the speed control knob (C). ➤ Click the knob (C) to confirm clockwise (CW) or counterclockwise (CCW) direction. Once confirmed, "CW/CCW" symbol becomes black. Figure 5. Set Rotation Direction
Torque	➤ The intensity of the torque is indicated on the main screen (9).
Lock	➤ Holding the Lock key (D) for 3 seconds, the instrument will lock its settings during operations. ➤ Unlock the control panel by holding the Lock key (D) for 3 seconds. ➤ If other buttons are clicked while the instrument is locked, the two LEDs aside lock button blink for many seconds.

5 EN - External Connections

USB	➤ for PC controlling, data logging using modbus protocol ➤ for software upgrading, using a dedicated program to requested by to service@velp.it and to be installed in the PC
RS232	➤ for interface with other instrument using RS232 protocol
Pt100	➤ This model has Pt100 connection in the back for fluid temperature measurement (measuring range from -200°C to + 550°C)

5.1 Modbus protocol

Controller register assignment

Register address	Register name	Write	Byte Number	Function/explanation
1	Instrument	NO	2	Instrument model
2	Serial number	NO	6	Instrument serial number
5	Product code	NO	12	Instrument product code
11	Main board sw	NO	8	Main board sw version
15	Display board sw	NO	8	Display board sw version
24	State	NO	2	Instrument state
25	Alarm	NO	2	The instrument notifies an alarm
26	Gear	NO	2	Gear set (only for OHS200)
27	Speed	NO	2	rpm measured
28	Torque	NO	2	Ncm measured
29	Timer	NO	4	Residual timer or counter
31	Timer set point	NO	4	Timer set value
33	Speed set point	NO	2	Speed set value
48	Pt100 connection	NO	2	Pt100 connected to the instrument
49	Temperature	NO	2	Temperature indication
50	Speed limit	NO	2	Speed maximum value
51	Torque limit	NO	2	Torque maximum value
52	Ramp	NO	2	Setting of acceleration
53	Method	NO	2	Reading of method number
54	Method n steps	NO	2	Total number of steps
55	Method step	NO	2	Current step visualized
56	Method n loops	NO	2	Total number of loops
57	Method loop	NO	2	Current loop visualized
102	Gear	SI	2	Setting of gear (only for OHS200)
103	Speed	SI	2	Setting of the speed
104	Timer	SI	4	Setting of the timer
106	Motor stop	SI	2	Motor stops stirring
107	Motor start	SI	2	Motor starts stirring

Sampling time: 1s or more

Examples (CRC 16 bit for standard MODBUS)

Reading single register: Speed (rpm)

Request		Reply		
Field	(Hex)	Field	(Hex)	Description
Address	0x64	Address	0x64	
Control command	0x03	Control command	0x03	
High start address	0x00	Number of bytes	0x02	
Low start address	0x1B	High register value	0x01	340 (RPM)
Number of High registers	0x00	Low register value	0x54	
Number of Low registers	0x01	High CRC	0xF4	
High CRC	0xFD	Low CRC	0x23	
Low CRC	0xF8			

Reading multiple register: Serial number

Request		Reply		
Field	(Hex)	Field	(Hex)	Description
Address	0x64	Address	0x64	

Control command	0x03	Control command	0x03	
High start address	0x00	Number of bytes	0x06	
Low start address	0x02	High register value	0x31	12345 0x00 0x31 0x32 0x33 0x34 0x35
Number of High registers	0x00	Low register value	0x00	
Number of Low registers	0x03	High register value	0x33	
High CRC	0xAD	Low register value	0x32	
Low CRC	0xFE	High register value	0x35	
		Low register value	0x34	
		High CRC	0x0A	
		Low CRC	0x0A	

Writing single register: STOP rotation

Request				
Field	(Hex)			Description
Address	0x64			
Control command	0x06			
High start address	0x00			
Low start address	0x6A			
High register value	0x00			Stop rotation
Low register value	0x01			
High CRC	0x61			
Low CRC	0xE3			

Virtual serial port

Baudrate	9600
Bits	8
Stop Bit	1
Parity	None
Maximum number of registers for single request	24

5.2 RS232 protocol

Command	Value (range)	Command Options *	Function
SS	0 (OFF) 35 – 500 (RUN)	Q/N	Sets/Reads the speed of the motor in RPM.
TQ	None	Q	Reads the actual torque in Ncm sensed by the motor
PI	None	Q	Returns a string with product information (software version and serial number). e.g: OHS v1.0.7-1.1.6 SN_609379
MS	None	Q	Returns the controller's current status: MS1 = Motor stopped MS2 = Motor accelerating MS3 = Motor decelerating MS4 = Motor running at set RPM MS5 = Motor overloaded or motor stalled or driver fault condition detected
SN	None	Q	Returns the serial number of the controller.

*

- Q Query command is allowed
N Command optionally accepts a value

Protocol Syntax

Commands are case sensitive and must be in upper case.

The end-of-line character is a carriage return <CR>, or Hex 13. Do not add a linefeed to the end of line.

Query Format: <COMMAND><CR>

This example queries the motors sensed torque.

Example:

Command = TQ<CR>

Reply = TQ###.#<CR> Where ### is the sensed torque on the motor.

Command Format: < COMMAND ><Value><CR>

This example sets the motor rotational speed.


Example:

Command = SS350<CR>

Reply = SS350<CR>

Set up serial port

Baudrate [bps]	9600
Data bits	8
Stop bit	1
Parity	None

Serial cable pinout

6 EN - VELP Ermes Configuration

VELP Ermes is a revolutionary cloud platform that transforms and improves your laboratory experience by creating an ecosystem of instruments, people and data. The VELP Ermes Cloud Platform is able to reduce distances and accelerate scientific processes in total safety. In order to access on Ermes, you need to enable your VELP account by selecting "Configure your VELP Ermes account" at <http://www.velp.com/en/login>.

To be able to communicate, the instrument needs to be in the operating range of laboratory Wi-Fi (2.4 GHz) and be configured as follows:

- Switch on the CONTROLLER Advance and select AP in the menu Set-up Wi-Fi (see chapter 7.5).
- Using the PC/Tablet/Mobile phone, select the RC_SERIAL NUMBER available on the Wi-Fi list, in order to connect directly to the instrument.
- Open a browser on the PC/Tablet/Mobile phone and insert the address 192.168.240.1 to reach the configuration page. Insert "admin" "admin" when requested as user name and password.
- Set the parameters required to connect to your Wi-Fi (network name, password, security, mac address, etc.) according to your internal procedure and save. If necessary, contact your IT administrator.
- Select menu Ermes from the Service menu (see chapter 7.6) and proceed with the product registration from the VELP Ermes Cloud Platform. For more information see FAQ on VELP website.

NOTE: To access to VELP Ermes it is necessary to have a VELP account.

7 EN - Menu

Clicking Menu button (13) the following figure appears



Figure 6. Main Menu

It's possible to move among submenus rotating the speed control knob (C).

Enter in a submenu with a single click of the speed control knob (C) once it's highlighted in blue.

Press Home to go back to the main screen.

7.1 Menu Structure

Timer

CW/CCW

7.2 Methods

In this submenu it is possible to set 10 different methods.
Once set, a method becomes white in the method list and a blue bar appears on the left.
Press START to begin the method.


Figure 7. Methods


Figure 8. Methods set


Figure 9. Stop Method

When a Method is working, on the main screen the current method (15) is shown.
It's not possible to modify speed, gear, or the timer.
It's possible to navigate in menu.
To stop a method before its end, enter in the method list and click STOP (Figure).

7.2.1 Method set-up

In each method it's possible to set parameters in Figure . Once method parameters are set, they are highlighted in white with a blue bar on the left as shown in Figure .


Figure 10. Method parameters


Figure 11. Method parameters set

- 5 steps can be set for each method.
- All steps are programmable with speed, timer, ramp, and intermittent mode.
- A Method is considered set when at least speed and timer are set for one step.

The first three screenshots show the 'STEPS' screen with 5 steps. In the first, all values are default (rpm: 0, Time: 00:00, Ramp Type: 0, Interm. Mode: No). In the second, Step 1 rpm is set to 500. In the third, Step 1 rpm is 500 and Step 1 Time is set to 00:00.

Step ID	rpm	Time hh:mm	Ramp Type	Interm. Mode
1	0	00:00	0	No
2	0	00:00	0	No
3	0	00:00	0	No
4	0	00:00	0	No
5	0	00:00	0	No

Figure 12

Steps

The next three screenshots show further configuration: Step 1 rpm is 500 and Time is 00:00; Step 1 rpm is 500 and Time is 01:30; Step 1 rpm is 500, Time is 01:30, and Ramp Type is 1. The final three screenshots show: Step 1 rpm is 500, Time is 01:30, Ramp Type is 1, and Interm. Mode is No; Step 2 rpm is set to 0 and Time to 00:00; and finally, Step 2 rpm is 0, Time is 00:00, Ramp Type is 1, and Interm. Mode is No.

Step ID	rpm	Time hh:mm	Ramp Type	Interm. Mode
1	500	00:00	0	No
2	0	00:00	0	No
3	0	00:00	0	No
4	0	00:00	0	No
5	0	00:00	0	No

Step ID	rpm	Time hh:mm	Ramp Type	Interm. Mode
1	500	01:30	0	No
2	0	00:00	0	No
3	0	00:00	0	No
4	0	00:00	0	No
5	0	00:00	0	No

Step ID	rpm	Time hh:mm	Ramp Type	Interm. Mode
1	500	01:30	1	No
2	0	00:00	0	No
3	0	00:00	0	No
4	0	00:00	0	No
5	0	00:00	0	No

Step ID	rpm	Time hh:mm	Ramp Type	Interm. Mode
1	500	01:30	1	No
2	0	00:00	0	No
3	0	00:00	0	No
4	0	00:00	0	No
5	0	00:00	0	No

Step ID	rpm	Time hh:mm	Ramp Type	Interm. Mode
1	500	01:30	1	No
2	0	00:00	1	No
3	0	00:00	0	No
4	0	00:00	0	No
5	0	00:00	0	No

- Rpm: set point speed
- Time: countdown performed for each step visualized also on the main screen
- Ramp: if is selected, Controller stirs for all the time at the speed set
if is selected, Controller reaches the speed set in time selected
- Interm. Mode: if YES is selected, but no intermittent mode is set for the method, Controller works in continuous mode.

Start Time	➤ It allows to set an hour at which the method starts. ➤ Enter in the menu, rotate speed control knob (C) till the desired hour. Click the knob to confirm. ➤ Even if a Start Time is set, the method begins if START button in method list is clicked.
Intermittent Mode	➤ It allows to set stirring period and stop period alternatively. ➤ Enter in the menu, rotate speed control knob (C) till the desired working time. Click to confirm. ➤ Rotate speed control knob (C) till the desired pause time. Click to confirm. ➤ Click again the knob to modify set values.
Loop	➤ It allows to repeat the whole method N times. ➤ Enter in the menu, rotate speed control knob (C) till the desired value. Click to confirm.

Click RESET to erase all method values.

7.3 Graph

In this submenu is possible to visualize the graph of the current analysis.

Click RESET to erase the graph. Only the last 60 minutes are shown. Click ZOOM to see the last 3 minutes.


Figure 13. Graph with method


Figure 14. Graph without method

When a method is set graph appears as in [Figure 13](#). When a method is not set, graph appears as in [Figure 14](#).

7.4 Safety

In this submenu it's possible to set all limits linked to the safety depending on the working conditions

Speed Limit	➤ It allows to set the speed full scale. ➤ The maximum speed value is set by default. ➤ Enter in the menu, rotate speed control knob (C) by 100 rpm steps till the desired value. Click to confirm.
Torque Limit	➤ It allows to set the torque limit. ➤ The maximum value is set by default. ➤ Enter in the menu, rotate speed control knob (C) by 10 Ncm steps till the desired value. Click to confirm.
Temperature Limit	➤ It allows to set a speed reduction if Pt100 detects a certain temperature increase (or decrease) in a defined time.


Figure 15. Temperature Limit

- Enter in the menu, rotate speed control knob (C) till the delta temperature desired value. Click to confirm.
- Rotate the knob (C) to set the interval time in which the delta temperature has to be considered. Click to confirm.
- Rotate the knob (C) to set the percentage of speed reduction desired if the selected delta temperature is detected in the set interval time. Click to confirm.
- Speed reduction can be set by 10% steps.

Acceleration

- It allows to choose among 3 different acceleration types when the instrument starts to stir or when a higher set point value is set.
- Enter in the menu, select Slow, Medium, or Fast depending on the customer application.
- Click speed control knob (C) to confirm.

Vibration sensor

- It allows to set a vibration sensitivity level.
- Enter in the menu, rotate speed control knob (C) to select OFF, Low, Medium or High. Click to confirm.
- When vibration sensor level is set different than OFF,  appears on the blue upper bar (3).

7.5 Set-up

Motor Type

- It allows to set the type of the motor connected to the instrument
- Enter in the menu, rotate speed control knob (C) to select the motor connected physically. Click to confirm.
 - **Motor 5:1 – 500 rpm:** 4 poles, 24V, 1/4HP, Max 8.3A, Max rpm 500
 - **Motor 6:1 – 417 rpm:** 4 poles, 24V, 1/11HP, Max 4.4A, Max rpm 417
 - **Motor 10:1 – 250 rpm:** 4 poles, 24V, 1/4HP, Max 8.3A, Max rpm 250
 - **Motor 5:1 – 500 rpm – 3/8:** 4 poles, 24V, 3/8HP, Max 15A, Max rpm 500

Start Mode


- It allows to set instrument's restart mode in case of blackout or power loss.
- Enter in the menu, rotate speed control knob (C) to select Stop or Run. Click to confirm.
 - **Stop:** when the instrument is switched on, it's requested a click of the control knob to start stirring.
 - **Run:** when the instrument is switched on, it restarts to work with the last set point set.
- If Stop is selected, A appears on the blue upper bar (2).
- If Run is selected, B appears on the blue upper bar (2).

Wi-Fi Set-up

- It allows to switch on wi-fi module for IoT transmission.
- Enter in the menu, rotate speed control knob (C) to select ON or OFF. Click to confirm.
- When wi-fi is ON,  appears on the blue upper bar (1).

Wi-Fi Information

- It describes all wi-fi information (MAC address, Power of signal, wi-fi name).

Display	➤ It allows for brightness to be set or if torque and temperature should be displayed on the main screen. ➤ Enter in the menu and scroll with speed control knob (C). Click to enter in submenus. ➤ Brightness: enter in this submenu, rotate speed control knob (C) to select the desired brightness value. Click to confirm. ➤ Brightness can be set by 10% steps. ➤ Torque: enter in this submenu, rotate speed control knob (C) to choose between ON or OFF if the display of torque on the main screen is desired or not. Click to confirm. ➤ Temperature: enter in this submenu, rotate speed control knob (C) to choose between ON or OFF if the display of temperature on the main screen is desired or not. Click to confirm ➤ If torque and temperature are displayed on the main screen, the display menu appears as in Figure 16. ➤ When temperature is displayed but Pt100 probe is not inserted in the instrument, the main screen appears as in Figure 17.  Figure 16. Torque and temperature shown  Figure 17. Main without temperature probe
Time & Date	➤ It allows to set hour and date. ➤ Enter in the menu and select Set Time. Rotate speed control knob (C) till the right time. Click to confirm. ➤ Return to the previous page clicking BACK button. ➤ Select Set Date. Rotate speed control knob (C) till the right day. Click to confirm ➤ Repeat the operation for month and year.
Temperature Unit	➤ It allows to choose the temperature and the torque unit that have to be visualized
Language	➤ It allows to select the interface language. ➤ Enter in the menu, rotate speed control knob (C) to select language. Click to confirm.
Protocol	➤ It allows to choose type of protocol between Modbus and RS232
7.6 Service	
Reset Torque	➤ It allows to reset the current torque. ➤ Click with the speed control knob (C) to reset torque.
Real Torque	➤ It allows to return to real torque value. ➤ Click with the speed control knob (C) to return to real torque value. ➤ Negative torques can't be visualized anymore.
Pt100 Alignment	➤ It allows for the alignment of the Pt100 probe to a reference thermometer. ➤ Enter in the menu, rotate speed control knob (C) to select the desired alignment value (from -10.0°C to 10.0°C). Click to confirm. ➤ Alignment sensitivity 0.1°C.

Reset Parameters	➤ It allows to return to default value for all functions. ➤ Click with the speed control knob (C) to reset parameters.
Update Software	➤ It allows to update the device with a new software version. ➤ Software version, dedicated program to be installed in the PC, and software installation guide have to be requested by e-mail to service@velp.it.
Counter	➤ It allows to see the number of instrument working hours.
Ermes	➤ It allows to connect the instrument to Ermes Cloud.

8 EN - Error messages

When the display shows an error message, the stirring function stops automatically.

Error code	Cause
AL1	Motor doesn't start stirring
AL2	High internal motor temperature
AL3	Motor overload
AL4	High driver temperature
AL5	Safety relay fault
AL10	Vibrations too high
AL11	Temperature too high (only with Pt100 inserted)
AL12	Temperature too low (only with Pt100 inserted)

To remove the error message, switch the instrument off and then switch on it.

If alarm persists on the display, please contact VELP Scientifica technical service department. service@velp.it

9 EN - Maintenance

Maintenance	➤ No routine or extraordinary maintenance is necessary; ➤ Repairs must be carried out by authorized Velp personnel only; ➤ Instrument must be transported in its original packaging any indications present on the original packaging must be followed (e.g., palletized);
Cleaning	➤ Disconnect the unit from the power supply and use a cloth dampened with a non-flammable non-aggressive detergent.

10 EN - Technical data

	Model	F201A0500
	Power supply	115-230V (+/-10%) – 50-60 Hz
	Dimensions (WxHxD)	90x180x170 mm (3.54x7.10x6.70 in)
	Weight	2,5 kg (5,51 lb)
	Power input	500 W
	Max power output	280 W (15A)
	Construction material (structure)	Stainless steel
General features	Working in continuous	Admitted
	Settable restart modality	Stop or work
	Noisiness	<< 60 dBa
	Environmental temperature admitted	+5...+40 °C
	Storage temperature admitted	-10...+60 °C
	Max humidity	80%
	Overvoltage category	II
	Pollution degree CEI EN61010-1	2
	Max altitude	2000 m
	Programmable speed range	35-500rpm (max speed depends to the motor selected)
Stirring	Speed selection	1 rpm step
	Rpm accuracy	± 1
	Stirring alarm	Motor fault
Torque	Max torque admitted	621 Ncm
Counters	Board counter	Working hours

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








We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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EN 61010-1	Safety requirements for electrical equipment for measurement, control, and laboratory use Part 1: General requirements
EN 62479	Assessment of the compliance of low-power electronic and electrical equipment with the basic restrictions related to human exposure to electromagnetic fields (10 MHz to 300 GHz)
ETSI EN 301 489-1 V2.2.0	Electromagnetic compatibility (EMC) standard for radio equipment and services; Part 1: Common technical requirements
ETSI EN 301 489-17 V3.2.0	ElectroMagnetic Compatibility (EMC) standard for radio equipment and services; Part 17: Specific conditions for Broadband Data Transmission Systems
EN 61326-1	Electrical equipment for measurement, control and laboratory use - EMC requirements - Part 1: General requirements
ETSI EN 300 328 V2.1.1	Wideband transmission systems; Data transmission equipment operating in the 2,4 GHz band

and satisfies the essential requirements of the following directives:

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2006/42/EC	Machinery directive
2014/53/EU	Radio Equipment directive (RED)
2015/863/EU (RoHS III)	Restriction of the use of certain hazardous substances
2012/19/EU (WEEE)	Waste of electric and electronic equipment

plus modifications / più modifiche / plus modifications / más sucesivas modificaciones / in der jeweils gültigen Fassung / 加修改.

15 Declaration of conformity

We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following regulations:

S.I. 2017/1206 Radio Equipment Regulations 2017

according to the relevant designated standards:

EN 61010-1 Safety requirements for electrical equipment for measurement, control, and laboratory use
Part 1: General requirements

EN 62479 Assessment of the compliance of low-power electronic and electrical equipment with the
basic restrictions related to human exposure to electromagnetic fields (10 MHz to 300 GHz)

ETSI EN 301 489-1 V2.2.0 Electromagnetic compatibility (EMC) standard for radio equipment and services;
Part 1: Common technical requirements

ETSI EN 301 489-17 V3.2.0 ElectroMagnetic Compatibility (EMC) standard for radio equipment and services;
Part 17: Specific conditions for Broadband Data Transmission Systems

EN 61326-1 Electrical equipment for measurement, control and laboratory use - EMC requirements -
Part 1: General requirements

ETSI EN 300 328 V2.1.1 Wideband transmission systems; Data transmission equipment operating in the 2,4 GHz
band

and satisfies the essential requirements of regulations:

S.I. 2008/1597 Supply of Machinery (Safety) Regulations 2008

S.I. 2012/3032 Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic
Equipment Regulations 2012

S.I. 2013/3113 Waste Electrical and Electronic Equipment Regulations 2013

plus modifications.

Thank you for having chosen VELP!

Established in 1983, VELP is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. VELP works according to **ISO 9001**, **ISO14001** and **ISO 45001** Quality System Certification. Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

Analytical instruments

Elemental Analyzers
Digestion Units
Distillation Units
Solvent Extractors
Fiber Analyzers
Dietary Fiber Analyzers
Oxidation Stability Reactor
Consumables

Laboratory Equipment

Magnetic Stirrers
Heating Magnetic Stirrers
Heating Plates
Overhead stirrers
Vortex mixers
Dispersers
COD Thermoreactors
BOD and Respirometers
Cooled Incubators
Flocculators
Overhead Shakers
Turbidimeter
Open Circulating Baths
Pumps

VELP Scientifica Srl
20865 Usmate (MB) ITALY
Via Stazione, 16
Tel. [+39 039 62 88 11](tel:+39039628811)
Fax. [+39 039 62 88 120](tel:+390396288120)

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