Instruction Manual
Manuale di istruzioni
Manuel d'instructions
Manual de instrucciones
Bedienungsanleitung

DLH - DLS Overhead stirrers

F201A0157, F201A0155

**General Information / Informazioni Generali / Informations Générales / Información General /
Allgemeine Hinweise**

Before using the unit, please read the following instruction manual carefully.

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Do not dispose of this equipment as urban waste, in accordance with EEC directive 2002/96/CE.

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This unit must be used for laboratory applications only.

The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions.

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This unit has been designed and manufactured in compliance with the following standards:

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Safety requirements for electrical equipment for measurement, control and for laboratory use

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Electrical equipment for laboratory use

General requirement - Canadian electrical code

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The unit is fitted with two fuses (2xT1,6 A L 250 V), found in the socket on the back. To replace one or more disconnect the mains cable and, using a screwdriver, lift up the small cover on the fuse box.

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It is responsibility of the user appropriately decontaminate the instrument in case of dangerous substances fall on or in it.

It is also responsibility of the user to use safety substances for cleaning or decontaminating, which do not react with internal parts of the instrument or with the material contained in it. In case of doubts on the compatibility of a cleaning solution, contact the manufacturer or local distributor.

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The solution may release toxic, dangerous or poisonous gases. Adequate safety measures must be taken, in accordance with the safety regulations in force, including the presence of hood and personal protective equipment (masks, gloves, goggles, etc.).

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The DLS and DLH overhead stirrers with electronic speed control and advanced safety features are able to satisfy the most diverse needs of the laboratory in terms of both viscosity and volume. The chuck is keyless and ensures a secure and easy-to-fit mount for stirring shafts of up to 10mm in diameter. The stirring speed is electronically controlled. The shaft bearings and belt drive ensure silent running. The stirring time is controlled electronically, values of up to 999 minutes and 59 seconds can be set and the time remaining displayed on the LCD. The relative torque value can also be displayed.


Fig. 1

2. Assembly and installation

Check the integrity of the unit after unpacking. The box includes:

- Overhead stirrer
- Instruction manual
- Allen key
- Allen head screw
- Shaft
- Power supply cable

2.1 Electrical connection

After unpacking the instrument, place the unit on the laboratory bench. Before connecting the instrument to the power supply, make sure that the values on the rating plate correspond to those of the power supply.

Ensure that the socket and the relative cut-off device conform to current safety norms and are easy to reach.

2.2 Start-up

Assemble the shaft using the Allen head screw. Fasten the instrument securely to the support rod (A00001300) with the double clamp (A00001301). Secure the receptacle using the ribbon clamp (A00001302) to avoid damage.

3. Operating controls

After a few seconds the instrument model and the software version will appear on the display. When the instrument is in use the graphic display shows the real chuck speed (rpm) and the set speed (rpm set). Speed and time functions are controlled by an encoder. Turn the encoder knob to select stirring speeds of from 50 to 2000 rpm. Press the encoder knob to activate the timer and turn the knob to select the stirring time. Stirring times of up to 999':59" can be selected. Press the encoder knob to start count-down. At the end of the set time stirring stops and "END" will appear on the display.

NOTE: the speed setting is saved when the instrument is turned off or if there is a power failure. When the instrument is turned on again, the value flashes for approximately 8 seconds before stirring starts. During this time the speed setting can be changed using the stirring speed knob.

3.1 Torque indicator

Keep the encoder knob pressed for 3 seconds to display the torque trend value, e.g. the value in Ncm of the effective torque on the stirring shaft. To reset the torque value and measure the relative torque, keep the knob pressed for 3 seconds.

NOTE: The torque value indicated is not an absolute value. The reading indicates the variation in the torque value compared to a set value defined by the operator. The most precise torque trend values are obtained once the instrument has reached the optimum working temperature. Optimum time 10-15 minutes.

3.2 Overload protection

The current and voltage are electronically limited in order to safeguard the motor and the electronic parts against electric surges. In the event of abnormal conditions arising the unit shuts down automatically.

Over Load Current surges

The stirrer can compensate for eventual current surges for a limited period of time; in this case  appears in the top rh corner of the display. If the situation persists, the speed of the motor is gradually reduced, when a speed of 0 is reached the unit shuts down.

Over Curr Rpm error

This error occurs if the chuck is accidentally blocked or if the number of revs exceeds the max. permitted value.

Over Temp Temperature error

This error occurs if the unit is used for heavy-duty work for prolonged periods of time and at an ambient temperature that exceeds the max. permitted temperature. Before resetting the instrument allow it to cool down.

Over Volt Voltage surge

The software limits the power supply of the motor to the nominal value of the motor itself (60V).

A high speed setting or a heavy load on the chuck may prompt a higher voltage. In this case the voltage is limited to the nominal value of the motor (60 Volt) and the number of revs is reduced accordingly. In this case, ▲ appears in the top right corner of the display. If for any reason the power supply of the motor exceeds the maximum threshold value of 70V, "OVERLOAD" the unit will shut down.

NOTE: to reset the unit remove the load from the chuck and turn the unit off and on again using the mains switch.

3.3 Set Start Mode

The instrument can start to work in two different modes:

- Stop: when the instrument is switched on through the main switch, 0000 rpm is the speed visualized on the display. The operator has to set the speed manually.
- Run: when the instrument is switched on, it restarts to work reaching the last set point without any operator's action.

The default Start Mode is "Run". It's possible to switch to "Stop" through the following steps:

- Switch on the instrument through the main switch. Set the speed to 0 rpm.
- Click the control knob to highlight the timer.
- Keep the control knob pressed for 3 seconds.
- The menu Start Mode is visualized. Rotate the control knob to switch from "Run" to "Stop".
- Click the knob to confirm the choice.

Repeat the previous operations to set "Run" Start Mode again.

4. Maintenance

No routine or extraordinary maintenance is necessary apart from periodically cleaning the unit as described in this manual. In compliance with the product guarantee law, repairs to our units must be carried out in our factory, unless previously agreed otherwise with local distributors. The instrument must be transported in its original packaging and any indications present on the original packaging must be followed (e.g. palletized).

It is the responsibility of the user, to properly decontaminate the unit in case of hazardous substances remaining on the surface or interior of the device. If in doubt about the compatibility of a cleaning or decontamination product, contact the manufacturer or distributor.

4.1 Cleaning

Disconnect the unit from the power supply and use a cloth dampened with a mild, non-flammable detergent.

5. Technical data

Code Model	F201A0157 DLH	F201A0155 DLS
Admitted power supply	110-230V, 50/60Hz (+/-10%)	110-230V, 50/60Hz (+/-10%)
Max. input / output power	190 W / 280 W	120 W / 180 W
Dimensions (WxHxD)	80x230x196 mm (3.1x9.0x7.7 in)	80x215x196 mm (3.1x8.5x7.7 in)
Weight	3 Kg (6,6 lbs)	2,5 Kg (5,5 lbs)
Speed range at nominal load	50 – 2000rpm	50 – 2000rpm
Speed indicator / resolution / accuracy	LCD / 1rpm / ± 10 rpm ≤ 300 rpm; $\pm 5\%$ > 300 rpm	LCD / 1rpm / ± 10 rpm ≤ 300 rpm; $\pm 5\%$ > 300 rpm
Max torque stirrer shaft	80Ncm	40Ncm
Max stirring quantity H ₂ O	40 litres	25 litres
Timer with automatic cut-off	From 000':00" to 999':59"	From 000':00" to 999':59"
Clamping chuck clamping range	From 1 to 10mm	From 1 to 10mm
Overload and overtemperature protection	Automatic cut-off	Automatic cut-off
Environmental ambient temperature	5 - 40°C (41 - 104°F)	5 - 40°C (41 - 104°F)
Storage temperature range	-10... + 60°C (14 - 140°F)	-10... + 60°C (14 - 140°F)
Max humidity	80%	80%
Noise Level	<<80dBa	<<80dBa
Operating mode	Continuous	Continuous
Overvoltage category CEI-EN 61010-1	II	II
Pollution degree / installation cat. CEI-EN61010-1	2	2
Max altitude installation	2000m (6561 ft) asl	2000m (6561 ft) asl
Fuse diam. 5x20mm	2xF1,6A L 250V	2xF1,6A L 250V

6. Accessories / Spare parts

A00001300	Support rod and base
A00001301	Double clamp
A00001302	Ribbon clamp
A00001304	Stirring shaft with floating blade
A00001305	Stirring shaft with folding blade
A00001306	Stirring shaft with fixed blade
A00001307	Stirring shaft with propeller
A00001308	Stirring shaft with paddle, 6 holes
A00001309	Stirring shaft with turbine
A00001310	Stirring shaft with turbo propeller
A00001311	Stirring shaft with anchor
10005266	Knob 26D blue
10000132	Gasket for overhead stirrer

			
A00001300	A00001301	A00001302	A00001304
			
A00001305	A00001306	A00001307	A00001308
			
	A00001309	A00001310	A00001311

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
We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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2002/96/CE (RAEE)

EN61010-2-051(2003)

and satisfies the essential requirements of the following directives:

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- Machinery directive 2006/42/EC / Macchine 2006/42/CE / Machines 2006/42/CE / Máquinas 2006/42/CE / Maschinen 2006/42/EG

- Low voltage directive 2014/35/EU / Bassa tensione 2014/35/EU / Basse tension 2014/35/EU / Baja tensión 2014/35/EU / Niederspannung 2014/35/EU

- Electromagnetic compatibility directive 2014/30/EU / Compatibilità elettromagnetica 2014/30/EU / Compatibilité électromagnétique 2014/30/EU / Compatibilidad electromagnética 2014/30/EU / Elektromagnetische Verträglichkeit 2014/30/EU

- plus modifications / più modifiche / plus modifications / más sucesivas modificaciones / in der jeweils gültigen Fassung.

Thank you for having chosen VELP!

Established in 1983, VELP is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. VELP works according to **ISO 9001**, **ISO14001** and **OHSAS 18001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

Analytical instruments

Elemental Analyzers
Digestion Units
Distillation Units
Solvent Extractors
Fiber Analyzers
Dietary Fiber Analyzers
Oxidation Stability Reactor
Consumables

Laboratory Equipment

Magnetic Stirrers
Heating Magnetic Stirrers
Heating Plates
Overhead stirrers
Vortex mixers
Homogenizers
COD Thermoreactors
BOD and Respirometers
Cooled Incubators
Flocculators
Overhead Shakers
Turbidimeter
Radiation Detector
Open Circulating Baths
Pumps

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VELP Scientifica Srl