Instruction Manual
Manuale di istruzioni
Manuel d'instructions
Manual de instrucciones
Bedienungsanleitung
指导手册

ECO 8 & ECO 25 Thermoreactors

F101A0127, F101A0125

**General Information / Informazioni Generali / Informations Générales / Información General /
Allgemeine Hinweise / 一般信息**

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This unit has been designed and manufactured in compliance with the following standards:

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Safety requirements for electrical equipment for measurement, control and for laboratory use

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IEC/EN 61010-1

Electrical equipment for laboratory use

UL 61010-1

General requirement - Canadian electrical code

CAN/CSA-C22.2 No.61010-1

VELP reserves the right to modify the characteristics of its products with the aim to constantly improving their quality.

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Thermo reaction can be carried out quickly and safely by placing the tubes into the thermoreactor heating block. The Start-Temp key offers the possibility to set 5 different work temperatures (70, 100, 120, 150 and 160 °C) whilst the Time key to set 4 different work times at the selected temperatures (30, 60, 120 or infinite). ECO 8 and ECO 25 offer the possibility to thermoreact respectively 8 and 25 tubes simultaneously (Ø 16 mm). ECO 8 model allows to hold also 1 test tube with Ø 22 mm.

2. Assembly and installation

Upon receipt and after having removed the packaging, please check the integrity of the instrument. The box includes:

- ECO 8 or ECO 25
- Power supply cord
- Safety cover
- Instruction manual

2.1 Electrical connection

Before connecting the instrument to the power supply, make sure that the values on the rating plate correspond to those of the power supply. Connect the unit to the power supply using the transformer supplied.

Ensure that the socket and the relative cut-off device conform to current safety norms and easy to reach.

2.2 Start-up

Use the on-off switch on the front of the unit to turn the unit on (position “–”). When switched on the unit all the leds blink for few seconds. When the leds are off the unit is ready to be programmed at different work temperatures and times.

3. Operating controls

ON - OFF SWITCH	It turns the unit on and off. If the switch is in the “0” position the unit is off; “–” means it is on.
START-TEMP	When Start-Temperature key is pressed, a work cycle will start with the pre-set values and the corresponding leds will light: Temperature 160 °C and Time 30 min. The next pressing permits to select a different temperature: 160, 70, 100, 120 and 150 °C. It's possible to select different work temperature until the heating block has not reached the set temperature and the count-down is not working. After that the key doesn't work anymore. If the selected temperature is lower than the running temperature of the heating block, (when the Start-Temp key is pressed) the work cycle will not start and all temperature leds will light intermittently, signalling an error. To correct the error, select a temperature that is higher than the running temperature or wait until the heating block cools down to a lower temperature than the temperature to be programmed.
TIME	When the Time key is pressed it's possible to select the following different time: 120, 60, 30 minutes or ∞ (infinite). When a time is selected, the corresponding led will light. After switching on the instrument, the Time key permits to set the time at the set temperature but not to start the reaction cycle. The run time will start when the heating block reaches the set temperature.
TEMP. LED	When a led lights it means that a work cycle is running at the selected temperature. If the leds are off a work cycle is not in progress.
TIME LED	When a pre-set time is selected the corresponding led will light. When the led blinks it means that the selected temperature has been reached and the count-down has started.

3.1 Heating block temperature verification

The temperature is controlled electronically which ensures stability without temperature oscillations around the set point. The Pt 100 probe provides high precision and does not require calibration because the internal software performs a self-calibration every time the instrument is turned on. Check the heating block temperature with a 3mm probe using the VELP Scientifica probe adaptor.

3.2 Power interruptions

The unit will register a black out or a lack of power supply as a manual stop. The work cycle will have to be restarted.

3.3 Safety features

The unit is protected against overheating by a thermostat that turns off the power supply when the temperature exceeds the safety temperature. Temperature setting is continuously controlled by electronic equipment that stops thermoregulation if anomalous measured are sensed. In case, all the leds start blinking and an acoustical warning signal sounds.

3.4 End of work cycle

At the end of the count-down the work cycle will stop automatically. To stop in every moment the work cycle turn off the general switch. When the work session is completed, should the unit remain connected to the power supply, switch it off using the on-off switch on the frontal panel.

4. Maintenance

No routine or extraordinary maintenance is necessary apart from periodically cleaning the unit as described in this manual. In compliance with the product guarantee law, repairs to our units must be carried out in our factory, unless previously agreed otherwise with local distributors. The instrument must be transported in its original packaging and any indications present on the original packaging must be followed (e.g. palletized).

It is the responsibility of the user, to properly decontaminate the unit in case of hazardous substances remaining on the surface or interior of the device. If in doubt about the compatibility of a cleaning or decontamination product, contact the manufacturer or distributor.

4.1 Cleaning

Disconnect the unit from the power supply and use a cloth dampened with an non-inflammable non-aggressive detergent.

5. Technical data

Code Model	F101A0127 ECO 8	F101A0125 ECO 25
Admitted power supply	115 V - 230 V / 50–60 Hz	115 V - 230 V / 50–60 Hz
Power	135 W	400 W
Dimensions (WxHxD)	135x95x230 mm	155x95x275 mm
Weight	2 Kg	3,6 Kg
Holes and diameter	8 (Ø16 mm) + 1 (Ø22 mm)	25 (Ø16 mm)
Construction material	Epoxy painted aluminum cast	
Thermoregulation	Pid microprocessor	
Selectable temperatures	70 - 100 - 120 - 150 - 160 °C	
Selectable times	30 - 60 - 120 min. - infinite	
Set temperature reached	Acoustic signal	
Count-down	Visual signal	
End of cycle	Acoustic signal	
Overtemperature	On the block at 200 °C	
Heating block temperature stability and homogeneity	± 0.5 °C	
Temperature precision	± 1.0 °C	
Time to reach 150 °C	12 min.	
Enviromental ambient temperature	5 - 40 °C (41 - 104°F)	
Storage temperature range	-10... + 60 °C (14 - 140°F)	
Max humidity	85%	
Pollution degree	2	

6. Accessories / Spare parts

CA0091666	Test tube for sample decomposition 22 mm external diameter, NS 19/26 cone with glass cap (ECO 8)
CA0091667	Condenser 200 mm type KS with 3 meters polyethylene tube (ECO 8)
10000002	Absorption attachment for condenser (ECO 8)
A00001039	Test tube extractor (ECO 25)
CM0091680	Set of 20 glass test tubes with internal Ø 14 mm and external Ø 16 mm (ECO 8 and ECO 25)
CA0091636	Holder for 12 round glass cells (ECO 8 and ECO 25)
A00000362	Adapter x probe Ø3mm for DK/DKL/ECO

10000667	ECO 8 Retarded fuse 5x20 3A
10000670	ECO 25 Retarded fuse 5x20 8A
10000232	Foot

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
We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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EN 61326-1
2015/863/EU (RoHS III)

2012/19/EU (WEEE)

and satisfies the essential requirements of the following directives:

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9. Declaration of conformity

We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following regulations:

S.I. 2016/1101 Electrical Equipment (Safety) Regulations 2016

S.I. 2016/1091 Electromagnetic Compatibility Regulations 2016

according to the relevant designated standards:

EN 61010-1 Safety requirements for electrical equipment for measurement, control, and laboratory use
Part 1: General requirements

EN 61010-2-010 Particular requirements for laboratory equipment for the heating of material

EN 61326-1 Electrical equipment for measurement, control and laboratory use - EMC requirements -
Part 1: General requirements

and satisfies the essential requirements of regulations:

S.I. 2008/1597 Supply of Machinery (Safety) Regulations 2008

S.I. 2012/3032 Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic
Equipment Regulations 2012

S.I. 2013/3113 Waste Electrical and Electronic Equipment Regulations 2013

plus modifications.

Thank you for having chosen VELP!

Established in 1983, VELP is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. VELP works according to **ISO 9001**, **ISO14001** and **ISO 45001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

Analytical instruments

Elemental Analyzers
Digestion Units
Distillation Units
Solvent Extractors
Fiber Analyzers
Dietary Fiber Analyzers
Oxidation Stability Reactor
Consumables

Laboratory Equipment

Magnetic Stirrers
Heating Magnetic Stirrers
Heating Plates
Overhead stirrers
Vortex mixers
Dispersers
COD Thermoreactors
BOD and Respirometers
Cooled Incubators
Flocculators
Overhead Shakers
Turbidimeter
Open Circulating Baths
Pumps

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VELP Scientifica Srl