Instruction Manual
Manuale di istruzioni
Manuel d'instructions
Manual de instrucciones
Bedienungsanleitung
指导手册

JLT4 - JLT6 Flocculators

F105A0108, F105A0109

General Information / Informazioni Generali / Informations Générales / Información General / Allgemeine Hinweise / 一般信息

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This unit has been designed and manufactured in compliance with the following standards:

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Safety requirements for electrical equipment for measurement, control and for laboratory use

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IEC/EN 61010-1

Electrical equipment for laboratory use

UL 61010-1

General requirement - Canadian electrical code

CAN/CSA-C22.2 No.61010-1

VELP reserves the right to modify the characteristics of its products with the aim to constantly improving their quality.

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Coagulation and flocculation of wastewaters

Chemical coagulants are added to wastewaters for the primary purpose of removing suspended solids and also phosphorus. The addition of chemicals reduces also heavy metal concentrations and improves disinfection efficiency. Chemical coagulation can be obtained by adding to wastewaters defined amounts of lime (calcium hydroxide), alum (aluminium sulphate) or iron salts (ferric or ferrous). The coagulation-sedimentation process typically involves:

- injection and mixing of the coagulant that neutralize the electric charge, prevalently negative, on suspended particles. Good and rapid mixing of the coagulant and wastewater is important to ensure efficient use of the chemical. Typical detention times in plant basins for turbulent mixing are 15 to 120 seconds.
- agglomeration of the coagulated particles into large settable flocs. This agglomeration is accomplished by stirring the water slowly so to allow the solids formed by the addition of chemicals to grow in size so that they can be removed by gravity settling. The stirring is usually performed by slowly rotating paddles during 10 to 30 minutes.
- sedimentation of the flocculated material by gravity separation. The solids are allowed to settle by gravity from the chemically treated standing wastewater. The effluent from the chemical settling tank is very clear and does not usually require further filtration.

The quantity of chemical coagulant required to achieve the test varies with time and from wastewater to wastewater. Typical coagulant doses are: 75 to 250 g/m³ for alum, 45 to 90 g/m³ for ferric chloride, 200 to 400 g/m³ for lime.

Reagents

The choice and dosing of the chemical coagulant to be adopted for the removal of suspended solids from wastewaters derive from the results of laboratory evaluations by the so-called Jar test. Multiple stirrers with reproducible stirring speeds allow to adopt standard conditions for the test, that are the basic requirement for reproducible results.

The standard conditions are related to:

- dimensions and shape of vessel and stirring paddle
- time and speed of stirring during turbulent mixing
- evaluation criteria for the obtained results.
- volume of wastewater
- time and speed of stirring during flocculation

The most diffused adopted conditions are:

- 1000 ml glass beakers (jars), tall form, Ø 105 mm
- stirring paddles 25 mm high, 75 mm wide, 1 mm thickness
- turbulent stirring after chemical addition: 120 rpm, 120 seconds
- first evaluation of results after 5 minutes of sedimentation
- 600 ml wastewater samples
- stirring height of paddles, middle height of sample
- slow speed flocculation: 30 rpm, 25 minutes

The results can be evaluated on the basis of different criteria:

- floc dimension evaluation with numerical degree (0=no flocs, 2=very small, 4=small, 6=medium, 8=large, 10=very large).
- time from the addition of chemicals to first appearing of flocs.
- evaluation of residual turbidity of supernatant, after a determined sedimentation time, by turbidimeter.
- measurement of electro-kinetic potential of suspended particles on a sample taken immediately after the addition and mixing of chemicals. Coagulations should be easier when zeta potential of particles is very low.
- evaluation of filterability of clarified water by standardized membrane filters under pressure. The reduction of water flow with time is related to the degree of clogging of filters due to residual unsettled suspended matter.

Other measurements are related to the pH of water after chemical treatment (lime raises the value while iron and aluminium salts lower it) and the temperature of water at which the experiments are performed.

References: American Society for Testing Materials Norm ASTM D-2035-G4 T; Passino R. and M.Beccari (1970) "Standardization of jar test in coagulation-flocculation processes for the removal of inorganic turbidities (in Italian). *Acqua e Aria*, Milano, 1-10, October 1970"; Degremont (1978) "Etude de la coagulation et de la flocculation des eaux. Methode n.703, 948-950. Memento technique de l'eau. 1200 pages, Paris"; U.S. Association of Environmental Engineering Professors (1972) "Environmental Engineering Unit Operations and Unit Processes Laboratory Manual. J.T. O'Connor Ed. 350 ppg".

Leaching test

Leaching tests are used to simulate the behaviour of a waste in a dump that is submitted to the leaching activity of rain water. There are two different methods that use leaching solutions of different type:

- acetic acid, proposed for dumps receiving both organic and inorganic wastes;
- carbon dioxide saturated water, for dumps receiving only inorganic wastes.

Waste sample to be examined: if a liquid phase is present, it must be removed by filtration or centrifugation and preserved in a refrigerator at 2-4 °C. The granulometry of the solid phase, must be lower than 9.5 mm (standard sieve).

Waste extraction: 20-25 g of the solid phase are weighed with a precision of ± 0.1 g.

Method a) A volume of distilled water corresponding to 16 times the weight of the sample is added.

The pH of the mixture is kept to 5±0.2 using 0.5 N acetic acid. The stirring is performed for 24 hours, controlling the pH every 15 minutes during the first hour and then every hour. The pH is corrected to 5±2 using again 0.5 N acetic acid. If after 24 hours of stirring the pH is higher than 5.2, a new correction is made and the stirring prolonged for other 4 hours.

Method b) The leaching solution is given by a volume corresponding to 20 times the sample weight of distilled water saturated by carbon dioxide for not less than 15 minutes. The stirring is performed during 6 hours without any pH correction.

Stirring speed: must be able to maintain continuously mixed the solid phase with the leaching solution.

Analysis of extract: the leaching solution and the liquid phase if present in the original sample are filtered by a membrane filter of 0.45 µm porosity and submitted to analysis for the single metals by the methods commonly used for waters.

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1. Introduction

During chemical-physical treatments of waters aimed to coagulate and precipitate pollutants, chemical coagulants are evaluated and their dosing optimized by jar tests. In relation to the variable composition of wastewaters the chemical dosing and the test frequency are to be adapted to local conditions. JLT Series is characterized by a settable speed of rotation, from 10 to 300 rpm, with the real speed showed on the relative display. The steel structure epoxy painted is studied and tested in order to give the instrument an excellent resistance against the attack of chemical agents and corrosion in general. The stainless-steel stirring rods are adjustable in height by a self-blocking chuck. The examined samples can be lighted by a back fluorescent lamp controlled by a switch located on front panel.

2. Assembly and installation

Upon receipt and after having removed the packaging, please check the integrity of the instrument. The box includes:

- Flocculator
- Plug for power supply
- Switching VELP 100-240V/12V without plug
- Instruction manual

2.1 Electrical connection

Before connecting the instrument to the power supply, make sure that the values on the rating plate correspond to those of the power supply. Connect the unit to the power supply using the transformer supplied.

Ensure that the socket and the relative cut-off device conform to current safety norms and easy to reach.

2.2 Start-up

Place the beakers containing the sample and then, set the speed and time using + or -.

3. Operating controls

The start-up is empowered by the keys + and -. It's possible to select the speed from 10 to 300 rpm. The progression speed is increased by keeping one of the two keys pushed, in order to reach the desired value faster.

It's possible to select two different type of running by the SET key on the "Electronic Timer Windows":

- Continuous (. . .); - With timer (minute/hour base);

In order to change the timer modality (minutes/hour base), maintain the set key pushed for a few seconds:

- With timer (minutes): "000" - with timer (hours): "h00"

If the time of running is not selected, the stirring is not active.

The display shows the count down and simultaneously the third decimal point turns on to indicate that the count is operative.

At the end, the instrument stops the stirring, the acoustic signal warns the end of the time and the display shows "End".

4. Maintenance

No routine or extraordinary maintenance is necessary apart from periodically cleaning the unit as described in this manual. In compliance with the product guarantee law, repairs to our units must be carried out in our factory, unless previously agreed otherwise with local distributors. The instrument must be transported in its original packaging and any indications present on the original packaging must be followed (e.g., palletized). It is the responsibility of the user, to properly decontaminate the unit in case of hazardous substances remaining on the surface or interior of the device. If in doubt about the compatibility of a cleaning or decontamination product, contact the manufacturer or distributor.

4.1 Cleaning

Disconnect the unit from the power supply and use a cloth dampened with an non-inflammable non-aggressive detergent.

5. Technical data

Power supply	100-240 V / 50-60 Hz
Power	19 W
Dimensions (WxHxD)	655x404x296 mm (JLT4) 935x404x296 mm (JLT6)
Weight (with plug)	13.0 Kg (JLT4) 17.0 Kg (JLT6)
Speed range	10 ÷ 300 rpm, 1 rpm step
Time range	0 ÷ 999 minutes / 0 ÷ 99 hours or continuous
Construction material	Epoxy painted metal structure
Environmental ambient temperature	5 - 40 °C (41 - 104°F)
Storage temperature range	-10... + 60 °C (14 - 140°F)
Max humidity	80%
Pollution degree	2

6. Accessories / Spare parts

A00001000	Plastic Becher 1000 ml	A00001001	Glass Becher 1000 ml
A00001002	Transparent plastic Imhoff cone	A00001003	Glass graduated Imhoff cone
A00001004	Stand for 2 Imhoff cones	10000235	Foot

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
We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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EN 61010-2-051

EN 61326-1

and satisfies the essential requirements of the following directives:

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9. Declaration of conformity ^{UK}CA

We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following regulations:

S.I. 2016/1101 Electrical Equipment (Safety) Regulations 2016

S.I. 2016/1091 Electromagnetic Compatibility Regulations 2016

according to the relevant designated standards:

EN 61010-1 Safety requirements for electrical equipment for measurement, control, and laboratory use
Part 1: General requirements

EN 61010-2-051 Particular requirements for laboratory equipment for mixing and stirring

EN 61326-1 Electrical equipment for measurement, control and laboratory use - EMC requirements -
Part 1: General requirements

and satisfies the essential requirements of regulations:

S.I. 2008/1597 Supply of Machinery (Safety) Regulations 2008

S.I. 2012/3032 Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic
Equipment Regulations 2012

S.I. 2013/3113 Waste Electrical and Electronic Equipment Regulations 2013

plus modifications.

Thank you for having chosen VELP!

Established in 1983, VELP is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. VELP works according to **ISO 9001**, **ISO14001** and **ISO 45001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

Analytical instruments

Elemental Analyzers
Digestion Units
Distillation Units
Solvent Extractors
Fiber Analyzers
Dietary Fiber Analyzers
Oxidation Stability Reactor
Consumables

Laboratory Equipment

Magnetic Stirrers
Heating Magnetic Stirrers
Heating Plates
Overhead stirrers
Vortex mixers
Dispersers
COD Thermoreactors
BOD and Respirometers
Cooled Incubators
Flocculators
Overhead Shakers
Turbidimeter
Open Circulating Baths
Pumps

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VELP Scientifica Srl