Instruction Manual
Manuale di istruzioni
Manuel d'instructions
Manual de instrucciones
Bedienungsanleitung
指导手册

MULTI-HS 6 Multiposition Digital Hotplate Stirrer

F20500530, F20510530

MULTI-HS 15 Multiposition Digital Hotplate Stirrer

F20500531, F20510531

General Information / Informazioni Generali / Informations Générales / Información General / Allgemeine Hinweise/ 基本信息

Before using the unit, please read the following instruction manual carefully.
Prima dell'utilizzo dello strumento si raccomanda di leggere attentamente il seguente manuale operativo.
Avant d'utiliser l'instrument, il est recommandé de lire attentivement le présent manuel d'instructions.
Antes de utilizar el instrumento, le recomendamos que lea con atención el siguiente manual de funcionamiento.
Bitte lesen Sie vor Inbetriebnahme des Geräts diese Bedienungsanleitung sorgfältig durch
使用前请仔细阅读以下使用说明书。

Caution, hot surface! / Attenzione, superficie calda! / Attention, surface chaude! / Prudencia, superficie caliente! / Vorsicht, heiße Oberfläche! / 注意：加热面高温

Do not dispose of this equipment as urban waste, in accordance with EEC directive 2002/96/CE.
Non smaltire l'apparecchiatura come rifiuto urbano, secondo quanto previsto dalla Direttiva 2002/96/CE.
Ne pas recycler l'appareil comme déchet solide urbain, conformément à la Directive 2002/96/CE.
No tirar el aparato en los desechos urbanos, como exige la Directiva 2002/96/CE.
Dieses Gerät unterliegt der Richtlinie 2002/96/EG und darf nicht mit dem normalen Hausmüll entsorgt werden.
按照 EEC 指令 2002/96/CE，不要将该设备作为城市垃圾处理。

The product can be used with flammable liquids / Il prodotto può essere utilizzato con liquidi infiammabili / Le produit peut être utilisé avec des liquides inflammables / El producto puede utilizarse con líquidos inflamables / Das Produkt kann mit brennbaren Flüssigkeiten verwendet werden. / 本产品可与易燃液体混合使用。

This unit must be used for laboratory applications indoor only. The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions. If the product is used in a not specified way by the manufacturer or with non specified accessories, product's safety may be compromised.

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This unit has been designed and manufactured in compliance with the following standards:

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Do not use with explosive or dangerous materials for which the equipment is not designed. The stirrer must not be used in explosive atmospheres, in bain-marie or to stir combustible liquids that have a low combustion temperature. The product is intended for use with very small quantities of flammable liquids or flammable liquids that have a fire point higher than 625°C and a flash point higher than 600°C.

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It is responsibility of the user appropriately decontaminate the instrument in case of dangerous substances fall on or in it accordingly to the safety datasheet of substances used and to the current laboratories' safety standards. It is not possible to decontaminate the product under steam.

It is also responsibility of the user to use substances for cleaning or decontaminating which do not react with internal parts of the instrument or with the material contained in it. In case of doubts on the compatibility of a cleaning solution, contact the manufacturer or local distributor.

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
1. Introduction

EN

Digital multi-position magnetic stirrer with heating available in 6 and 15 position. The MULTI-HS feature exclusive feature such as timer, safety lock, brushless motor to guarantee maximum performance and safety.

- 1 Aluminum alloy plate with ceramic coating
- 2 Main switch
- 3 Speed control knob
- 4 Temperature control knob
- 5 Display

- A Heating ON led
- B Temperature set point led
- C Timer led
- D Autoreverse led
- E Timer/autoreverse/lock button

The Aluminum alloy heating plate with ceramic coating (1) ensures:

- Optimum heat distribution and a high specific power thanks to the circular configuration
- Temperature homogeneity
- High resistance to thermal stress and thermal shock
- High resistance to chemical agents and surface abrasion

Magnetic stirring is generated by neodymium magnets, driven by brushless motors which offer a virtually unlimited duration.

Size and shape of the magnetic stir bar determines the stirring efficiency at any given speed.

The stirring bar code which satisfies most applications for MULTI-HS 6 is the stirring bar A00001056 (6x35 mm); for MULTI-HS 15 is the A00001057 (6x20 mm)

NOTE: The vessel must be made of a suitable material to withstand the foreseen temperature.

NOTE: Using the heating plate at high temperatures may cause discoloring. This does not alter the thermal, mechanical and chemical resistance of the plate in any way.

2. Assembly and installation

- Unpacking
 - Check the integrity of the unit after unpacking.
- The box includes
 - MULTI-HS 6/15 Multiposition Digital Hotplate Stirrer
 - Power supply cord
 - Instruction manual
 - Magnetic stir bar Ø 6x35mm (A00001056)
- First installation
 - Place the unit on non-flammable surface
 - Make sure that the values on the rating plate, correspond to those of the power supply
 - Move the main switch (2) to the OFF position
 - Ensure that the socket provided with grounding is accordant to current safety norms and easy to reach. Use only the cable provided with the instrument.
 - Insert the mains power cable into the socket

NOTE: the main cable must remain far away from the hot plate. It can be substituted only by main cables with same features (T=90°C, connector C15).

3. Operating controls

Commissioning	➤ Switch on the instrument through the main switch (2) ➤ Display shows software version, the last set point values and OFF on both displays (if "Start mode" is set on Stop, see chapter 5)
Stirring	➤ Adjust motors speed set point by turning the speed control knob (3) ➤ Click the knob to start stirring ➤ The speed increases until set point achievement ➤ A microprocessor ensures constant speed even when the viscosity changes (counter-reaction) ➤ Switch off the stirring by clicking the knob
Heating	➤ Adjust heating plate temperature set point by turning the temperature control knob (4) ➤ During all the time of set point adjusting, the temperature set point led (B) is switched on ➤ Click the knob to start heating ➤ Display shows set point temperature. ➤ While the instrument is heating, the Heating ON led (A) turns on ➤ Switch off the heating by clicking the knob ➤ When heating function is off and the heating plate temperature exceeds 50 °C, the message "Hot" is shown on the temperature display until temperature falls below 50 °C NOTE: this warning is not active if the instrument is not powered. NOTE: in case of black-out, once the power is back the device restarts in the set restarting mode (see chapter 5 "Strt mode")
Timer	➤ To access the timer function, click the Timer/autoreverse/lock button (E) ➤ On the displays t l n E h h.n n are visualized ➤ Set the timer turning the speed control knob until the desired value and click to confirm ➤ Timer led (C) turns on ➤ Timer countdown starts when one of the main functions is on or when temperature reaches set point value (see chapter 5, "Time Strt") ➤ When countdown is active, timer led flashes ➤ Esc timer menu: a) click twice timer/autoreverse button b) wait 5 seconds without touching any knob ➤ When timer countdown ends, E n d t l n e is shown until operator action. It's possible stirring continues if the parameter "time Sped" is ON (see chapter 5)
Autoreverse	➤ To access autoreverse function click twice the Timer/autoreverse/lock button (D) ➤ On the displays A.rEu n n. s are visualized ➤ Set the autoreverse turning the speed control knob until the desired value and click to confirm ➤ The autoreverse led (D) turns on ➤ Autoreverse starts when stirring is on
Lock	➤ Holding the Timer/autoreverse/lock button (E) for 3 seconds, the instrument will lock it's settings during operations. ➤ Unlock the control panel by holding the Timer/autoreverse/lock button (E) for 3 seconds. ➤ If other buttons are clicked while the instrument is locked, the LEDs around the lock button will blink for few seconds.

4. Setting mode

Press both knobs for 3 seconds to enter into the setting mode when heating and stirring are off.

Clicking the Timer/autoreverse button (F) it is possible to pass from one of the following parameter to the next ones:

Parameter shown		Default value	Range	Description
Display 1	Display 2			
TENP	LINt	120	OFF - 50 - 120	Set full scale temperature: ➤ It limits the maximum value of temperature set point for the heating plate* (from 50 to 120°C step 10°C) ➤ If "OFF" is set, it is not possible to perform any action with the temperature control knob
SPEd	LINt	1000	100 - 1000	Set full scale speed: ➤ It limits the maximum value of speed set point for the motor (from 100 to 1000 rpm, step 100 rpm)
Strt	NOdE	StOP	StOP - run	Set starting mode: ➤ Stop: when the instrument is switched on through the main switch, OFF is visualized on both displays ➤ Run: when the instrument is switched on, it restarts to work with the last set points set
tINE	Strt	SnAP	SnAP - SEtP	Set start timer countdown: ➤ Snap: timer countdown starts as soon as one main function (Temperature/Speed) is on ➤ SetP: timer countdown starts when temperature set point is reached
tINE	SPEd	OFF	On - OFF	Continue stirring after timer end: ➤ If "on", stirring continues after timer end
ANb.t	CAL	0.0	-10.0 ÷ 10.0	Heating plate probe calibration: ➤ It allows plate ambient temperature alignment to a reference thermometer.
COun	NOtO			Running time. Days (24H): ➤ Motors and plate working time are displayed in hours until 9999 hours (around 416 days). After, they are shown in days starting from 417 and a dot turns on as following, to differentiate days from hours . 4 1 7
COun	HEAt			
rSEt		nO	YES - nO	Reset: ➤ It allows to restore default setting parameters

Click the stirring control knob to visualize the set value. Turn the knob to set the new value and click to confirm.

To esc setting mode wait 10 seconds without touching any knob or push both knobs at the same time.

5. Error messages

When the display shows an error message, the stirring and heating functions stop automatically. To remove the error message, disconnect the instrument from the power supply.

Error code	Cause
AL1	Overtemperature
AL2	Excessive heating time

If an error message appears on the display, please contact VELP Scientifica's technical service department.

6. Maintenance and cleaning

Maintenance	➤ No routine or extraordinary maintenance is necessary; ➤ Repairs must be carried out by authorized Velp personnel only; ➤ Instrument must be transported in its original packaging any indications present on the original packaging must be followed (e.g. palletized); ➤ It is the responsibility of the user, to properly decontaminate the unit in case of hazardous substances remaining on the surface or interior of the device. If in doubt about the compatibility of a cleaning or decontamination product, contact the manufacturer or distributor.
Cleaning	➤ Disconnect the unit from the power supply and use a cloth dampened with a non-inflammable non-aggressive detergent.

7. Technical data

	Models	F20500530, F20500531	F20510530, F20510531
General features	Power supply	230 V / 50-60 Hz (+/-10%)	115V / 60 Hz (+/-10%)
	Dimensions (WxHxD)	240x65x550 mm (9.4x2.6x21.7 in)	240x65x550 mm (9.4x2.6x21.7 in)
	Weight	4.8 Kg (10.5 lb)	4.8 Kg (10.5 lb)
	Power input	620W, 2,7A	620W, 5,4A
	Construction material (body)	Aluminium/stain and steel	Aluminium/stain and steel
	Working in continuous	Admitted	Admitted
	Settable restart modality	Stop or work	Stop or work
	Maximum load on the plate	25kg	25kg
	Noisiness	<< 80 dBa	<< 80 dBa
	Environmental temperature admitted	+5...+40 °C	+5...+40 °C
	Storage temperature admitted	-10...+60 °C	-10...+60 °C
	Max humidity	80%	80%
	Overvoltage category	II	II
	Pollution degree CEI EN61010-1	2	2
	Max altitude	2000 m	2000 m
Heating plate	Heating output	600 W	600 W
	Heating plate dimensions	432 x 223 mm	432 x 223 mm
	Programmable temperature range	Room temperature - 120 °C	Room temperature - 120 °C
	Type of temperature control	Digital	Digital
	Temperature selection	1 °C	1 °C
	Hot plate alarm	Over 50°C	Over 50°C
	Overtemperature alarm	Over 140 °C	Over 140 °C
	Construction material (plate)	Aluminium	Aluminium
	Full scale temperature range	From 50°C to 120°C	From 50°C to 120°C
	Safety circuit	Thermostat	Thermostat
Stir	Stirring capacity	6x400 ml or 15x250 ml	6x400 ml or 15x250 ml
	Programmable speed range	100 – 1000 rpm	100 – 1000 rpm
	Motors type	6 or 15 BLDC	6 or 15BLDC
	Speed selection	5 rpm step	5 rpm step
	Autoreverse	From 5s to 99m:59s	From 5s to 99m:59s
	Autoreverse selection	1 s	1 s
	Full scale speed range	From 100 to 1000 rpm	From 100 to 1000 rpm
Timer	Timer	From 1min to 99h:59min	From 1min to 99h:59min
	Timer selection	1 min	1 min
	2 modalities	Immediate start or at temperature set point	Immediate start or at temperature set point
Counters	Motor counter	Hours of work	Hours of work
	Heating plate counter	Hours of work	Hours of work

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We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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DIN EN 61326-1	Electrical equipment for measurement, control and laboratory use. EMC requirements. General requirements
2015/863/EU (RoHS III)	Restriction of the use of certain hazardous substances
2012/19/EU (WEEE)	On waste electrical and electronic equipment
DIN EN 61010-1	Safety Requirements for Electrical Equipment for Measurement, Control, and Laboratory Use - Part 1: General Requirements.
DIN EN61010-2-051	Safety requirements for electrical equipment for measurement, control and laboratory use - Part 2-051: Particular requirements for laboratory equipment for mixing and stirring
DIN EN61010-2-010	Safety requirements for electrical equipment for measurement, control and laboratory use - Part 2-010: Particular requirements for laboratory equipment for the heating of materials

and satisfies the essential requirements of the following directives:

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10. Declaration of conformity ^{UK}CA

We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following regulations:

S.I. 2016/1101 Electrical Equipment (Safety) Regulations 2016

S.I. 2016/1091 Electromagnetic Compatibility Regulations 2016

according to the relevant designated standards:

EN 61010-1 Safety requirements for electrical equipment for measurement, control, and laboratory use
Part 1: General requirements

EN 61010-2-010 Particular requirements for laboratory equipment for the heating of material

EN 61010-2-051 Particular requirements for laboratory equipment for mixing and stirring

EN 61326-1 Electrical equipment for measurement, control and laboratory use - EMC requirements -
Part 1: General requirements

and satisfies the essential requirements of regulations:

S.I. 2008/1597 Supply of Machinery (Safety) Regulations 2008

S.I. 2012/3032 Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic
Equipment Regulations 2012

S.I. 2013/3113 Waste Electrical and Electronic Equipment Regulations 2013

plus modifications.

Thank you for having chosen VELP!

Established in 1983, VELP is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. VELP works according to **ISO 9001**, **ISO14001** and **ISO 45001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

Analytical instruments

Elemental Analyzers
Digestion Units
Distillation Units
Solvent Extractors
Fiber Analyzers
Dietary Fiber Analyzers
Oxidation Stability Reactor
Consumables

Laboratory Equipment

Magnetic Stirrers
Heating Magnetic Stirrers
Heating Plates
Overhead stirrers
Vortex mixers
Dispersers
COD Thermoreactors
BOD and Respirometers
Cooled Incubators
Flocculators
Overhead Shakers
Turbidimeter
Open Circulating Baths
Pumps

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