Instruction Manual
Manuale di istruzioni
Manuel d'instructions
Manual de instrucciones
Bedienungsanleitung
指导手册

OCB Open Circulating Bath

General Information / Informazioni Generali / Informations générales / Información general / Allgemeine Informationen / 一般信息

Before using the unit, please read the following instruction manual carefully.
Prima dell'utilizzo dello strumento si raccomanda di leggere attentamente il seguente manuale operativo.
Avant d'utiliser l'instrument, il est recommandé de lire attentivement le présent manuel d'instructions.
Antes de utilizar el instrumento, le recomendamos que lea con atención el siguiente manual de funcionamiento.
Bitte lesen Sie vor Inbetriebnahme des Geräts diese Bedienungsanleitung sorgfältig durch.
在使用本装置之前，请仔细阅读以下使用说明书。

Warning! Hot surface.
Attenzione superficie calda
Attention surface chaude
Atención: superficie caliente
Achtung, heiße Oberfläche
注意表面高温！

Do not dispose of this equipment as urban waste, in accordance with EEC directive 2002/96/CE.
Non smaltire l'apparecchiatura come rifiuto urbano, secondo quanto previsto dalla Direttiva 2002/96/CE.
Ne pas recycler l'appareil comme déchet solide urbain, conformément à la Directive 2002/96/CE.
No tirar el aparato en los desechos urbanos, como exige la Directiva 2002/96/CE.
Dieses Gerät unterliegt der Richtlinie 2002/96/EG und darf nicht mit dem normalen Hausmüll entsorgt werden.
根据 EEC 指令 2002/96/CE，请不要将本设备作为城市垃圾处理。

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This unit has been designed and manufactured in compliance with the following standards:

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IEC/EN 61010-1

Electrical equipment for laboratory use

UL 61010-1

General requirement - Canadian electrical code

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Contents / Indice / Index / Índice / Inhalt / 内容

1.	INTRODUCTION	8
1.1	FUNCTIONS	8
2.	ASSEMBLY AND INSTALLATION	8
2.1	ELECTRICAL CONNECTION	8
2.2	START-UP	9
2.3	INFORMATION ABOUT CONSTRUCTION MATERIALS	9
3.	OPERATING CONTROLS	9
3.1	SETTING THE OPERATING TEMPERATURE	9
3.2	SETTING THE TIMER	9
3.3	TEMPERATURE ALIGNMENT	9
3.4	UNITS OF MEASUREMENT	10
3.5	SAFETY FEATURES AND WARNING ERRORS	10
4.	MAINTENANCE	10
4.1	CLEANING	10
5.	TECHNICAL DATA	10
6.	ACCESSORIES	10

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10. DECLARATION OF CONFORMITY **UK**
CA 30

The OCB Open Circulating Bath consists of an immersion heating head with electronic temperature control; a circulation pump and a stainless-steel coated heating coil. The temperature can be set to a max. value of 105°C at 1°C intervals. The heating unit must be positioned above the polycarbonate tank containing the heated liquid.

Different types of test tube racks can be used simultaneously: the 5-liter polycarbonate tank holds 2 test tube rack, the 7-liter polycarbonate tank holds 3 test tube racks. The heating fluid is mixed continuously by the heating unit and circulates in the tank containing the test tube racks and samples.

1.1 Functions

Turn the heating unit on using the mains switch (1) on the front panel.

The display (2) shows the temperature of the heating liquid measured using an immersion probe type Pt100 (3).

The three buttons: ↑, SET, ↓, on the front panel control the operating temperature of the heating liquid, the operating time and other optional functions such as the unit of measure for temperature readings (°C or °F); the alignment of the temperature reading against the reference temperature and the display of error messages.


The stainless-steel coated heating coil (4) heats the liquid in the thermostatic tank until the set temperature is reached.

The led (5) lights up to indicate that heating is in progress.

The heating liquid is kept circulating by the pump (6). If, due to evaporation, the minimum-level sensor detects that the heating liquid in the tank has reached the minimum level, the unit stops heating and the error message “H2O” appears on the display. When this error message appears the liquid in the thermostatic tank must be topped-up to the correct level before resetting the alarm by turning the unit off and on again using the mains switch.

The safety thermostat knob on the rear panel can be used to select the maximum safety temperature (°C) of the heating bath. The temperature probe (8) checks the temperature of the liquid constantly and heating stops if the temperature exceeds the set value.

NOTE: The values shown on the safety thermostat knob are approximate since this thermostat has a safety function only.

2. Assembly and installation

Check the integrity of the unit after unpacking. The box includes:

- OCB Heating unit
- Instruction manual
- Power supply cable

After unpacking the instrument place the polycarbonate tank (code A00000090 or A00000091) on the laboratory bench on a flat, non-flammable surface and within reach of an electric socket. Make sure that access to the electric socket is not obstructed and that it is easy to reach. Fill the tank with deionized or distilled water until the optimum level is reached depending on the type of test tube rack or test tubes to be heated. Place the test tube rack/s in the tank.

NOTE: Make sure that no water enters the heating unit. To avoid overfilling it is advisable to leave a 25mm gap between the water level and the top of the tank. The use of distilled water, de-ionized water or special heating liquids is recommended in order to keep the inside of the tank and the components of the heating unit clean. Normal tap water can however be used although this will result in the build-up of deposits and lime-scale which could affect the efficiency of the instrument.

2.1 Electrical connection

Before connecting the instrument to the power supply, make sure that the values on the rating plate correspond to those of the power supply. Connect the heating unit to the power supply using the cable supplied. The unit is earthed.

The electric socket must be earthed. Position the power cable so that it is well away from hot surfaces.

NOTE: The mains switch disconnects the instrument from the mains supply, reducing unnecessary consumption.

2.2 Start-up

Turn the heating unit on using the luminous mains switch on the front panel. The temperature set-point programmed the last time the unit was used remains active; the default setting on new instruments is 40.0°C.

2.3 Information about construction materials

Structure Stainless Steel
Heating element Stainless Steel

Paint PET
Thermostatting Bath PC

Frontal panel PET
Test tube rack PC

3. Operating controls

When the heating unit is turned-on the display shows the software version installed followed by the last temperature set-point which flashes for approx. 5 seconds.

The display then shows the main window which indicates the temperature reading of the Pt100 probe.

NOTE: The unit automatically sets the heating bath temperature to the value of the last temperature set-point used. The water-circulation pump remains active for the whole time that the instrument is in use.

3.1 Setting the operating temperature

From the main window press the middle key  to display the last temperature set-point used.

With the set-point flashing use the arrow keys  and , to increase or decrease the value to the temperature required. Operating temperatures can be set within a range from room temperature up to 105°C. The setting is recorded after approx. 5" after which the display returns to the main window showing the current temperature of the heating bath.

The new operating temperature will remain in use even after the instrument has been turned off and will be the temperature set-point the next time the instrument is used.

NOTE: The maximum temperature reached during the current work-cycle can be viewed by keeping the key  pressed. The maximum temperature value is deleted from memory each time the instrument is turned off.

3.2 Setting the timer

From the main window press the key  twice to view the operating time. Operating times can be set from 00:00 to 24:59 (h:min.). The default value is 0:00.

With the default time flashing use the keys  and  to increase or decrease the hour value.

To set the minutes value press the key  or wait for 4 seconds.

With the value flashing use the keys  and  to increase or decrease the minutes value. The setting will automatically be memorized after approx. 5 seconds.

If no keys are pressed for 5 seconds the display automatically records the time setting and returns to the main window.

Count-down starts immediately and is indicated on the display by the flashing arrows.

At the end of the operating time the instrument automatically turns the heating unit off and the message "End" appears on the display. The water circulation pump remains active.


To start a new work-cycle press the key . The operating time returns to the default value

(0:00) and the maximum temperature reached can be viewed by pressing .

The operating time can be changed during the work-cycle by following the same procedure indicated above even if an operating time has already be programmed.

With an operating time setting of 0:00 count-down does not start and the instrument runs in continuous mode.

NOTE: To view the operating time remaining press the key  from the main window. The remaining time automatically returns to zero when the instrument is turned off.

3.3 Temperature alignment


The temperature control microprocessor on the main board automatically calibrates the Pt100 probe therefore no calibration procedure is necessary. It is however possible to enter an offset value of from -9.9 to +9.9°C.

Temperature alignment from -9.9 to +9.9°C Default value 0.0°C
from -17 to +17°F Default value 0.0°F

Turn the unit on keeping the two keys  and  pressed. The display reads:

Use the arrow keys  and , to increase or decrease the offset value for "Temperature alignment":

If no keys are pressed for 5 seconds the display shows the software version followed by the set-point and, after a further 5 seconds, the temperature reading taken by the probe. The offset value entered is recorded automatically.

3.4 Units of measurement

The measurement units for temperature readings can be set to °C or °F.

Turn the unit on keeping the two keys  and  pressed. The display shows the unit of measure currently in use.

Use the arrow keys  and , to switch from °C to °F and vice versa.

If no keys are pressed for 5 seconds the display shows the software version followed by the set-point and, after a further 5 seconds, the temperature reading of the probe in the unit of measure selected is shown.

°F
58.0

3.5 Safety features and warning errors

The instrument is fitted with safety features and controls capable of detecting the following errors:

- Probe off-line;
- Faulty probe (failure or short circuit);
- Temperature not within the permitted range (5...+105°C).
In this case heating stops immediately and the message "Err" appears on the display.
- Low water level. Heating stops and the message "H2O" appears on the display

In the case of "Err" or "H2O", remove the fault that caused the alarm and reset the instrument by turning it off and on again from the mains switch.

4. Maintenance

No routine or extraordinary maintenance is necessary apart from periodically cleaning the unit as described in this manual. In compliance with the product guarantee law, repairs to our units must be carried out in our factory, unless previously agreed otherwise with local distributors.

It is the responsibility of the user, to properly decontaminate the unit in case of hazardous substances remaining on the surface or interior of the device. If in doubt about the compatibility of a cleaning or decontamination product, contact the manufacturer or distributor.

4.1 Cleaning


Disconnect the unit from the power supply and use a cloth dampened with a mild, non-flammable detergent.

5. Technical data

Code	F40300240	F40310240
Admitted power supply	230 V / 50–60 Hz	115 V / 50–60 Hz
Power	900 W	
Dimensions (WxHxD) without tank	136x275x180 mm (5,35x10,8x7 in)	
Weight	2.85 Kg (6.28 lb)	
Temperature setting range	5 ...105°C / 41...221°F	
Display resolution	0,1°C	
Temperature stability	0,2°C	
Temperature probe	Pt100	
Safety thermostat	programmable from 0 to 110°C	
Construction material	Steel	
Environmental ambient temperature	5 - 40°C (41 - 104°F)	
Storage temperature range	-10... + 60°C (14 - 140°F)	
Max humidity	85%	
Pollution degree	2	

6. Accessories

- A00000090 Polycarbonate tank 5 liters
(100°C, 2 racks can be supported)
- A00000091 Polycarbonate tank 7 liters
(100°C, 3 racks can be supported)
- A00000092 20-position Test tube rack Ø 18mm
- A00000093 36-position Test tube rack for Eppendorf Ø 11mm
- A00000094 7-position Test tube rack Ø 30mm

		
A00000090	A00000091	
		
A00000092	A00000093	A00000094

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
8. Warranty / Garanzia / Garantie / Garantía / Garantie / 担保

The unit is guaranteed against production defects for **25 months** from our invoice date. In accordance with this guarantee VELP SCIENTIFICA undertakes to repair any instruments resulting as faulty due to the quality of the materials used or poor workmanship. Faults arising due to inexpert handling/use or carelessness will not be replaced or repaired under warranty.

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address / indirizzo / adresse / dirección / Adresse / 地址

Via Stazione, 16, 20865 USMATE (MB), Italy

under our responsibility declare that the product is manufactured in conformity with the following standards:

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2015/863/EU (RoHS III)

2012/19/EU (WEEE)

plus modifications and that the documents listed in annex I are available at VELP's offices as foreseen by the machinery directive.

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10. Declaration of conformity ^{UK} _{CA}

We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following regulations:

S.I. 2016/1101 Electrical Equipment (Safety) Regulations 2016

S.I. 2016/1091 Electromagnetic Compatibility Regulations 2016

according to the relevant designated standards:

EN 61010-1 Safety requirements for electrical equipment for measurement, control, and laboratory use
Part 1: General requirements

EN 61010-2-010 Particular requirements for laboratory equipment for the heating of material

EN 61326-1 Electrical equipment for measurement, control and laboratory use - EMC requirements -
Part 1: General requirements

and satisfies the essential requirements of regulations:

S.I. 2008/1597 Supply of Machinery (Safety) Regulations 2008

S.I. 2012/3032 Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic
Equipment Regulations 2012

S.I. 2013/3113 Waste Electrical and Electronic Equipment Regulations 2013

plus modifications.

Thank you for having chosen VELP!

Established in 1983, VELP is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. VELP works according to **ISO 9001**, **ISO14001** and **OHSAS 18001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

Analytical instruments

Elemental Analyzers
Digestion Units
Distillation Units
Solvent Extractors
Fiber Analyzers
Dietary Fiber Analyzers
Oxidation Stability Reactor
Consumables

Laboratory Equipment

Magnetic Stirrers
Heating Magnetic Stirrers
Heating Plates
Overhead stirrers
Vortex mixers
Homogenizers
COD Thermoreactors
BOD and Respirometers
Cooled Incubators
Flocculators
Overhead Shakers
Turbidimeter
Radiation Detector
Open Circulating Baths
Pumps

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