**Instruction Manual**  
**Manuale di istruzioni**  
**Manuel d'instructions**  
**Manual de instrucciones**  
**Bedienungsanleitung**  
**说明书**

## RC2 Heating Plate, 2-Position

F20700430, F20710430

**General Information / Informazioni Generali / Informations Générales / Información General /  
Allgemeine Hinweise / 一般信息**

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**This unit has been designed and manufactured in compliance with the following standards:**

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Safety requirements for electrical equipment for measurement, control and for laboratory use

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**IEC/EN 61010-1  
IEC/EN 61010-2-051**

Electrical equipment for laboratory use

**UL 61010-1**

General requirement - Canadian electrical code

**CAN/CSA-C22.2 No.61010-1**

VELP reserves the right to modify the characteristics of its products with the aim to constantly improving their quality.

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The plug disconnects the instrument. Therefore, place the instrument where it can be quickly disconnected.

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The heated solution may release toxic, dangerous or poisonous gases. Adequate safety measures must be taken, in accordance with the safety regulations in force, including the presence of hood and personal protective equipment (masks, gloves, goggles, etc.).

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Do not use with explosive or dangerous materials for which the equipment is not designed. The stirrer must not be used in explosive atmospheres, in bain-marie or to stir combustible liquids that have a low combustion temperature. The product is intended for use with very small quantities of flammable liquids or flammable liquids that have a fire point higher than 625°C and a flash point higher than 600°C.

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It is responsibility of the user appropriately decontaminate the instrument in case of dangerous substances fall on or in it.

It is also responsibility of the user to use safety substances for cleaning or decontaminating, which do not react with internal parts of the instrument or with the material contained in it. In case of doubts on the compatibility of a cleaning solution, contact the manufacturer or local distributor.

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## 1. Introduction

The RC2 heating plate is used in the laboratory for general use and for all those applications that require a precise regulation for heating liquids.

The Aluminum alloy heating plate ensures:

- Optimum heat distribution and a high specific power thanks to the circular configuration
- Temperature homogeneity
- High resistance to thermal stress and thermal shock

## 2. Assembly and installation

Check the integrity of the unit after unpacking. The box includes:

- RC2 Heating Plate
- Instruction manual

### 2.1 Electrical connections

After having unpacked the instrument, place the unit on the laboratory bench.

Before connecting the instrument to the power supply, make sure that the values on the rating plate correspond to those of the power supply.

Ensure that the socket and the relative cut-off device conform to current safety norms and are easy to reach.

### 2.2 Start-up

Rotate the temperature knobs completely to the left. Then, set the temperature by turning the dedicated knobs.

## 3. Operating controls

### TEMPERATURE KNOB

To select the temperature of the heating plate turn the corresponding knob (Heating). The temperature can be programmed between 50 and 370 °C.

### ON-OFF SWITCH

The on-off switch turns the unit on and off. If the switch is in the "OFF" position the unit is off; if the switch is in the "ON" position the unit is powered on. The main switch completely disconnect the instrument from the power supply when not used, thus reducing energy waste.

### 3.1 Error message

The unit is fitted with safety devices which cut-off the power supply to the heating plate in the case of malfunctions.

#### Display

Flashing Heating LED (once every second)

Flashing Heating LED (2 times every second)

Flashing Heating LED (once every 3 seconds)

#### Cause

Overheated heating plate or thermocouple open

Fault in the temperature reading circuit

Thermocouple wires reversed or potentiometer board not connected correctly

Should any of the above occur, please contact your nearest VELP Scientifica service centre.

## 4. Maintenance

No routine or extraordinary maintenance is necessary apart from periodically cleaning the unit as described in this manual. In compliance with the product guarantee law, replace the fuses and repairs to our units must be carried out in our factory, unless previously agreed otherwise with local distributors and using Velp spare parts only.

The instrument must be transported in its original packaging and any indication present on the original packaging must be followed (e.g. palletized).

It is the responsibility of the user, to properly decontaminate the unit in case of hazardous substances remaining on the surface or interior of the device. If in doubt about the compatibility of a cleaning or decontamination product, contact the manufacturer or distributor.

### 4.1 Cleaning

Disconnect the unit from the power supply and use a cloth dampened with an non-inflammable non-aggressive detergent.

## 5. Technical data

Velp Code	<b>F20700430</b>	<b>F20710430</b>
Power supply	230V/50-60Hz +/-10%	115V/50-60Hz +/-10%
Dimensions WxHxD	340x107x246 mm	340x107x246 mm
Weight	3.1 Kg	3.1 Kg
Overall power	600W x 2 (1200W)	600W x 2 (1200W)
Construction material	Painted aluminum	Painted aluminum
Diameter of the heating plate	135mm x 2	135mm x 2
Programmable temp. range	50 – 370 °C	50 – 370 °C
Type of temperature control	Analog	Analog
Overtemperature protection	Yes	Yes
Room temperature range	+5...+40 °C	+5...+40 °C
Storage temperature range	-10...+60 °C	-10...+60 °C
Max humidity	80%	80%
Level of electrical protection	IP 42	IP 42
Pollution degree	2	2
Overvoltage category	II	II
Noise level	<< 80 dBa	<< 80 dBa
Max altitude	4000 m	4000 m

## 6. Accessories / Spare parts

A00001069	Support rod
A00000351	Handle for AluBlock removal
A00000382	Extension for support rod ARE – AREX
A00000375	Hemispheric bowl 5L flasks, plate 135
A00000374	Hemispheric bowl 3L flasks, plate 135
A00000373	Hemispheric bowl 25ml flasks, plate 135
A00000346	PTFE Safety cover for bowl 1000 ml
A00000345	PTFE Safety cover for bowl 500 ml
A00000344	PTFE Safety cover for bowl 250 ml
A00000343	PTFE Safety cover for bowl 100 ml
A00000342	PTFE Safety cover for bowl 50 ml
A00000341	MonoAluBlock, 40 pos. Ø12 x h 14 mm
A00000340	MonoAluBlock, 17 pos. Ø28 x h 43 mm
A00000339	MonoAluBlock, 17 pos. Ø28 x h 30 mm
A00000338	MonoAluBlock, 17 pos. Ø28 x h 24 mm
A00000337	MultiAluBlockTM, 11 pos. Ø12 x h 14 mm
A00000334	Hemispheric bowl 1000ml flasks, plate 135
A00000333	Hemispheric bowl 500ml flasks, plate 135
A00000332	Hemispheric bowl 250ml flasks, plate 135
A00000331	Hemispheric bowl 100ml flasks, plate 135
A00000330	Hemispheric bowl 50ml flasks, plate 135
A00000327	MultiAluBlockTM, 4 pos. Ø21 x h31 mm
A00000325	MultiAluBlockTM, 4 pos. Ø28 x h 30 mm
A00000324	MultiAluBlockTM, 4 pos. Ø28 x h 43 mm
10000232	Foot 14Dx8H
10005213	Knob 24D blue
10006341	Ceramic retarded fuse 6,3x32mm 8A (for F20700430)
10006304	Ceramic retarded fuse 6,3x32mm 15A (for F20710430)

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
We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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2015/863/EU (RoHS III)

2012/19/EU (WEEE)

and satisfies the essential requirements of the following directives:

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## 9. Declaration of conformity <sup>UK</sup> <sub>CA</sub>

We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following regulations:

S.I. 2016/1101 Electrical Equipment (Safety) Regulations 2016

S.I. 2016/1091 Electromagnetic Compatibility Regulations 2016

according to the relevant designated standards:

EN 61010-1 Safety requirements for electrical equipment for measurement, control, and laboratory use  
Part 1: General requirements

EN 61010-2-010 Particular requirements for laboratory equipment for the heating of material

EN 61010-2-051 Particular requirements for laboratory equipment for mixing and stirring

EN 61326-1 Electrical equipment for measurement, control and laboratory use - EMC requirements -  
Part 1: General requirements

and satisfies the essential requirements of regulations:

S.I. 2008/1597 Supply of Machinery (Safety) Regulations 2008

S.I. 2012/3032 Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic  
Equipment Regulations 2012

S.I. 2013/3113 Waste Electrical and Electronic Equipment Regulations 2013

plus modifications.

## Thank you for having chosen VELP!

Established in 1983, VELP is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. VELP works according to **ISO 9001**, **ISO14001** and **ISO 45001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

### Analytical instruments

Elemental Analyzers  
Digestion Units  
Distillation Units  
Solvent Extractors  
Fiber Analyzers  
Dietary Fiber Analyzers  
Oxidation Stability Reactor  
Consumables

### Laboratory Equipment

Magnetic Stirrers  
Heating Magnetic Stirrers  
Heating Plates  
Overhead stirrers  
Vortex mixers  
Dispersers  
COD Thermoreactors  
BOD and Respirometers  
Cooled Incubators  
Flocculators  
Overhead Shakers  
Turbidimeter  
Open Circulating Baths  
Pumps

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