Instruction Manual
Manuale di istruzioni
Manuel d'instructions
Manual de instrucciones
Bedienungsanleitung
指导手册

RESPIROMETRIC Sensor

F102B0153

General Information / Informazioni Generali / Informations Générales / Información General / Allgemeine Hinweise / 一般信息

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This unit has been designed and manufactured in compliance with the following standards:

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Safety requirements for electrical equipment for measurement, control and for laboratory use

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Declaration of conformity

We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following regulations:
S.I. 2016/1101 Electrical Equipment (Safety) Regulations 2016
S.I. 2016/1091 Electromagnetic Compatibility Regulations 2016
according to the relevant designated standards:
EN 61010-1 Safety requirements for electrical equipment for measurement, control, and laboratory use
Part 1: General requirements
EN 61326-1 Electrical equipment for measurement, control and laboratory use - EMC requirements -
Part 1: General requirements
and satisfies the essential requirements of regulations:
S.I. 2008/1597 Supply of Machinery (Safety) Regulations 2008
S.I. 2012/3032 Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment
Regulations 2012
S.I. 2013/3113 Waste Electrical and Electronic Equipment Regulations 2013
plus
modifications.

Contents / Indice / Index / Índice / Inhalt

DECLARATION OF CONFORMITY / DICHIARAZIONE DI CONFORMITÀ / DÉCLARATION DE CONFORMITÉ
/ DECLARACIÓN DE CONFORMIDAD / KONFORMITÄTSERKLÄRUNG / 符合性声明  2

DECLARATION OF CONFORMITY  2

1. INTRODUCTION 5

2. ASSEMBLY AND INSTALLATION 5

3. OPERATING CONTROLS 5

 3.1 WIRELESS OR NON-WIRELESS MODE 5

 3.2 DESCRIPTION OF KEYS (WIRELESS MODE) 5

 3.3 DESCRIPTION OF KEYS (NON-WIRELESS MODE) 5

 3.4 LAST DATA RECOVERY 5

 3.5 SCALE SELECTION 6

 3.6 ANALYSIS DURATION AND SAMPLING TIME 6

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1. Introduction

BOD measurement with RESPIROMETRIC Sensor is based on the RESPIROMETRIC method that provides a direct measurement of the oxygen consumed by microorganisms from an air or oxygen-enriched environment in a closed atmosphere, at a constant temperature ($20\text{ °C} \pm 0.5$) and with a soft sample stirring.

The measurement of BMP (Biochemical Methane Potential or Biochemical Methanization Potential) represents the anaerobic correspondent of the BOD tests in fact an adequate anaerobic microbial inoculation is used to degrade the test substrate anaerobically. During the test the production of biogas and/or methane is proportional to the quantity of organic substance gradually degraded.

2. Assembly and installation

Upon receipt and after having removed the packaging, please check the integrity of the instrument. The box includes:

- RESPIROMETRIC Sensor
- Alkali holder
- Instruction manual

The RESPIROMETRIC Sensor is delivered with 2 lithium batteries installed: remove the insulating strip from the rear without opening it.

3. Operating controls

The display of the electronic measuring unit is normally turned off reducing the consumption of internal batteries. The touch of a key activates the display which automatically is turned off after three seconds. The key activates different functions if the unit is performing a measurement or if it is in memory code and if the unit is in wireless mode. When the display is turned off and no analysis is performed, RESPIROMETRIC Sensor shows the symbol $\square\square$ by pressing any keys, in wireless mode.

3.1 Wireless or non-wireless mode

RESPIROMETRIC Sensor is already set on the wireless mode. To switch to a non-wireless mode, keep pressed SET when the display is off until "S.E.T." appears. When the three dots disappear:

- Key **SET** Pass to wireless mode ON
 Key **START** Pass to wireless mode OFF

3.2 Description of keys (wireless mode)

Out of a measuring cycle (before starting)

- | | |
|-----------------------------------|---|
| Key SET with display off | Wireless symbol $\square\square$ is shown |
| Key SET with display on | Wireless symbol $\square\square$ is shown |
| Key START with display off | Wireless symbol $\square\square$ is shown |
| Key START with display on | The cycle of measurement is started |

During a measuring cycle (after starting)

- | | |
|-------------------------------|--|
| Key SET | The scale in use is shown |
| Key START | The current measure is shown (also after the 5-day period if Reset is not performed) |
| Key START kept pressed | The memory mode is entered (also after the 5-day period if Reset is not performed) |

Memory mode (it is active during a measuring cycle and also after if Reset is not performed)

3.3 Description of keys (non-wireless mode)

Out of a measuring cycle (before starting)

- | | |
|-----------------------------------|-------------------------------------|
| Key SET with display off | The scale is shown |
| Key SET with display on | The scale is changed |
| Key START with display off | The scale is shown |
| Key START with display on | The cycle of measurement is started |

During a measuring cycle (after starting)

- | | |
|-------------------------------|--|
| Key SET | The scale is shown |
| Key START | The current measure is shown (also after the 5-day period if Reset is not performed) |
| Key START kept pressed | The memory mode is entered (also after the 5-day period if Reset is not performed) |

Memory mode (it is active during a measuring cycle and also after if Reset is not performed)

- | | |
|------------------|---|
| Key SET | The stored value is shown |
| Key START | The order number of stored value is increased |

If no values are stored, the display shows three points after a prolonged pressing of key **START**.

Reset

Keys **SET+START** kept pressed Measurement cycle stopped, stored values deleted and the scale is shown. Before starting a new measuring cycle, it is requested to reset the stored values.

NOTE: before beginning an additional RESPIROMETRIC Sensor, it's necessary to wait that the display shows "An".

3.4 Last data recovery

If during an analysis the Wireless DataBox™ loses some data, it's possible to retrieve the data for the last seven samplings by pressing SET twice: in this way data will be automatically sent to the Wireless DataBox™.

3.5 Scale selection

Scale	Volume of sample
0 ÷ 4000 mg O ₂ /l	25ml
0 ÷ 999 mg O ₂ /l	100 ml
0 ÷ 600 mg O ₂ /l	150 ml
0 ÷ 250 mg O ₂ /l	250 ml
0 ÷ 90 mg O ₂ /l	400 ml
PRE	Refer to the method
BMP	Refer to the method

The scale is related to the expected value. The possibility of measuring diluted samples extends the range. Results are shown as mg O₂/l (ppm) for BOD test and values in mbar for PRE scale.

Wireless mode

Before the analysis set the scale on the RESPIROSoft™

Non-wireless mode (after the reset)

Key **SET** (display on)

The scale is changed

Key **SET** or **START**

The scale is shown

3.6 Analysis duration and sampling time

Wireless mode

for the Wireless Mode please refer to the Wireless DataBox Operating Manual

Non-wireless mode

Respirometric Sensor will perform only 5-day test with 24-hour sampling time for BOD analysis or 30 days using PRE scale.

Wireless mode

Before the analysis set the scale on the RESPIROSoft™

3.7 Alarms

AL1 Calibration error. Restart the analysis. If the problem persists contact the service.

AL2 Analysis start denied. Check on RESPIROSoft™.

AL3 Max. sensors q.ty reached. 50 analyses are already in progress. Wait the end.

AL4 One parameter or more are wrong

AL5 No PC response. Check the distance between the RESPIROMETRIC Sensor and the Wireless DataBox™.

AL6 The analysis is stopped from the PC

AL7 Unknown command result

AL8 Data sending failed. Restart the analysis. If the problem persists contact the service.

AL9 No PC response. Check if the RESPIROSoft™ is active.

AL10 No automatic measure is done

AL11 Wireless DataBox™ not found. Check if the Wireless DataBox™ is connected.

AL12 Wireless DataBox™ connection isn't strong enough

AL13 Analysis start is denied due to low battery level

bAt Low battery (wireless mode OFF). Remove the batteries with the fingers and replace them with new ones paying attention to the polarity indicated inside the instrument.

Err Calibration error (wireless mode OFF).

--- The measured value overcomes the maximum value of the chosen scale.

4. Maintenance

RESPIROMETRIC Sensor warns about the required change of batteries when the unit is resetted and the charge is lower than a minimum threshold. Use a coin lever in the rear part of RESPIROMETRIC Sensor and open it. Batteries are removed using the fingers. Introduce two fresh batteries at the same time paying attention to the polarity. Battery replacement deletes any stored value.

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Thank you for having chosen VELP!

Established in 1983, VELP is today one of the world's leading manufacturers of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. VELP works according to **ISO 9001**, **ISO14001** and **ISO 45001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

Analytical instruments

Elemental Analyzers
Digestion Units
Distillation Units
Solvent Extractors
Fiber Analyzers
Dietary Fiber Analyzers
Oxidation Stability Reactor
Consumables

Laboratory Equipment

Magnetic Stirrers
Heating Magnetic Stirrers
Heating Plates
Overhead stirrers
Vortexers and Shakers
Homogenizers
COD Thermoreactors
BOD and Respirometers
Cooled Incubators
Flocculators
Overhead Shakers
Turbidimeter
Radiation Detector
Open Circulating Baths
Pumps

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