# **Instruction Manual**

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### **TEMPSoft™**

A00000404

#### **General Information / Informazioni Generali / Informations Générales / Información General / Allgemeine Hinweise**

Before using the unit, please read the following instruction manual carefully.

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Do not dispose of this equipment as urban waste, in accordance with EEC directive 2012/19/UE.

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#### **This unit must be used for laboratory applications indoor only.**

The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions. If the product is used in a not specified way by the manufacturer or with non specified accessories, product's safety may be compromised.

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VELP reserves the right to modify the characteristics of its products with the aim to constantly improving their quality.

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The VELP TEMPSoft™ is an extremely intuitive software able to manage the FOC Connect internal temperature directly from the PC.

This software allows to:

- control up to 10 FOC Connect,
- set working ramps (different temperatures and times) with the possibility to set minimum and maximum temperature alarm thresholds,
- monitor the internal temperature with a dedicated graph,
- export data in multiple formats,
- store test data in conformity with GLP (Good Laboratory Practice)

The communication between the PC and the instrument is managed by RS485 cable.

## 2. Installation

### Unpacking

Check the integrity of the unit after unpacking.

The box includes

- USB pen drive with TEMPSoft™ installer
- Connection cable FOC Connect USB-RS485
- Instruction manual

### Installation

- Connect the USB pen drive to one of the USB ports of the PC
- Run the **setup.exe** from the USB pen drive to start the installation
- The installation guide is displayed. Click **Next** and follow the instructions. The TEMPSoft™ icon will be created on the desktop automatically.
- At the end of the procedure, restart the PC to complete the software installation
- Connect the USB RS485 cable to the PC. The PC recognizes the presence of a new hardware device and it starts the driver installation. If requested, specify the path of the **Driver** folder on the USB pen drive

**Note:** in order to install the TEMPSoft™ Software correctly, administrator rights for the PC must be available. The PC on which the TEMPSoft™ Software will be installed must not be overloaded with other software. This is to prevent unstable operation of the TEMPSoft™ Software.

### Starting the program

- Start the TEMPSoft™  
The faster way is double-click the shortcut icon on the desktop.
- Enter the access credentials:  
User            User  
Password    velp
- The main window will appear

**Note:** the User level cannot delete rows from the database and modify methods. To perform these operations, it is necessary to use the follow access credentials:

User            Superuser  
Password    velp01

### 3. Quick start guide

<i>Menu</i>	<i>Submenu</i>	<i>Description</i>
File	New	Create a new analysis database and "Laboratory information" file.
	Open	Load an already existing database.
	Save as	Save an analysis database (including "Laboratory information" file) assigning a file name. If the file name already exists, it can be overwritten.
	Laboratory information	Enter the details of the laboratory concerned. This information will be included in the test report.
	Exit	Exit from the program.
Report	Create report	Create a standard report of the selected analyses.
	Export data	Export data of the selected rows using different format file.
Modify	Filter	Create database filters.
	Standard column	Reset the order of the database columns to the default setting.
	Configure	Configure the sequence of the database columns, show/hide columns.
	Remove	Remove one or more rows from the database.
	Select all	Select all the rows of the database containing a value.
System	Languages	Change the language of the program.
Analytics	Methods	Display, modify, create or delete methods.
Service	Enable communication	Enable the PC-FOC Connect communication
	Disable communication	Disable the PC-FOC Connect communication
	Register New FOC	Register a new FOC Connect on the TEMPSoft
	ID Table	Permits to see the instrument ID/ S/N list
	Acoustic signals	Enable/Disable acoustic signals.
	Login	Modify the software entry level.
	Diagnostic	Perform a diagnosis of the status of the instrument.
	Control panel	For VELP or authorized technical staff only.
Help	About	Show software version and VELP contacts.

### 4. PC-Incubator connection

- Connect the USB RS485 cable to the FOC Connect; then turn on the FOC Connect
- On the TEMPSoft, press **Service / Register New FOC**  
On the showed popup, select "Automatic search" (or the specific port to which the USB RS485 cable is connected); then press **Connect**
- Assign a new numeric ID (different to 100) to the incubator

**Note:** to connect other incubators, it is necessary to use the accessory "Interconnection cable FOC Connect" (part number A00000405). Use this cable to connect the incubators in series, one incubator at the time; then repeat the above procedure.

- Press **Service / Enable communication** to activate the PC-Incubators communication  
On the showed popup, select "Automatic search" (or the specific port to which the USB RS485 cable is connected); then press **Connect**  
The available incubators will be showed in the table on the right part of the software window.

## 5. Work flow

### Method creation

Press **Analytcs / Methods**.

Press **Create** to create a new method, **Modify** to modify an existing method or **Copy** to copy an existing method.

For each method, enter the method name, the ramps number and the temperature and the time for each ramp.

Press **Ok** to confirm.

### Filling the database

Select **File / New** to create a new database or select **File / Open** to load an already existing database.

Fill the database with the following information:

<b>Sample name</b>	The name of the sample to be processed
<b>Instrument S/N</b>	The serial number of the incubator to be used for the analysis
<b>Method</b>	The method to be used for the analysis
<b>Sampling time</b>	The time interval between successive samples
<b>Alarm temperature Min. (optional)</b>	The software will advise if the incubator temperature goes below this value
<b>Alarm temperature Max. (optional)</b>	The software will advise if the incubator temperature goes above this value
<b>Notes (optional)</b>	Notes about the analysis

**Note:** each database can store max 1000 analyses (rows). For each database, it is possible to fill the laboratory information.

### Starting the analysis

Select a filled row on the database. Then, using the right button of the mouse, press **Start**.

In the first column, a green arrow will appear.

Once the analysis has been started, the database columns cannot be modified.

The other columns will be filled automatically by the software.

The graph on the lower part of the software window refers to the selected row.

### End of the analysis

Once the last ramp of the method is completed, the analysis automatically ends.

Furthermore, it is possible to manually stop the analysis. In order to do it, select the row of the analysis in progress; then, using the right button of the mouse, press **Stop**.

At the end of the analysis, the green arrow in the first column will be replaced by:

Analysis ended correctly

Analysis stopped manually

Analysis ended with warnings

### Modify the showed database columns

Press **Modify / Configure**.

It is possible to modify the order of the columns and to hide some columns.

To restore the default database columns, press **Modify / Standard columns**.

### Filtering the database

Press **Modify / Filter**.

It is possible to apply some filters on the database in order to show only the information complying the filters.

Once the filter parameters have been filled, press **Save New** to save the set filter.

Then, press **Activate Filter** to apply it on the database.

**Note:** the filter status ("Filter ON" or "Filter OFF") is showed in the lower part of the software window.

In the same menu, it is possible to update or delete a previously created filter by pressing the relative buttons.

To remove the filter, press **Deactivate Filter**.

### Exporting data

The information about the analysis can be exported in different formats.

Select the desired rows.

Press **Report / Create Report** to create a standard method.

Press **Report / Export Data** to select which columns you want to export and to select the format file (TXT, LIMS, EXCEL).

Press **Generate Report** to confirm.

## 6. Accessories

A00000405 Interconnection cable FOC Connect

## 7. Software installation requirements

<b>Operating system</b>	Windows 7 (SP 1), Windows 8 or Windows 10
<b>Hardware</b>	Core 2 Duo, RAM 4 GB
<b>Available space on the hard disk</b>	100 MB
<b>Recommended graphic resolution</b>	1024 x 768

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**MULTI**

We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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2011/65/EU (RoHS)    2012/19/EU (RAEE)

and satisfies the essential requirements of the following directives:

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## Thank you for having chosen VELP!

Established in 1983, VELP is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. VELP works according to **ISO 9001**, **ISO 14001** and **OHSAS 18001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

### Analytical instruments

Elemental Analyzers  
Digestion Units  
Distillation Units  
Solvent Extractors  
Fiber Analyzers  
Dietary Fiber Analyzers  
Oxidation Stability Reactor  
Consumables

### Laboratory Equipment

Magnetic Stirrers  
Heating Magnetic Stirrers  
Heating Plates  
Overhead stirrers  
Vortex mixers  
Homogenizers  
COD Thermoreactors  
BOD and Respirometers  
Cooled Incubators  
Flocculators  
Overhead Shakers  
Turbidimeter  
Radiation Detector  
Open Circulating Baths  
Pumps

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VELP Scientifica Srl