Instruction Manual Manuale di istruzioni Manuel d'instructions Manual de instrucciones Bedienungsanleitung 指导手册

TX4 Digital IR Vortex Mixer

F202A0270

General Information / Informazioni Generali / Informations Générales / Información General / Allgemeine Hinweise / 一般信息

Before using the unit, please read the following instruction manual carefully.
Prima dell'utilizzo dello strumento si raccomanda di leggere attentamente il seguente manuale operativo.
Avant d'utiliser l'instrument, il est recommandé de lire attentivement le présent manuel d'instructions.
Antes de utilizar el instrumento, le recomendamos que lea con atención el siguiente manual de funcionamiento.
Bitte lesen Sie vor Inbetriebnahme des Geräts diese Bedienungsanleitung sorgfältig durch
在使用本装置之前，请仔细阅读以下使用说明书。

Do not dispose of this equipment as urban waste, in accordance with EEC directive 2002/96/CE.
Non smaltire l'apparecchiatura come rifiuto urbano, secondo quanto previsto dalla Direttiva 2002/96/CE.
Ne pas recycler l'appareil comme déchet solide urbain, conformément à la Directive 2002/96/CE.
No tirar el aparato en los desechos urbanos, como exige la Directiva 2002/96/CE.
Dieses Gerät unterliegt der Richtlinie 2002/96/EG und darf nicht mit dem normalen Hausmüll entsorgt werden.
根据 EEC 指令 2002/96/CE，请不要将本设备作为城市垃圾处理。

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This unit has been designed and manufactured in compliance with the following standards:

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Safety requirements for electrical equipment for measurement, control and for laboratory use

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**IEC/EN 61010-1
IEC/EN61010-2-051**

Electrical equipment for laboratory use

UL 61010-1

General requirement - Canadian electrical code

CAN/CSA-C22.2 No.61010-1

VELP reserves the right to modify the characteristics of its products with the aim to constantly improving their quality.

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Connect the instrument using only the supplied external power supply. The plug disconnects the instrument. Therefore, place the instrument where it can be quickly disconnected.

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Position the instrument on a flat surface, with a distance from the wall of 30 cm (at least).

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Before placing the tube on the stirring cup, make sure that its cap is tightly closed to prevent sample leakage.

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The solution may release toxic, dangerous or poisonous gases. Adequate safety measures must be taken, in accordance with the safety regulations in force, including the presence of hood and personal protective equipment (masks, gloves, goggles, etc.).

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Do not use with explosive and dangerous materials for which the equipment is not designed. The stirrer must not be used in explosive atmospheres.

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Keep hand/arm vibration to a minimum (Directive 2002/44/EC) when operated manually:

- Keep the sample container as vertical as possible.
- Only exert the minimum pressure necessary for the mixing process.
- Use attachments for multiple samples when there is a large number of samples.

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9. DECLARATION OF CONFORMITY  27

Mixing by spinning can be carried out quickly and safely by placing the tube containing the sample on the mixing cup head. The Vortex Mixer with LCD display can be used in automatic (SENSOR) or continuous mode (TIMER): in both cases the mixing speed can be adjusted using the speed selector knob (1).

When the instrument is switch on, it is always set in SENSOR mode. Mixing starts when the test-tube crosses the infrared The instrument is factory set to continuous mode (Timer), mixing starts automatically when the speed selector / settings knob is pressed. To select Sensor mode, rotate the knob to highlight "Sensor" on the display, press the knob to activate the operating mode selected. Mixing starts when the test-tube crosses the infrared sensor field (5) near the mixing cup head (4). The press-on fitting means that the mixing cup head can easily be replaced by any of the various accessories available for test-tubes of different sizes.


Fig. 1

2. Assembly and installation

Check the integrity of the unit after unpacking. The box includes:

- TX4 Digital IR Vortex Mixer with mixing cup head
- EU adapter plug
- Power supply 100-240V/12V no plug
- Instruction manual

2.1 Electrical connections

After having unpacked the instrument, place the unit on the laboratory bench.

Before connecting the instrument to the power supply, make sure that the values on the rating plate correspond to those of the power supply. Connect the unit to the power supply using the transformer supplied.

Ensure that the socket and the relative cut-off device conform to current safety norms and are easy to reach.

2.2 Start-up

Check that the mixing cup head is correctly positioned. Also check that the speed regulation knob is turned to "0 rpm". Use the on-off switch on the front of the unit to turn the unit on (position "ON"). The start-up and the control of the speed are effected with the knob. When switched on the unit will automatically start in the last operating mode set (Timer / Sensor).

2.3 Information about construction materials

Housing	Zinc alloy/Polymer	Feet	PVC
Paint	PET	Platform	PP
Frontal control panel	PET	Foamed part for platforms	PE
Cup head	SEBS		

3. Operating controls

Connect the unit to mains and turn it on using the on-off switch. Select the operating mode using the mode select button on the front of the unit.


Fig. 2


When **Timer** mode is selected, the timer can be set and the mixing cup head runs continuously at the set speed.

When the **Touch** mode is selected, mixing starts automatically when the test-tube crosses the sensor field (infrared optical system) shown in Figure 2.


In both operating modes there is no need for any pressure on the mixing cup head. The sensor field is purpose-designed to detect the presence of any kind of test-tube.

NOTE: in case of black-out, the device will not automatically restart after the break and reconnect the power supply.

SETTINGS KNOB

With the instrument in stand-by, turn the knob to select and modify the operating parameters (rpm, timer, Start and Sensor) (Fig. 3). With the instrument running, turn the knob to increase or decrease the speed from 0 to 3000rpm or from 0 to 800rpm depending on the full-scale setting. Press the speed selector / settings knob to start or stop mixing in Timer mode (Start or Stop highlighted - Fig 4), to select the operating mode (Sensor or Timer highlighted - Fig 5) and/or to modify speed and time parameters (rpm or timer highlighted - Fig 6 and 7). With the timer highlighted (Fig. 7), turn the knob to set the time from 1 second to 999min. and 59 seconds in Timer mode. During use the display will show the countdown and then END (Fig. 9). To deactivate the time keep the speed selector / settings knob pressed for 3 seconds. The display will show --- : --: mixing will be active for an infinite time. When a parameter is highlighted (Fig. 3), keep the knob pressed for 3 seconds to set the max speed (Fig. 8). This function can be used in Timer mode only. Use the 800rpm full scale setting if you are using a platform. When Timer is highlighted (Fig. 5), keep the knob pressed for 3 seconds to set the sensor time from 1 second to 999min. and 59 seconds (Fig. 10), this function can be used in Sensor mode only. The display shows the countdown and mixing stops automatically.

ON-OFF SWITCH

The on-off switch turns the unit on and off. If the switch is in the "OFF" position the unit is off; if the switch is in the "ON" position the unit is on. Always turn the unit off after use.

4. Maintenance

No routine or extraordinary maintenance is necessary apart from periodically cleaning the unit as described in this manual. In compliance with the product guarantee law, repairs to our units must be carried out in our factory, unless previously agreed otherwise with local distributors.

To replace the cup head, pull it upwards and off. The other push-fit accessories can then be inserted.

The instrument must be transported in its original packaging and any indications present on the original packaging must be followed (e.g. palletized).

It is the responsibility of the user, to properly decontaminate the unit in case of hazardous substances remaining on the surface or interior of the device. If in doubt about the compatibility of a cleaning or decontamination product, contact the manufacturer or distributor.

4.1 Cleaning

Disconnect the unit from the power supply and use a cloth dampened with an non-inflammable non-aggressive detergent.

5. Technical data

Voltage	AC 100-240V (±10%); 50-60 Hz
Power input	12 VDC ; 9.2W
Current consumption	0.77A
Dimensions (WxHxD)	150x130x165 mm (5.9x5.1x6.5 in)
Weight	2.7 Kg (6.2 lb)
Speed range	0-3000 rpm or 0-800 rpm
Timer range	1 sec. - 999min. 59sec.
Spin diameter	4.5 mm (0.177 in)
Max. weight on the mixing cup head / Accessories	0.5 kg (1.1 lb)
Permitted operation	Continuous
Operation mode	Timer or Sensor
Level of electrical protection IEC/EN60529+A1	IP 42
Noise level	<< 80 dBa
Admitted temperature	Storage: -10...+60 °C (+14...+140 F°) Working: 5°...+40 °C (+41...+104 °F)
Admitted humidity	Max. 80%
Overvoltage category	II
Pollution degree IEC/EN61010-1	2
Max. altitude	2000 m

6. Accessories / Spare parts

A00000012	Foam stand for 19 microvials 1.5 ml - Eppendorf®
A00000013	Customizable soft foam top
A00000014	Foam stand 5 test tubes Ø16mm
A00000015	Foam stand for microtiter
A00000016	Small rubber supporting plate Ø 50mm
A00000019	Foam stand 4 test tubes Ø29mm

Do not exceed 800 rpm when using one of the above accessories. High speeds may cause the sample to spill.

		
A00000012	A00000013	A00000014
		
A00000015	A00000016	A00000019

10005195	Mixing cup head TX, ZX, RX
10005213	Knob 24D blue

10005880	Foot 28Dx5H antivibration
40001099	Switching VELP 100-240V/12V without plug

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
We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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EN 61010-2-051

EN 61326-1

and satisfies the essential requirements of the following directives:

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9. Declaration of conformity ^{UK}CA

We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following regulations:

S.I. 2016/1101 Electrical Equipment (Safety) Regulations 2016

S.I. 2016/1091 Electromagnetic Compatibility Regulations 2016

according to the relevant designated standards:

EN 61010-1 Safety requirements for electrical equipment for measurement, control, and laboratory use
Part 1: General requirements

EN 61010-2-051 Particular requirements for laboratory equipment for mixing and stirring

EN 61326-1 Electrical equipment for measurement, control and laboratory use - EMC requirements -
Part 1: General requirements

and satisfies the essential requirements of regulations:

S.I. 2008/1597 Supply of Machinery (Safety) Regulations 2008

S.I. 2012/3032 Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic
Equipment Regulations 2012

S.I. 2013/3113 Waste Electrical and Electronic Equipment Regulations 2013

plus modifications.

Thank you for having chosen VELP!

Established in 1983, VELP is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. VELP works according to **ISO 9001**, **ISO14001** and **OHSAS 18001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

Analytical instruments

Elemental Analyzers
Digestion Units
Distillation Units
Solvent Extractors
Fiber Analyzers
Dietary Fiber Analyzers
Oxidation Stability Reactor
Consumables

Laboratory Equipment

Magnetic Stirrers
Heating Magnetic Stirrers
Heating Plates
Overhead stirrers
Vortex mixers
Dispersers
COD Thermoreactors
BOD and Respirometers
Cooled Incubators
Flocculators
Overhead Shakers
Turbidimeter
Open Circulating Baths
Pumps

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