Instruction Manual
Manuale di istruzioni
Manuel d'instructions
Manual de instrucciones
Bedienungsanleitung
指导手册

VTF Digital thermoregulator Fuzzy logic

F208B0063

General Information / Informazioni Generali / Informations Générales / Información General / Allgemeine Hinweise / 一般信息

Before using the unit, please read the following instruction manual carefully.

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Do not dispose of this equipment as urban waste, in accordance with EEC directive 2002/96/CE.

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This unit must be used for laboratory applications only.

The manufacturer declines all responsibility for any use of the unit that does not comply with these instructions.

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This unit has been designed and manufactured in compliance with the following standards:

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Electrical equipment for laboratory use	UL 61010-1
---	-------------------

General requirement - Canadian electrical code	CAN/CSA-C22.2 No.61010-1
--	---------------------------------

VELP reserves the right to modify the characteristics of its products with the aim to constantly improving their quality.
Nell'impegno di migliorare costantemente la qualità dei prodotti, VELP si riserva la facoltà di variarne le caratteristiche.
Dans le but d'améliorer constamment la qualité de ses produits, VELP se réserve le droit d'apporter des modifications aux caractéristiques de ceux-ci. VELP se reserva el derecho de modificar las características de sus productos con el objetivo de mejorar constantemente su calidad. VELP behält sich zum Zwecke der ständigen Verbesserung der Produktqualität das Recht auf Änderung der Geräteeigenschaften vor.
VELP 保留修改其产品特性的权利，目的是不断提高其质量。

Safety Regulations / Norme di Sicurezza / Consignes de Sécurité / Advertencias de Seguridad / Sicherheitshinweise / 安全条例

The plug disconnects the instrument. Therefore, place the instrument where it can be quickly disconnected.
La spina è il mezzo di disconnessione dell'apparecchio. Pertanto, non posizionare l'apparecchio in modo che sia difficile azionare il mezzo di disconnessione. Le bouchon est le moyen de déconnexion de l'appareil. Par conséquent, placer l'appareil où il peut être rapidement débranché. El tapón es el medio de desconexión del dispositivo. No coloque el dispositivo en una forma que es difícil de desconectar. Der Stecker trennt das Gerät. Daher Stellen Sie das Instrument, wo es schnell getrennt werden kann. 该插头可以断开仪器的连接。因此，要把仪器放在可以快速断开的地方。

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1. Introduction

EN

The Thermoregulator is ideal for many applications and meets the most demanding requirements in terms of precision, reliability and flexibility of use thanks to the application of Fuzzy Logic technology. The Fuzzy Logic electronic automatically adapts thermoregulation to the varying factors such as power, load and thermal dispersion specific to each application, optimizing both overheating and oscillations around the temperature set-point. The accuracy and precision of thermoregulation at each end of the scale and whatever the volume being processed, is a fundamental characteristic of the Vertex.

WARNING: if the Thermoregulator is connected to an instrument different than the VELP Scientifica AREX, AREX DIGITAL PRO, AREC.X or PW10, please check the compatibility of the DIN 5 pole connector as illustrated in wiring diagram in this manual.

NOTE: the most precise results are obtained measuring aqueous solutions temperature.

The microprocessor offers various other functions:

- a working time of up to 24 hours and 59 minutes can be set with automatic switch off;
- the maximum sample temperatures reached during the test can be recorded.

The structure is made of a non-scratch technopolymer resistant to chemical agents and offers a high level of IP54.

The instrument has an in-built safety and control circuit which is constantly active and shuts-down the thermoregulator immediately in the following situations: the temperature probe is not connected; the temperature probe is faulty (cut-off or short-circuited), the temperature is out-of-range and/or the probe is not inserted inside the sample.


Fig. 1

- | | |
|--------------------|--------------------|
| 1 Display | 2 Set |
| 3 Minus | 4 Plus |
| 5 Measurement unit | 6 Timer display on |
| 7 Temp. probe | 8 Green "ON" led |
| 9 Clamps | 10 Plug |
| 11 Sliding clamp | |

2. Assembly and installation

Check the integrity of the unit after unpacking.

The box includes:

- Thermoregulator, with power cable for connection to VELP Scientifica AREX, AREX DIGITAL PRO, AREC.X and PW10
- Pt100 temperature probe
- Instruction manual

Connecting the Pt100 probe

Connect the probe to the Thermoregulator as shown on Fig. 2.

NOTE: the mains cable must remain far away from the hot plate.

Installing the VTF

The VTF has an innovative integrated system to simplify installation on the support rod and facilitate the positioning of the Pt100 probe in the most commonly used containers.

The two clamps allow height-regulation whilst a sliding clamp allows horizontal regulation (Fig. 3).


Fig. 2


Fig. 3

Using the VTF with magnetic stirrers

The power cable is suitable for connection to the most common heating magnetic stirrers with a dedicated socket for the remote control of the heating plate such as the VELP Scientifica models AREX, AREX DIGITAL PRO and AREC.X. The socket must supply a tension of between 9 and 15V DC. The integrated installation system allows the instrument to be installed on support rods with a diameter of between 10 and 13 mm.

Other applications

When using the thermoregulator for other applications bear in mind that the transistor output can take a maximum current of 50mA. Where power loads of up to 2200W are required the use of an external power relay is necessary (code A00000001, VELP PW10).

For the thermoregulation of liquids that are not compatible with the construction material of the probe (stainless steel AISI 316), the use of a glass-coated probe is recommended (code A00000003).

When it is necessary to thermoregulate a liquid that is not in close proximity to the VTF, a 1m probe extension cable is available (code A00000002).

3. Operating controls

Install the thermoregulator on the support rod and place the probe in the liquid to be processed. To power the VTF plug the spiral cable into the dedicated socket on the heating magnetic stirrer and turn the stirrer on. The display shows the software version after which it flashes for approximately five seconds showing last set point (Fig. 4).

Fig. 4

The current Pt100 probe temperature reading is then displayed on the main window. Thermoregulation starts automatically based on temperature Set Point selected.

NOTE: the magnetic stirrer always exercises primary control of the heating plate temperature. To control the heating plate temperature using the VTF thermoregulator, set the magnetic stirrer temperature to maximum. The magnetic stirrer temperature control can also be used as a safety thermostat since the heating plate will not exceed the temperature set on the stirrer. In this case, longer heating time will be necessary in order for the sample to reach the VTF set point temperature.

Setting the working temperature

From the main window, press the central  button. The set-up window is displayed (Fig. 4).

Set Point Temperature From -10 to +300 °C Default value: 10 °C, for VTF

When the display is flashing use the  and  keys to select the temperature required. The temperature setting is saved when no keys are touched for approx. 5 seconds (Fig. 4).

If no key is pressed for approx. 5 seconds when the set-up window is flashing, the value shown on the display is automatically saved. After 5 seconds the VTF will evaluate the temperature reading of the probe immersed in the liquid and will proceed to thermoregulate the liquid to the temperature selected.

NOTE: to display the maximum temperature reached since the instrument was turned on, with the main window displayed keep the  key pressed. The maximum temperature reached is deleted from the memory every time the instrument is turned off.

Setting the timer

From the main window press the  key twice to display the working time.

Working time From 00:00 to 24:59 (h:min.) Default value: 0:00 (Fig. 6)

When the display is flashing, use the  and  keys to select the working time required.

To switch from hours to minutes press the  key once or simply wait 4 seconds (Fig. 7).


Fig. 6

Fig. 7

The value is saved when no keys are touched for approx. 5 seconds. If no key is pressed for 5 seconds the display returns the main window and the previously displayed value is saved.

Count-down starts the moment the working time has been saved. The alternate flashing of the two arrows to the left of the display indicates that the VTF is running in timer-mode (Fig. 8).

Fig. 8

At the end of the working time the instrument automatically turns off the thermoregulation, and the display shows "End" (Fig. 9).


Fig. 9

To return to the main window in order to start new work-cycle, press . Working time will return to the default value of 0:00 and the maximum temperature reached can be viewed by pressing . Working time can be changed when VTF is in use even if the instrument is running in timer-mode. If working time is set to 0:00 the instrument will run in continuous mode.

NOTE: to display the time left, with the main window displayed keep the  key pressed.

The time left is automatically reset to zero every time the instrument is turned off.

Temperature alignment

The micro-chip carries out automatic probe calibration therefore a probe calibration procedure is unnecessary. Nevertheless, an offset value of from - 9.9 to + 9.9 °C can be set.

Temperature alignment	From - 9.9 to + 9.9 °C	Default value	0.0 °C
	From - 17 to + 17 °F		0.0 °F

Turn the instrument on by keeping the  and  keys pressed.

Use the  and  keys to modify the “Temperature alignment” parameters (Fig. 10).


Fig. 10

If no key is pressed for 5 seconds the display shows the software version and the default set point value.

If no key is pressed for a further 5 seconds the display shows the probe temperature reading and saves the value.

Unit of measure

Turn the instrument on by keeping the  and  keys pressed.

Use the  and  keys to modify the “Unit of measure” (Fig. 11).

Fig. 11

If no key is pressed for 5 seconds the display shows the software version and default set point value. If no key is pressed for a further 5 seconds the display shows the probe temperature reading in the new unit of measure and saves the value.

Safety devices and alarm signals

The instrument has an in-built safety and control circuit which is constantly active and shuts-down the VTF when an alarm message appears on the display:

- AL .1 = Unconnected probe or faulty probe (cut-off or short circuited) or temperature out of range.
- AL .2 = Slow temperature increase read by the probe. This function is active in case of:
 - Probe temperature < 50°C
 - T set point – T probe > 5°C
- AL .3 = Fast temperature decrease read by the probe.

To reset the thermoregulator turn it off and on again after having found and removed the cause of the alarm.

In order to bypass the AL .2 and AL .3, turn the instrument on by keeping the  and  keys pressed and select AL.oF using  and  keys. Press  to confirm.

4. Start-up

- Install the VTF on the support rod. Make sure that the probe is immersed in the liquid for at least 15 mm and that it isn't in contact with the bottom of the container;
- Make sure the magnetic stirrer is turned on;
- Power the VTF by plugging the spiral cable into the dedicated socket on the heating magnetic stirrer;
- Set the magnetic stirrer temperature to maximum;
- Select the temperature and the working time required;
- To optimize thermoregulation and improve temperature homogeneity, stir the sample gently but continuously.

WARNING

- When working with set points near the temperature of evaporation, check that the probe remains immersed in the sample by at least 10mm throughout the work-cycle.
- The sample being processed is subject to evaporation and for this reason it may not reach the temperature selected.

5. End-of-work operations

To interrupt thermoregulation and leave the VTF thermoregulator powered, turn off heating plate using the knob on the magnetic stirrer. In this way heating stops but the VTF continues to display the temperature.

If the timer has been set, the VTF will automatically stop at the end of the set time and the display will show “End”.

To reactivate thermoregulation press .

At the end of the process, if the VTF has to be left connected to the magnetic stirrer, turn the stirrer off.

6. Maintenance

No routine or extraordinary maintenance is necessary apart from periodically cleaning the unit as described in this manual. In compliance with the product guarantee law, repairs to our units must be carried out in our factory, unless previously agreed otherwise with local distributors. The instrument must be transported in its original packaging and any indications present on the original packaging must be followed (e.g. palletized).

It is the responsibility of the user, to properly decontaminate the unit in case of hazardous substances remaining on the surface or interior of the device. If in doubt about the compatibility of a cleaning or decontamination product, contact the manufacturer or distributor.

Cleaning

Disconnect the unit from the power supply and use a cloth dampened with a mild, non-flammable detergent.

7. Technical data

Power supply	12 V (9÷15)
Power	1.2 W
Dimensions (WxHxD)	75x145x120 mm
Weight	0.250 kg
Display	4 digit LCD, Multi function
Type of thermoregulation	Fuzzy Logic
Connector	5 pole 270° DIN
Level of electrical protection (CEI EN60529)	IP54
Temperature range	0...+ 50 °C
Storage temperature range	- 10...+ 60 °C
Max. humidity	90 %
Construction material	Technopolymer
Temperature setting (interval)	- 10...+ 300 °C 14...572 °F
Resolution (precision)	0.2 °C (± 0.5 °C) 1 °F (± 0.9 °F)
Working time setting	From 0:00 to 24:59 h:min
Continuous mode Remaining time reading	Possible
Type of probe	Pt100, 250 mm lenght
Probe speed τ 90%	5 s
Immersion depth	15 mm minimum

8. Accessories / Spare parts

A00000001	Derivation element PW10
A00000002	Probe extension cable, length 1m
A00000003	Glass probe
A00000280	Clamp with probe support
40000579	Pt100 probe
40000580	Friction assembled

			
A00000001	A00000002	A00000003	A00000280

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
We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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2015/863/EU (RoHS III)

2012/19/EU (WEEE)

and satisfies the essential requirements of the following directives:

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2006/42/EC

2014/30/EU

Plus modifications and the documents required in Annex I of the Machinery Directive are available at our premises.

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plus modifications and that the documents listed in the annex I are available at Velp's offices as per the machinery directive.

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11. Declaration of conformity ^{UK}CA

We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following regulations:

S.I. 2016/1101 Electrical Equipment (Safety) Regulations 2016

S.I. 2016/1091 Electromagnetic Compatibility Regulations 2016

according to the relevant designated standards:

EN 61010-1 Safety requirements for electrical equipment for measurement, control, and laboratory use
Part 1: General requirements

EN 61326-1 Electrical equipment for measurement, control and laboratory use - EMC requirements -
Part 1: General requirements

and satisfies the essential requirements of regulations:

S.I. 2008/1597 Supply of Machinery (Safety) Regulations 2008

S.I. 2012/3032 Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic
Equipment Regulations 2012

S.I. 2013/3113 Waste Electrical and Electronic Equipment Regulations 2013

plus modifications.

Thank you for having chosen VELP!

Established in 1983, VELP is today one of the world's leading manufacturer of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. VELP works according to **ISO 9001**, **ISO14001** and **ISO 45001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

Analytical instruments

Elemental Analyzers
Digestion Units
Distillation Units
Solvent Extractors
Fiber Analyzers
Dietary Fiber Analyzers
Oxidation Stability Reactor
Consumables

Laboratory Equipment

Magnetic Stirrers
Heating Magnetic Stirrers
Heating Plates
Overhead stirrers
Vortex mixers
Dispersers
COD Thermoreactors
BOD and Respirometers
Cooled Incubators
Flocculators
Overhead Shakers
Turbidimeter
Open Circulating Baths
Pumps

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