



Effect of the Addition of Chili Pepper Powder on Vegetable Oils Oxidative Stability

By: Antonella Cavazza¹, Stefania Corti², Claudia Mancinelli², Chiara Bignardi¹, Claudio Corradini¹.

¹ Dipartimento di Chimica, Università degli Studi di Parma, Parco Area delle Scienze 17/A, 43124 Parma, Italy.

² VELP Scientifica srl, via Stazione 16, 20865 Usmate, Monza and Brianza, Italy.

ARTICLE INFO

Article history:

Received: 5 June 2015

Revised: 6 October 2015

Accepted: 9 October 2015

Published online: 20 October 2015

Keywords:

Virgin olive oils

Sunflower oil

Chili Pepper

Capsaicinoids


Pungency

OXITEST

ABSTRACT

Red chili pepper represents an excellent source of functional bioactive compounds, and its use as food additive could be useful for protecting other products from spoilage. Olive oils show a limited shelf-life due to a progressive quality deterioration linked mostly to oxidation reactions providing a loss of quality in terms of organoleptic and healthy properties. In this work, the addition of red chili pepper powder to virgin, extra-virgin olive and sunflower oils was proved to be effective in improving their stability during 12 months of shelf-life. Oil stability was monitored by employing the Oxitest system, based on accelerating the oxidation process. Three pepper powders of different pungency were added to the virgin olive oil samples with the aim of investigating whether capsaicinoids could be responsible for the observed effects. The strongest effect was recorded when employing the less pungent pepper showing that the protective effect could not be attributed to the capsaicinoids content. Besides, a screening on different oils demonstrated that the less stable samples received a stronger protective effect. The evaluation of a possible dose-effect relation was also carried out pointing out that a threshold dose of about 1 % was required to record an effect. The results suggest potential perspectives in the field of food technology, as addition of a small amount of sweet pepper could be proposed to prolong food shelf-life.

© AOCS 2015



VELP Scientifica, via Stazione 16, 20865 Usmate (MB), Italy

Tel: +39 039 628 811 Fax: +39 039 628 8120 www.velp.com