

# Comparison table

## Glass Crucibles



**Glass crucible P0**  
6pcs/box

**Glass crucible P1**  
6pcs/box

**Glass crucible P2**  
6pcs/box

**Glass crucible P3**  
6pcs/box

Code	A00000138	A00000139	A00000140	A00000137
Filter	P0	P1	P2	P3
Nominal pore size	160-250 micron	100-160 micron	40-100 micron	16-40 micron
Equipment	FIWE Advance - FIWE - CSF6	FIWE Advance - FIWE - CSF6	FIWE Advance - FIWE - CSF6	FIWE Advance - FIWE - CSF6
Example of application fields	Special crucible: filtration of very coarse precipitates	Special crucible: filtration of coarse precipitates	Standard crucible: preparatory work with crystalline precipitates	Special crucible: analytical work with medium-fine precipitates