

VELP Scientifica Dispersers



OV 625 Digital



OV 725 Digital

Brushless Motor	No	Yes
Volume range min. (H2O)	0.2 ml	0.2 ml
Volume range max. (H2O)	Up to 2,5 l	Up to 2,5 l
Viscosity max.	5.000 mPa s	5.000 mPa s
Speed min.	3.000 rpm	2.000 rpm
Speed max.	25.000 rpm	25.000 rpm
Speed control	Digital	Digital
Speed display	LED 3 Digit	LED 4 Digit
Timer	No	Yes
Noise without element	75 dB(A)	70 dB(A)
Width mm	92 mm	100 mm
Depth mm	130 mm	235 mm
Height mm	280 mm	280 mm
Weight kg	2,4 kg	3,7 kg
Permissible ambient temperature min.	5 °C	5 °C
Permissible ambient temperature max.	40 °C	40 °C
Permissible relative humidity	80%	80%
Voltage	230 V / 115V 50-60 Hz (+/-10%)	230 V / 115V 60 Hz (+/-10%)
Power consumption	800 W	500 W
Motor rating output	500 W	400 W