



ESP

Magnetic Stirrers

Ultraflat magnetic stirrer with no moving mechanical components



Reliable and Easy-to-use

Flat and light weight magnetic stirrer with no moving mechanical components: it is therefore maintenance-free. The clear surface is ideal for visualization of particles and color changes.

- Speed regulation up to 1100 rpm
- Accepts up to 5 liter flasks
- Ideal for microtitration, microbiology and biochemistry



INSTRUMENT - CODE

ESP

100-240 V / 50-60 Hz

F206A0179

TECHNICAL DATA

CONSTRUCTION MATERIAL	Technopolymer
AMBIENT OPERATING TEMPERATURE	5 - 40 °C
STIRRING SPEED	Up to 1100 rpm
STIRRING VOLUME	Up to 5 L
DIMENSION (WXHXD)	160x33x230 mm (6.3x1.3x9 in)
POWER	5 W
WEIGHT	0.9 Kg (2.0 lb)



We reserve the right to make technical alterations
We do not assume liability for errors in printing, typing or transmission

VELP Official Partner

DESIGNED AND MANUFACTURED IN ITALY



ITALY - HQ
Via Stazione 16
20865 Usmate (MB) Italy
Tel. +39 039 628811
velpitalia@velp.com

INDIA
velpindia@velp.com

USA
155 Keyland Court, Bohemia
NY 11716 - U.S.
Tel. +1 631 573 6002
velpusa@velp.com

CHINA
Xinlong Rd Building 28, Lane 1333
Shanghai city - China
Tel. +86 180 1755 7329
velpchina@velp.com