



MICROSTIRRER

Magnetic Stirrers

Compact magnetic stirrer specially designed for microtitrations



Reliable and Easy to use

Compact and light weight magnetic stirrer with metal structure for visualization of particles and color changes. Specially designed for microtitrations.

- Speed regulation up to 1100 rpm
- Accepts up to 5 liter flasks
- Excellent speed control, even at low speed
- Epoxy-coated steel casing
- Remains cold, even after days of continuous use



INSTRUMENT - CODE

MICROSTIRRER

100-240 V / 50-60 Hz

F203A0161

TECHNICAL DATA

CONSTRUCTION MATERIAL	Epoxy painted metal structure
SPEED CONTROL	Excellent speed control even at low revolutions
STIRRING SPEED	Up to 1100 rpm
STIRRING VOLUME	Up to 5 L
DIMENSION (WXHxD)	120x48x128 mm (4.7x1.9x5 in)
POWER	0.6 W
WEIGHT	0.55 Kg (1.1 lb)



We reserve the right to make technical alterations
We do not assume liability for errors in printing, typing or transmission

VELP Official Partner

DESIGNED AND MANUFACTURED IN ITALY



ITALY - HQ
Via Stazione 16
20865 Usmate (MB) Italy
Tel. +39 039 628811
velpitalia@velp.com

INDIA
velpindia@velp.com

USA
155 Keyland Court, Bohemia
NY 11716 - U.S.
Tel. +1 631 573 6002
velpusa@velp.com

CHINA
Xinlong Rd Building 28, Lane 1333
Shanghai city - China
Tel. +86 180 1755 7329
velpchina@velp.com