



RC 2

Heating Plates

Two-position heating plates with aluminum top



Practical and Easy to use

Heating plates with aluminum top and temperature regulation. Designed for the longest lifespan and to guarantee excellent resistance to chemicals. It accepts different combinations of AluBlocks™ that can maximize temperature control and save handling time.

- Maximum temperature 370 °C
- Accepts flasks of up to 15 L volume
- Uniform heat distribution with coated aluminum alloy heating plate



INSTRUMENT - CODE

RC2	230 V / 50-60 Hz	F20700430
RC2	115 V / 50-60 Hz	F20710430

TECHNICAL DATA

	RC2
HOUSING MATERIAL	Epoxy painted aluminum structure
PLATE MATERIAL	Aluminum alloy coated with a special protection
PLATE DIMENSION	135 mm
MODELS	2-position
TEMPERATURE RANGE	Room temp. - 370 °C
POWER INPUT	1200 W
WEIGHT	3.3 Kg (7.3 lb)
DIMENSIONS (WxHxD)	340x90x246 mm (13.4x3.5x9.7 in)



We reserve the right to make technical alterations
We do not assume liability for errors in printing, typing or transmission

VELP Official Partner

DESIGNED AND MANUFACTURED IN ITALY



ITALY - HQ
Via Stazione 16
20865 Usmate (MB) Italy
Tel. +39 039 628811
velpitalia@velp.com

INDIA
velpindia@velp.com

USA
40, Burt Drive, Unit #1, Deer Park
NY 11729 - U.S.
Tel. +1 631 573 6002
velpusa@velp.com

LATAM
velplatam@velp.com

CHINA
Xinlong Rd Building 28, Lane 1333
Shanghai city - China
Tel. +8621 34500630
velpchina@velp.com